R<sub>x</sub> Only

8086098(11)

2020-04

English

# Veritor™ System

## For Rapid Detection of Respiratory Syncytial Virus (RSV)

CLIA-waived kit configured for testing nasopharyngeal swab samples freshly collected and dispensed directly onto assay test device.

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## For Rapid Detection of Respiratory Syncytial Virus (RSV)

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BD Veritor System RSV Devices	30 devices	Foil pouched device containing one reactive strip. Each strip has one test line of monoclonal antibody specific to RSV viral antigen and murine monoclonal control line antibodies.
RV Reagent D	30 tubes with 400 µL reagent	Detergent with <0.1% sodium azide (preservative).
Flexible minitip flocked swab	30 each	Swab for nasopharyngeal collection.
RSV Positive Control Swab	1 each	RSV Positive Control Swab, RSV antigen (noninfectious cell lysate) with <0.1% sodium azide (preservative).
RSV Negative Control Swab	1 each	RSV Negative Control Swab, (detergent-treated non-infected cells) with <0.1% sodium azide (preservative).

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**Storage and Handling: Kits may be stored at 2–30 °C. DO NOT FREEZE. Reagents and devices must be at room temperature (15–30 °C) when used for testing.**

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

1. The BD Veritor System RSV Kit includes swabs with a flocked tip for nasopharyngeal specimen collection.
2. Insert the swab into one nostril of the patient, reaching the surface of the posterior nasopharynx.
3. Rotate the swab over the surface of the posterior nasopharynx.
4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the BD Veritor System Kit.

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To use Walk Away mode - Connect the AC power adapter to the Analyzer and a power source

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**If a printer is not connected, record the result before removing the assay device.**

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

- Turn the Analyzer on by pressing the blue power button once.  
The display window will briefly show “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.
- When the display window reads:  
“INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”
  - **Double-click the blue power button.**

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

- Insert the test device into the slot on the right side of the Analyzer.

**The test device must remain horizontal to prevent spilling the specimen out of the sample well.**

- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

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**Step 10D: Remove the test device**

- Remove and then discard the test device appropriately. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to Analyze Now mode at the conclusion of each read sequence.
- ✉ If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

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#### **REPORTING OF RESULTS**

- Positive Test** Positive for the presence of RSV antigen. A positive result may occur in the absence of viable virus.
- Negative Test** Negative for the presence of RSV antigen. Infection due to RSV cannot be ruled out because the antigen present in the sample may be below the detection limit of the test. In the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral cell culture or an FDA-cleared RSV molecular assay.
- Control Invalid** Do not report results. Repeat the test.

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1. Remove and discard the cap from the RV Reagent D tube corresponding to the sample to be tested.
  2. Insert the control swab into the tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
  3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
  4. Continue processing the swab according to the Test Procedure for Nasopharyngeal Swabs above, beginning at Step 5.
- If the kit controls do not perform as expected, do not test patient specimens. Contact BD Technical Services at 1.800.638.8663.

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### REFERENCES

- Hall CB, Weinberg GA, Iwane MK, et al., The Burden of Respiratory Syncytial Virus Infection in Young Children. *N Engl J Med* 2009;360:588–98.
- Nair H, Nokes DJ, Gessner BD, et al., Global burden of acute lower respiratory infections due to respiratory syncytial virus in young children: a systematic review and meta-analysis. *Lancet* 2010;375:1545–55.
- Falsey AR and EE Walsh, Respiratory syncytial virus infection in adults. *Clin Microbiol Rev* 2000;13:371–84.
- Murata Y and AR Falsey, RSV Infection in Elderly Adults, In: Patricia Cane, Editor(s), *Perspectives in Medical Virology*, 2006, Elsevier, Volume 14, Pages 163–82.
- Crowcroft NS, Cutts F and Zambon MC, Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. *Commun Dis Public Health*. 1999;2:234–41.
- Thompson WW, Shay DK, Weintraub E, et al., Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA*. 2003;289:179–86.
- Henrickson KJ and CB Hall, Diagnostic assays for respiratory syncytial virus disease. *Pediatr Infect Dis J* 2007;26 (Suppl):S36–40.
- Popow-Kraupp, T and JH Aberle, Diagnosis of Respiratory Syncytial Virus Infection. *Open Microbiol J* 2011;5:128–34. Published online 2011 December 30.
- Barenfanger J, Drake C, Leon N, et al., Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. *J Clin Microbiol* 2000;38:2824–8.

10. Clinical and Laboratory Standards Institute. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI 2005. Wayne, PA.
11. Garner JS, Hospital Infection Control Practices Advisory Committee, U. S. Department of Health and Human Services, Center for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect Control Hospit Epidemiol 1996;17:53-80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

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Français

## For Rapid Detection of Respiratory Syncytial Virus (RSV)

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##### Best.- Nr. Beschreibung

256038	BD Veritor System for Rapid Detection of Respiratory Syncytial Virus (RSV), CLIA-befreites Kit, 30 Tests
256042	BD Veritor System for Rapid Detection of Respiratory Syncytial Virus (RSV), Labor-Testkit, 30 Tests
256045	BD Veritor System for Rapid Detection of Flu A+B, CLIA-befreites Kit, 30 Tests
256041	BD Veritor System for Rapid Detection of Flu A+B, Labor-Testkit, 30 Tests
256061	BD Veritor System RSV Control Swab Set, 10 Tupferpaare
220252	COPAN Flexible Minitip Flocked Swab, 100 Tupfer
256066	BD Veritor Plus Analyzer
256068	BD Veritor InfoScan-Modul
443907	USB-Drucker kabel für BD Veritor Analyzer

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## Change History

Revision	Date	Change Summary
(10)	2018-12	<p>Changed name of BD Synapsys™ Microbiology Informatics Solution to BD Synapsys™ Informatics Solution. Replaced figure in "Specimen Collection and Handling" section to illustrate proper nasopharyngeal sampling technique.</p> <p>Made grammatical and content changes to provide more clarity in usage instructions.</p> <p>Clarified brand name references of BD Veritor Plus Analyzer.</p> <p>Replaced items in Availability list.</p>
(11)	2020-04	<p>Additional result documentation capabilities made possible with the implementation of the BD Synapsys Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect (replacing the BD Veritor System Reader and the BD InfoSync Module).</p> <p>In sections C and D, revised Steps 10C and 10D regarding the use of the BD Veritor Plus Analyzer with the BD Veritor InfoScan.</p> <p>Added access information to obtain the document from <a href="http://bd.com/e-labeling">bd.com/e-labeling</a>.</p>

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US Customers only: For symbol glossary, refer to [bd.com/symbols-glossary](http://bd.com/symbols-glossary).

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Fabricant / Proizvođač / Gyártó / Fabricante / Аткарушы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производител / Výrobca / Proizvođač / Tillverkare / Üretici / Виробник / 生产厂商

Use by / Используйте до / Spotføjebute do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption / 사용 기한 / Upotrijebiti do / Felhasználhatóság dátuma / Usare entro / Дейін пайдаланура / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použítí do / Upotrebiti do / Använd före / Son kullanna tarihi / Використати до/line / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)  
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)  
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
AAAA-MM-DD / AAAA-MM (MM = slutning af måned)  
JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)  
EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)  
AAAA-MM-DD / AAAA-MM (MM = fin del mes)  
AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)  
AAAA-MM-JJ / AAAA-MM (MM = fin du mois)  
GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)  
ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja)  
AAAA-MM-GG / AAAA-MM (MM = fine mese)  
ЖӨЖӨЖ-АА-КК / ЖӨЖӨЖ-АА / (АА = айдың соңы)  
YYYY-MM-DD/YYYY-MM (MM = 월말)  
MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)  
GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)  
JJJJ-MM-DD / JJJJ-MM (MM = einde maand)  
AAAA-MM-DD / AAAA-MM (MM = slutten av måneden)  
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
AAAA-MM-DD / AAAA-MM (MM = fim do mês)  
AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)  
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца)  
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)  
GGGG-MM-DD / GGGG-MM (MM = kraj meseca)  
AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)  
YYYY-AA-GG / YYYY-AA (AA = айн сопу)  
PPPP-MM-DD / PPPP-MM (MM = кінець місяця)  
YYYY-MM-DD / YYYY-MM (MM = 月末)

REF Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Katalooginumber / Numéro catalogue / Kataloški broj / Katalógusszám / Numero di catalogo / Каталог номери / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs / Catalogus number / Numer katalogowy / Număr de catalog / Номер по каталору / Katalogové číslo / Kataloški broj / Katalog numarasi / 目錄号

EC REP Authorized Representative in the European Community / Оторизирани представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriseret repræsentant i De Europæiske Fællesskaber / Autoriseret Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogus / Représentant autorisé pour la Communauté européenne / Autorizirani predstavnik u Europskoj uniji / Meghatalmazott képviselő az Európai Közösségekben / Rappresentante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / 유럽 공동체의 위임 대표 / Igaltais atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevogde vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Reprezentant autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom spoločenstve / Autorizovano predstavníctvo u Europskoj uniji / Auktoriserad representant i Europeiska gemenskapen / Avrupa Topluluğu Yetkili Temsilcisi / Уповноважений представник у країнах ЄС / 欧洲共同体授权代表

**IVD**

In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostika meditsiinaparatuur / Dispositif médical de diagnostic in vitro / Medicinska pomagaia za In Vitro Dijagnostiku / In vitro diagnostikai orvosi eszköz / Dispositivo medicale per diagnostica in vitro / Жасанды жағдайда жұргізілетін медициналық диагностика аспабы / In Vitro Diagnostik 의료 기기 / In vitro diagnostikos prietaisais / Medicinas ierīces, ko lieto in vitro diagnostikā / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medicinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispositiv medical pentru diagnostic in vitro / Медицинский прибор для диагностики ин витро / Medicinska pomôcka na diagnostiku in vitro / Medicinali uređaj za in vitro dijagnostiku / Mediceinteknis produkt för in vitro-diagnostik / In Vitro Diagnostik Tibbi Cihaz / Медицинский прибор для диагностики ин витро / 体外诊断医疗设备

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi piirang / Limites de température / Dozvoljena temperatura / Hömersékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperaturas ierobežojumi / Temperatuurilimiet / Temperaturbegrænsning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenje teploty / Ograničenje temperature / Temperaturgräns / Sicaklık sınırlaması / Обмеження температури / 温度限制

**LOT**

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Κωδικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kood / Numéro de lot / Lot (kod) / Tétel száma (lot) / Codice batch (lotto) / Топтама коды / 배치 코드 (코트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot nummer / Batch-code (parti) / Kod partii (serie) / Código do lote / Cod de serie (lot) / Код партии (LOT) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Partii Kodu (lot) / Код партії / 批号 (亚批)

Contains sufficient for <n> tests / Содержит достаточно в достаточном для <n> теста / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Contenido suficiente para <n> pruebas / Kullaldane <n> testide jaoks / Contenu suffisant pour <n> tests / Sadržaj za <n> testova / <n> testzhez elegendő / Contenido suficiente para <n> teste / <n> тесттери үшін жеткілікті / <n> 테스트가 충분히 포함됨 / Pakankamas kiekis atlikti <n> testu / Satur pietiekami <n> pārbaudēm / Inhoud voldoende voor "n" testen / Inholder tilstrækkelig til <n> tester / Zawiera ilość wystarczającą do <n> testów / Conteúdo suficiente para <n> testes / Conținut suficient pentru <n> teste / Достаточно для <n> тестов(a) / Obsah vystačí na <n> testov / Sadržaj dovoljan za <n> testova / Innehåller tillräckligt för <n> analyser / <n> test için yeterli mizmeze içerir / Вистачить для аналізів: <n> / 足够进行 <n> 次检测

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prostudujte pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Пайдалану нұсқаулығымен танысып алыңыз / 사용 지침 참조 / Skaitykite naudojimo instrukcijas / Skaitīt lietošanas pamācību / Raadpleeg de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se brugsanvisningen / Kullanim Talimatları'na başvurun / Див. інструкцію за використання / 请参阅使用说明

Do not reuse / Не използвайте отново / Neponoživajte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Мрън етанαχρησιμοποίητέ / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 사용하지 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbrug / Nie stosować powtórnie / Não reutilize / Nu refolositi / Не использовать повторно / Neponoživajte opakovane / Ne upotrebljavajte ponovo / Får ej återanvändas / Tekrar kullannamain / Не використовувати повторно / 请勿重复使用

**SN**

Serial number / Серийни номер / Sériové číslo / Seriennummer / Seriennummer / Σειριακός αριθμός / Nº de serie / Seerianumber / Numéro de série / Serijski broj / Rozsazatszám / Numero di serie / Топтамалық нөмірі / 일련 번호 / Serijos numeris / Serijas numurs / Série nummer / Numer serijny / Número de série / Număr de serie / Серийный номер / Seri numaras / Номер серії / 序列号

For IVD Performance evaluation only / Само за оценка качеството на работа на IVD / Pouze pro vyhodnocení výkonu IVD / Kun til evaluering af IVD ydelse / Nur für IVD-Leistungsbewertungszwecke / Móno για αξιολόγηση απόδοσης IVD / Sólo para la evaluación del rendimiento en diagnóstico in vitro / Ainult IVD seadme hindamiseks / Réserve à l'évaluation des performances IVD / Samo u znanstvene svrhe za In Vitro Dijagnostiku / Kizárolag in vitro diagnosztikához / Solo per valutazione delle prestazioni IVD / Жасанды жағдайда «пробирка ішінде» диагностикада тек жұмысты бағалау үшін / IVD 성능 평가에 대해서만 사용 / Tik IVD prietais veikimo charakteristikoms tikrinti / Vienigi IVD darbības novērtēšanai / Uitsluitend voor doeltreffendheidsonderzoek / Kun for evaluering av IVD-ytelse / Tyklo ko oceny wydajności IVD / Uso exclusivo para evaluación de IVD / Numai pentru evaluarea performanței IVD / Только для оценки качества диагностики ин витро / Určené iba na diagnostiku in vitro / Samo za procenu učinka u in vitro dijagnostici / Endast för utvärdering av diagnostisk användning in vitro / Yalnızca IVD Performans değerlendirmesi için / Тільки для оцінювання якості діагностики ин витро / 仅限 IVD 性能评估

For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Κατώτερο όριο θερμοκρασίας / Limite inferior de temperatura / Alumine temperatuuripiir / Limite inférieure de température / Najniža dozvoljena temperatura / Alsó hőmérsékleti határ / Limite inferiore di temperatura / Температураның төменгі рұқсат шегі / 하한 온도 / Žemiausia laikymo temperatūra / Temperaturas zemākā robeža / Laagste temperatuurilimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite minimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodnja hranica teploty / Donja granica temperature / Nedre temperaturgräns / Sicaklık alt sınırı / Мінімальна температура / 温度下限

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This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

R<sub>x</sub> Only

bd.com/e-labeling

Europe, CH, GB, NO: +800 135 79 135	
International: +31 20 794 7071	
AR +800 135 79 135	LT 8800 30728
AU +800 135 79 135	MT +31 20 796 5693
BR 0800 591 1055	NZ +800 135 79 135
CA +1 855 805 8539	RO 0800 895 084
CO +800 135 79 135	RU +800 135 79 135
EE 0800 0100567	SG 800 101 3366
GR 00800 161 22015 7799	SK 0800 606 287
HR 0800 804 804	TR 00800 142 064 866
IL +800 135 79 135	US +1 855 236 0910
IS 800 8996	UY +800 135 79 135
LI +31 20 796 5692	VN 122 80297

Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, Maryland 21152 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

BD, the BD logo, Synapsys, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2020 BD. All rights reserved.

# Insert ends here

DO NOT PRINT THIS PAGE

**CLSI document should not be included in print file to vendor.**

**NOTES:**

The CLSI document is to be updated at the same time as the package insert and should be released as a separate PDF document to the web.

The package insert should also be released as a separate PDF document to the web.

The package insert and CLSI document should be released as one document to SAP.

Final art provided to Buyer/print vendor should not include this page or the CLSI portion of this document.

Proofreader: please proof the CLSI document along with the package insert.

**Procedure\* BD Veritor™ System For Rapid Detection of Respiratory Syncytial Virus (RSV)**

For use with nasopharyngeal swab specimens.

Prepared by	Date Adopted	Supersedes Procedure #

Review Date	Revision Date	Signature

Distributed to	# of Copies	Distributed to	# of Copies

---

**\*Any modifications to this document are the sole responsibility of the facility. This “Sample Procedure” is not intended as a substitute for your facility procedure manual, instrument manual, or reagent labeling/package insert. This “Sample Procedure” is intended as a model for use by your facility to meet the needs of your laboratory.**

# For Rapid Detection of Respiratory Syncytial Virus (RSV)

## CLIA Complexity-WAIVED

For use with nasopharyngeal swab specimens.

For *in vitro* diagnostic use only.

**A Certificate of Waiver is required to perform this test in a CLIA waived setting. To obtain a Certificate of Waiver, please contact your state health department.**

**Additional CLIA waiver information is available at the Centers for Medicare and Medicaid website at [www.cms.hhs.gov/CLIA](http://www.cms.hhs.gov/CLIA) or from your state health department.**

**Failure to follow the instructions or modification to the test system instructions will result in the test no longer meeting the requirements for waived category.**

## INTENDED USE

The BD Veritor™ System for Rapid Detection of Respiratory Syncytial Virus (RSV) is a chromatographic immunoassay with an instrumented read for the direct and qualitative detection of RSV fusion protein from a direct nasopharyngeal swab from patients suspected of having a viral respiratory infection. This test is intended for *in vitro* diagnostic use to aid in the diagnosis of RSV infections in infants and pediatric patients under the age of 6 years. Negative results do not preclude RSV infection and should not be used as the sole basis for treatment or for other management decisions. A negative test is presumptive. It is recommended that negative test results be confirmed by viral cell culture or an alternative method, such as a FDA-cleared molecular assay. The test is intended for professional and laboratory use. It is to be used in conjunction with the BD Veritor System Instrument.

## SUMMARY AND EXPLANATION

Viral respiratory tract infections are responsible for widespread disease. Respiratory syncytial virus is a leading cause of lower respiratory tract infections (LRI) in young children in both the developed and developing worlds. Worldwide, it is estimated that RSV is responsible for greater than 30 million cases of LRI in children under 5 years of age each year.<sup>1,2</sup>

Diagnostic methods for detection of respiratory viruses include viral cell culture, direct fluorescent antibody (DFA), rapid immunoassays, and nucleic acid amplification assays such as the polymerase chain reaction (PCR).<sup>3,4</sup> Each has been demonstrated to have clinical utility for the detection of respiratory viruses including RSV. Rapid immunoassays available for specific viruses such as influenza A/B and RSV allow a quick diagnosis so that patients may be appropriately isolated and treated to prevent the nosocomial spread of infections to fellow patients with compromised cardiac, respiratory or immune functions.<sup>5</sup> In addition, rapid tests assist with the selection of appropriate antiviral therapy. The most common specimen types collected for RSV testing include nasopharyngeal washes, nasopharyngeal aspirates, nasal swabs and nasopharyngeal swabs.

The BD Veritor System for Rapid Detection of RSV (also referred to as the BD Veritor System and BD Veritor System RSV) is a chromatographic immunoassay to detect RSV fusion protein extracted from nasopharyngeal swab specimens from symptomatic patients. All BD Veritor System RSV test devices are interpreted by a BD Veritor System Instrument, either a BD Veritor Reader or BD Veritor Plus Analyzer (the "Analyzer"). When using the BD Veritor Plus Analyzer, workflow steps depend on the selected operational mode and the Analyzer configuration settings. In Analyze Now mode, the instrument evaluates assay devices after manual timing of their development. In Walk Away mode, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Connection of the Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the implementation of the BD Synapsys™ Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect. Please refer to the Analyzer *Instructions for Use* for details on these features and contact BD Technical Support for more information.

## PRINCIPLES OF THE PROCEDURE

The BD Veritor System for Rapid Detection of RSV is a qualitative, digital immunoassay for the detection of RSV fusion protein in samples processed from nasopharyngeal specimens. When specimens are processed and added to the test device, RSV antigen binds to anti-RSV antibodies conjugated to detector particles in the RSV test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of RSV antibody on the membrane. A positive result for RSV is determined by the BD Veritor System Instrument (purchased separately) when antigen-conjugate is deposited at the Test "T" position and the Control "C" position on the BD Veritor System RSV assay device. The Instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

## REAGENTS

The following components are included in the BD Veritor System for Rapid Detection of RSV kit:

BD Veritor System RSV Devices	30 devices	Foil pouched device containing one reactive strip. Each strip has one test line of monoclonal antibody specific to RSV viral antigen and murine monoclonal control line antibodies.
<b>RV Reagent D</b>	30 tubes with 400 µL reagent	Detergent with <0.1% sodium azide (preservative).
Flexible minitip flocked swab	30 each	Swab for nasopharyngeal collection
RSV Positive Control Swab	1 each	RSV Positive Control Swab, RSV antigen (noninfectious cell lysate) with <0.1% sodium azide (preservative)
RSV Negative Control Swab	1 each	RSV Negative Control Swab, (detergent-treated non-infected cells) with <0.1% sodium azide (preservative)

**Materials Required But Not Provided:** BD Veritor™ Plus Analyzer (Cat. No. 256066), Timer, Tube Rack for specimen testing.

**Optional Equipment:** BD Veritor™ InfoScan Module (Cat. No. 256068), USB Printer cable for BD Veritor™ Analyzer (Cat. No. 443907), Epson Printer model TM-T20 II. BD Veritor Plus Connect (contact BD Technical Services for details).

### Warnings and Precautions:

#### Warning

**H302** Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

**P273** Avoid release to the environment. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. The RSV Positive Control Swab and the positive control line on the BD Veritor System for Rapid Detection of RSV device have been prepared from RSV-infected tissue culture cells which have been inactivated by detergent treatment and sonication then subsequently tested by bioassay procedures.
4. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>6-9</sup> and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used BD Veritor System test devices as biohazardous waste in accordance with federal, state and local requirements.
6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. Other than the flocked swabs that are used for specimen collection, kit components should not come in contact with the patient.
8. Do not use kit components beyond the expiration date.
9. Do not reuse the BD Veritor System test device.
10. Do not use the kit if the Control RSV Positive Swab and Control RSV Negative Swab do not yield appropriate results.
11. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
12. To avoid erroneous results, specimens must be processed as indicated in the assay procedure section. The addition of excess sample may give invalid test results.
13. Proper specimen collection, storage and transport are critical to the performance of this test.
14. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.

**Storage and Handling:** Kits may be stored at 2–30 °C. **DO NOT FREEZE.** Reagents and devices must be at room temperature (15–30 °C) when used for testing.

## SPECIMEN COLLECTION AND HANDLING

**Specimen Collection and Preparation:** The acceptable specimen for testing with the BD Veritor System for Rapid Detection of RSV is a nasopharyngeal (NP) swab. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers.

Inadequate specimen collection or improper specimen handling and/or transport may yield a false negative result; therefore, training in specimen collection is highly recommended due to the importance of specimen quality for generating accurate test results.

### Specimen Transport and Storage:

Freshly collected specimens should be processed and tested within one hour. It is essential that correct specimen collection and preparation methods be followed.


1. The BD Veritor System RSV Kit includes swabs with a flocked tip for nasopharyngeal specimen collection.

2. Insert the swab into one nostril of the patient, reaching the surface of the posterior nasopharynx.

3. Rotate the swab over the surface of the posterior nasopharynx.

4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the BD Veritor System Kit.

### Dos and Don'ts of Sample Collection

- Do collect sample as soon as possible after onset of symptoms
- Do test sample immediately
- BD recommends flocked swabs which are provided in the BD Veritor System RSV Kit
- Do not use cotton tips and wood shafts
- Do not use calcium alginate swabs

## TEST PROCEDURE FOR NASOPHARYNGEAL SWABS

### NOTES:


- **Reagents, specimens and devices must be at room temperature (15–30 °C) for testing.**
- The CLIA-waived BD Veritor System for Rapid detection of RSV kit is only intended for nasopharyngeal swab specimens that are collected and tested directly (i.e dry swabs that have **NOT** been placed in transport media). The kit includes a pre-diluted processing reagent in a ready to use "unitized" tube. This CLIA-waived kit is **NOT INTENDED** for testing liquid samples such as wash or aspirate samples or swabs in transport media as results can be compromised by over dilution.

### Prepare for testing

The following steps assume that users of a BD Veritor Plus Analyzer have chosen and set all configuration options, and that the Analyzer is ready to use. To choose or change these settings, see the BD Veritor Plus Analyzer *Instructions for Use*, section 4.7. A printer is not necessary to display results. However, if your facility has chosen to connect the BD Veritor Plus Analyzer to a printer, check that the Analyzer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

#### Step 1: For each patient specimen:


- Remove one RV Reagent D tube/tip and one BD Veritor System RSV device from its foil pouch immediately before testing.
- Label with patient's name or ID number.
- Place the labeled RV Reagent D tube(s) in the designated area of the tube rack.

### Prepare the Sample


#### Step 2:

- Remove and discard the cap from the RV Reagent D tube corresponding to the sample to be tested.


#### Step 3:

- Insert the patient sample swab all the way into the RV Reagent D tube and swirl it against the inside wall three (3) times.

#### Step 4:

- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Properly dispose of the swab.

#### Step 5:

- Press the attached tip firmly onto the RV Reagent D tube containing the processed sample (threading/twisting not required).
- Vortex or mix thoroughly by swirling or flicking the bottom of the tube.
- Do not use tips from any other product, including other products from BD or other manufacturers.


After step 5, choose from the model and workflow option below before continuing to step 6:

	BD Veritor Reader or Analyzer in <b>Analyze Now</b> mode	BD Veritor Plus Analyzer in <b>Walk Away</b> mode	BD Veritor Plus Analyzer with the BD Veritor InfoScan module in <b>Analyze Now</b> mode—or— <b>Walk Away</b> mode	
Instructions in section:	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>

**Step 6A: Adding the specimen**

- Invert the RV Reagent D tube and hold the tube vertically (approximately one inch above the labeled BD Veritor System RSV device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System RSV device.


**NOTE:** Squeezing the tube too close to the tip may cause leakage.

**Step 7A: Timing development**

- After adding the sample, allow the test to run for 10 minutes before inserting into the BD Veritor Instrument.
- **CAUTION: incorrect results may occur if development time is less than 10 minutes.** Some lines may appear on the device sooner. Do not read device visually.
- **NOTE:** If running test under laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8A: Using the BD Veritor Instrument:**

- During incubation time, turn the BD Veritor Instrument on by pressing the power button once.
- Insert assay device when the 10 minute assay development time is complete. Follow the on-screen prompt to complete the procedure.
- The status of the assay analysis process appears in the display window, followed by the result display.

**Step 9A: Record the Result**


- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately.

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

**To use Walk Away mode - Connect the AC power adapter to the Analyzer and a power source**


### Step 6B: Starting Walk Away mode

- Turn the Analyzer on by pressing the blue power button once.
- When the display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE,”
  - **Double-click** the blue power button.

### Step 7B: Adding the specimen:

- When the display window reads “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”:
  - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
  - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System RSV device.


**NOTE:** Squeezing the tube close to the tip may cause leakage.

**CAUTION:** A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.

### Step 8B: Start the development and reading sequence

- Insert the test device into the slot on the right side of the Analyzer.
  - The test device must remain horizontal to prevent spilling the specimen out of the sample well.**
- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

**Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**


### Step 9B: Record the Result

- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately.

**ATTENTION:** TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).

### Step 6C: Adding the specimen

- Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System RSV device.  
**NOTE: Squeezing the tube close to the tip may cause leakage.**

### Step 7C: Timing development

- Allow the test to develop for **10 minutes**.
- **CAUTION: incorrect results may occur if development time is less than 10 minutes.** Some lines may appear on the device sooner. Do not read device visually.
- If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.


### Step 8C: Using the Analyzer

During the incubation time, turn the BD Veritor Plus Analyzer on by pressing the blue button once.


The display window briefly shows “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When assay development time is complete and the Analyzer display window reads: “INSERT TEST DEVICE OR DOUBLE CLICK FOR WALK AWAY MODE”:
  - Insert the BD Veritor System RSV device into the BD Veritor Plus Analyzer.

### Step 9C: Using the Bar Code scanner

- Follow the prompts on the display window to complete any required barcode scans of:
  - OPERATOR ID
  - SPECIMEN ID and/or
  - KIT LOT NUMBER.

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new test sequence.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the Kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

After required scans are completed, the Analyzer displays a countdown timer and test analysis begins.

- **Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**
- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

If a printer is not connected, record the result before removing the assay device.

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

### Step 10C: Remove the test device


- Remove and then discard the test device appropriately. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test.

If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

**D****Using the BD Veritor Plus Analyzer In “Walk Away” mode:  
with the BD Veritor InfoScan module installed****To use Walk Away mode - Connect the AC power adapter to the Analyzer and a power source****Step 6D: Starting Walk Away mode**


- Turn the Analyzer on by pressing the blue power button once.  
The display window will briefly show “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When the display window reads:  
“INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”  
– **Double-click the blue power button.**

**Step 7D: Using the Bar Code scanner**


- Follow the prompts on the display window to complete any required barcode scans of:
  - OPERATOR ID
  - SPECIMEN ID and/or
  - KIT LOT NUMBER


- **Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 6D. To restart this step, double-click the power button.**
- **Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.**
- **The Analyzer can record the Kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.**

**Step 8D: Add the specimen to the test device:**

- When the display window reads: “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”:
  - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System RSV device sample well).
  - Gently squeeze the ridged portion of the tube away from the tip, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System RSV device. **NOTE: Squeezing the tube close to the tip may cause leakage.**
  - **CAUTION: A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.**

**Step 9D: Start the development and reading sequence**

- Insert the test device into the slot on the right side of the Analyzer.
- The test device must remain horizontal to prevent spilling the specimen out of the sample well.**
- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
  - The display window shows the remaining analysis time.

**Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**

- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).**

**Step 10C: Remove the test device**

- Remove and then discard the test device appropriately. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to Analyze Now mode at the conclusion of each read sequence.

If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

**OPTIONAL TEST PROCEDURE: Testing for RSV and INFLUENZA A+B using a single NP swab. For patients under the age of 6 years.**

**Note: The BD Veritor™ System for Rapid Detection of Flu A+B (Cat. No. 256045) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of RSV (Cat. No. 256038).**

This procedure allows for use of the remaining processed sample from Step 5 above to test additionally for Influenza A+B. When using this optional test procedure, the sample may be used up to 15 minutes after initial processing.

1. Collect NP swab from the patient and follow Steps 1–5 of the test procedure above as instructed for RSV.
2. Using the sample from Step 5, continue the test procedure using the test device for Flu A+B.
3. Refer to the product insert for BD Veritor System for Rapid Detection of Flu A+B, (Cat. No. 256045) for the test procedure and full description of the BD Veritor Flu A+B test. Follow the Instrument on-screen prompts to complete the procedure and obtain test results. Refer to the product insert for the BD Veritor™ System Flu A+B CLIA waived Kit (Cat. No. 256045) for result interpretation.

#### **INTERPRETATION OF RESULTS:**

The BD Veritor System Instrument (purchased separately) must be used for interpretation of all test results. Operators should not attempt to interpret assay results directly from the test strip contained within the BD Veritor System RSV assay device.

<b>Display</b>	<b>Interpretation</b>
RSV: +	Positive Test for RSV (RSV antigen present)
RSV: -	Negative Test for RSV (no RSV antigen detected)
CONTROL INVALID	Control line error. Repeat the test.

**Invalid Test** – If the test is invalid, the BD Veritor System Instrument will display “CONTROL INVALID” and the test or control must then be repeated. If the “CONTROL INVALID” reading recurs, contact BD Technical Support.

#### **REPORTING OF RESULTS**

**Positive Test** Positive for the presence of RSV antigen. A positive result may occur in the absence of viable virus.

**Negative Test** Negative for the presence of RSV antigen. Infection due to RSV cannot be ruled out because the antigen present in the sample may be below the detection limit of the test. In the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral cell culture or an FDA-cleared RSV molecular assay.

**Control Invalid** Do not report results. Repeat the test.

#### **QUALITY CONTROL:**

**To utilize the Analyzer’s QC documentation capability, specimen barcode scanning must be enabled on an Analyzer equipped with the BD Veritor InfoScan module. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration.**

Each BD Veritor RSV device contains both positive and negative internal/procedural controls:

1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
2. The membrane area surrounding test lines functions as a background check on the assay device.

**The BD Veritor System Instrument evaluates the positive and negative internal/procedural controls after insertion of each BD Veritor System test device. The BD Veritor System Instrument prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result.**

**NOTE: The internal controls do not assess proper sample collection technique**

#### **External Positive and Negative Controls:**

RSV + and RSV - control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the BD Veritor System Instrument perform as expected. Prepare kit control swabs and test using the same procedure (either Analyze Now or Walk Away mode) as used for patient specimen swabs. When using the barcode scanning feature to document QC procedures, scan the barcode on the control swab packaging when prompted for a Specimen ID.

Your laboratory’s standard Quality Control procedures and applicable local, state and/or federal regulations or accreditation requirements dictate the performance of external quality control procedures.

BD recommends controls be run once for:

- each new kit lot,
- each new operator,
- each new shipment of test kits,
- as required by internal quality control procedures and in accordance with local, state and federal regulations or accreditation requirements.

#### **Test Procedure For Kit Swab Controls:**

1. Remove and discard the cap from the RV Reagent D tube corresponding to the sample to be tested.
2. Insert the control swab into the tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
4. Continue processing the swab according to the Test Procedure for Nasopharyngeal Swabs above, beginning at Step 5.

If the kit controls do not perform as expected, do not test patient specimens. Contact BD Technical Services at 1.800.638.8663.

#### **LIMITATIONS OF THE PROCEDURE**

- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of RSV antigens from nasopharyngeal swabs.
- The BD Veritor System for Rapid Detection of RSV is capable of detecting both viable and non-viable RSV particles. The BD Veritor System for Rapid Detection of RSV performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Results from the BD Veritor System for Rapid Detection of RSV test should be correlated with the clinical history, epidemiological data and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of RSV infection, and should be confirmed by viral cell culture or in the U.S., an FDA-cleared RSV molecular assay.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other non-RSV viral or bacterial infections.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no RSV activity when disease prevalence is low. False negative test results are more likely during peak RSV activity when prevalence of disease is high.
- This device has been evaluated for use with human specimen material only.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, RSV viruses that have undergone minor amino acid changes in the target epitope region.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection.
- The validity of the BD Veritor System for Rapid Detection of RSV test has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
- Therapeutic anti-RSV monoclonal antibodies may interfere with the BD Veritor System for Rapid Detection of RSV.
- Performance characteristics have not been established for use with patients older than 5 years of age or for immunocompromised patients.

#### **EXPECTED VALUES**

The rate of positivity observed in RSV testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, time of year, age of the patient, geographic location and most importantly, local disease prevalence. In the 2012/2013 clinical trial the overall prevalence of RSV as determined by viral cell culture positive specimens was 25.6% (range of 7.7 to 65.2%). Overall prevalence of RSV as determined by PCR positive specimens tested was 34.2% (range of 15.4 to 69.6%).

#### **PERFORMANCE CHARACTERISTICS**

##### **Clinical Performance**

Performance characteristics for the BD Veritor System for Rapid Detection of RSV test were established in a prospective multi-center clinical study conducted at eight Point of Care (POC) U.S. testing sites during the 2012–2013 respiratory season. The performance of the BD Veritor RSV test was compared to a commercially available PCR method as well as viral culture. A total of 540 specimens were enrolled in the clinical trial. A total of 523 specimens had evaluable results for all three testing methods: PCR, viral cell culture and BD Veritor RSV. The study population was 42.6% females and 57.4% males. The table below shows age distribution of the study population.

Demographics Summary – Age Group		
Age Group	Number	Percentage
<2	305	58.3
2–5	218	41.7
<i>Total</i>	<i>523</i>	<i>100</i>

The table below summarizes the performance obtained with clinical specimens using the BD Veritor System RSV test in comparison to a commercially available PCR. The overall positive percent agreement (PPA) and negative percent agreement (NPA) of the BD Veritor System RSV with a PCR comparator, based on these 523 specimens, are 81.6% (146/179) and 99.1% (341/344), respectively.

BD Veritor RSV Compared to PCR			
BD Veritor RSV	PCR		
	P	N	Total
<b>P</b>	146	3	149
<b>N</b>	33	341	374
<i>Total</i>	<i>179</i>	<i>344</i>	<i>523</i>

Reference Method: PCR  
PPA: 81.6% (95% CI: 75.2%, 86.6%)  
NPA: 99.1% (95% CI: 97.5%, 99.7%)

BD Veritor System RSV test performance compared to viral cell culture was also evaluated in this study. For the same 523 specimens 91.8% (123/134) were positive by both BD Veritor RSV and cultures, 93.3% (363/389) were negative by both BD Veritor RSV and culture. There were 26 BD Veritor RSV positive, viral cell culture negative specimens of which 23 were demonstrated to be RSV positive by an FDA cleared molecular assay.

Invalid rates for the BD Veritor System for Rapid Detection of RSV were calculated by dividing the number of invalids by the total number of evaluable specimens tested by the BD Veritor System. The overall invalid rate for the BD Veritor System for RSV based on the 523 specimens was determined to be 0.2% (1/523, 95% CI: 0.0%, 1.1%).

#### Reproducibility

The reproducibility of the BD Veritor System for Rapid Detection of RSV test was evaluated at two POC sites (P-1, P-2) and one clinical laboratory site (S-1). The reproducibility panel was composed of 12 simulated RSV samples. These included moderate positive samples, low positive samples (near the assay limit of detection), high negative samples (i.e., containing very low concentrations of virus) and negative samples. The panel was tested by two operators at each site for five consecutive days. The results are summarized below.

BD Veritor RSV Reproducibility (% RSV Positive Results)								
Sample	P-1		P-2		S-1		Total	
	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI
High negative RSV	6.7% (2/30)	1.8%, 21.3%	6.7% (2/30)	1.8%, 21.3%	13.3% (4/30)	5.3%, 29.7%	8.9% (8/90)	4.6%, 16.6%
Low positive RSV	90.0% (27/30)	74.4%, 96.5%	76.7% (23/30)	59.1%, 88.2%	80.0% (24/30)	62.7%, 90.5%	82.2% (74/90)	73.1%, 88.8%
Moderate positive RSV	100% (30/30)	88.6%, 100%	100% (30/30)	88.6%, 100%	100% (30/30)	88.6%, 100%	100% (90/90)	95.9%, 100%
Negative	0.0% (0/30)	0.0%, 11.3%	0.0% (0/30)	0.0%, 11.3%	0.0% (0/30)	0.0%, 11.3%	0.0% (0/90)	0.0%, 4.1%

## Analytical Studies

### Analytical Sensitivity (Limit of Detection)

The limit of detection (LOD) for the BD Veritor System for Rapid Detection of RSV test was established for the following RSV strains. The LOD for each strain represents the lowest concentration producing a positivity rate of  $\geq 95\%$  based on testing 60 to 80 replicates.

RSV Viral Strain	Calculated LOD (TCID <sub>50</sub> /mL)	No. Positive / Total	% Positive
VR-26 (Long subgroup A)	1.43 x 10 <sup>5</sup>	57/60	95.0
VR-955 (9320 subgroup B)	3.98 x 10 <sup>4</sup>	57/60	95.0
VR-1540 (A-2)	1.94 x 10 <sup>3</sup>	59/60	98.3
VR-1580 (Washington subgroup B)	1.08 x 10 <sup>4</sup>	58/60	96.7
VR-1400 (Wild Type subgroup B)	2.96 x 10 <sup>3</sup>	76/80	95.0

TCID<sub>50</sub>/mL = 50% Tissue Culture Infectious Dose

### Analytical Specificity (Cross Reactivity)

The BD Veritor System for Rapid Detection of RSV test was evaluated with bacteria and yeast at a target concentration of approximately 5 x 10<sup>6</sup> CFU/mL (CFU – Colony Forming Units). The viruses were evaluated at concentrations of 10<sup>4</sup>TCID<sub>50</sub>/mL or greater. Of the microorganisms tested, none showed cross-reactivity in the RSV test.

<i>Bacteroides fragilis</i>	<i>Neisseria mucosa</i>	<i>Veillonella parvula</i>
<i>Bordetella pertussis</i>	<i>Neisseria</i> sp. ( <i>Neisseria perflaus</i> )	Adenovirus, Type 1
<i>Candida albicans</i>	<i>Neisseria subflava</i>	Adenovirus, Type 7
<i>Chlamydia pneumoniae</i>	<i>Peptostreptococcus anaerobius</i>	Cytomegalovirus
<i>Corynebacterium diphtherium</i>	<i>Porphyromonas asaccharolyticus</i>	Enterovirus
<i>Escherichia coli</i>	<i>Prevotella oralis</i>	HSV Type 1
<i>Fusobacterium nucleatum</i>	<i>Propionibacterium acnes</i>	Human Coronavirus OC43
<i>Haemophilus influenzae</i>	<i>Proteus mirabilis</i>	Human metapneumovirus (hMPV-27 A2)
<i>Haemophilus parainfluenzae</i>	<i>Pseudomonas aeruginosa</i>	Human Parainfluenza
<i>Kingella kingae</i>	<i>Serratia marcescens</i>	Influenza A/California/7/2009 H1N1
<i>Klebsiella pneumoniae</i>	<i>Staphylococcus aureus</i>	Influenza A/Brisbane/10/2007 H3N2
<i>Lactobacillus</i> sp.	<i>Staphylococcus epidermidis</i>	Influenza A/Victoria/3/75 H3N2
<i>Legionella</i> sp.	<i>Streptococcus mutans</i>	Influenza B/Brisbane/60/2008
<i>Moraxella catarrhalis</i>	<i>Streptococcus pneumoniae</i>	Influenza B/Florida/4/2006
<i>Mycobacterium tuberculosis</i>	<i>Streptococcus pyogenes</i>	Influenza B/Lee/40
<i>Mycoplasma pneumoniae</i>	<i>Streptococcus</i> sp. Group C	Measles virus
<i>Neisseria gonorrhoeae</i>	<i>Streptococcus</i> sp. Group G	Mumps virus
<i>Neisseria meningitidis</i>	<i>Streptococcus salivarius</i>	Rhinovirus

## Interfering Substances

Various substances were evaluated with the BD Veritor System for Rapid Detection of RSV test. These substances included whole blood (2%) and various medications. No interference was noted with this assay for any of the substances at the concentrations tested.

Substance	Concentration
Ayr Saline Nasal Gel	10 mg/mL
4-Acetamidophenol	10 mg/mL
Acetylsalicylic acid	20 mg/mL
Albuterol	0.083 mg/mL
Amantadine Hydrochloride	500 ng/mL
Beclomethasone	500 ng/mL
Budesonide	500 ng/mL
Chlorpheniramine maleate	5 mg/mL
Dexamethasone	10 mg/mL
Dextromethorphan	10 mg/mL
Diphenhydramine HCl	5 mg/mL
Fexofenadine	500 ng/mL
FluMist™	1%
Flunisolide	500 ng/mL
Fluticasone	500 ng/mL
Guaiaacol Glyceryl Ether	20 mg/mL
Homeopathic Allergy Medicine	10 mg/mL
Ibuprofen	10 mg/mL
Loratidine	100 ng/mL

Substance	Concentration
Menthol Throat Lozenges	10 mg/mL
Mometasone	500 ng/mL
Mupirocin	500 ng/mL
Oseltamivir	500 ng/mL
Oxymetazoline	0.05 mg/mL
Phenylephrine	1 mg/mL
Pseudoephedrine HCl	20 mg/mL
Purified Mucin Protein	1 mg/mL
Ribavirin	500 ng/mL
Rimantadine	500 ng/mL
Synagis	4 µg/mL
Tobramycin	500 ng/mL
Triamcinolone	500 ng/mL
Zanamivir	1 mg/mL
Four OTC nasal sprays	10%
Four OTC throat drops	12.5%
Two OTC mouthwashes	5%
Whole Blood	2%

## CLIA WAIVER STUDY

The accuracy of the BD Veritor System for Rapid Detection of RSV test was evaluated at eight non-laboratory intended use Point of Care (POC) sites. A total of 22 operators representative of CLIA waived site personnel (intended users) participated in the study. No training on the use of the test was provided. The BD Veritor RSV test results obtained by the intended users were compared with results obtained by a commercially available PCR method. There were 523 prospectively collected specimens evaluated in this study. Among those there were 179 specimens that tested positive and 344 that tested negative by the PCR method. The positive percent agreement (PPA) of the BD Veritor RSV assay for specimens that tested positive by the comparator PCR method was 81.6% (146/179) with the 95% confidence interval (CI) of 75.2–86.6%, and the negative percent agreement (NPA) was 99.1% (341/344) with the 95% CI of 97.5–99.7%.

The results of the BD Veritor RSV test obtained on the same specimens by the intended users were also compared with the results obtained by viral cell culture. Among the 523 specimens tested, there were 134 specimens that tested positive and 389 that tested negative by the viral culture method. The sensitivity of the BD Veritor RSV assay for specimens that tested positive by the comparator viral cell culture was 91.8% (123/134) with the 95% CI of 85.9–95.4%, and the specificity was 93.3% (363/389) with the 95% CI of 90.4–94.4%. Of the 26 BD Veritor RSV positive, viral cell culture negative specimens there were 23 specimens that tested positive by the PCR comparator method.

The summary of the results is presented below.

BD Veritor RSV Compared to PCR				BD Veritor RSV Compared to Viral Culture		
BD Veritor RSV	PCR			Viral Culture		
	P	N	Total	P	N	Total
<b>P</b>	146	3	149	123	26*	149
<b>N</b>	33	341	374	11	363	374
<b>Total</b>	179	344	523	134	389	523
PPA: 81.6% (95% CI: 75.2%, 86.6%) NPA: 99.1% (95% CI: 97.5%, 99.7%)				Sensitivity: 91.8% (95% CI: 85.9%, 95.4%) Specificity: 93.3% (95% CI: 90.4%, 94.4%)		
*Of the 26 BD Veritor RSV positive, viral cell culture negative specimens, there were 23 specimens that tested positive by the PCR comparator method.						

The overall invalid rate for the BD Veritor System for RSV based on the 523 specimens tested in the study was 0.2% (1/523) with the 95% CI of 0.0–1.1%.

Another study was designed to assess the capability of untrained users to test weakly reactive samples and deliver results with accuracy. This study was conducted at three CLIA waived intended use sites using simulated swab samples. The samples were spiked with RSV virus targeting three concentrations (high negative ~5% positivity, low positive ~95% positivity and moderate positive ~100% positivity). The swabs were provided to the operators in panels that were masked and randomized prior to shipment to the sites. The samples were also tested by trained laboratorians at one clinical site. Each site had two operators who each tested the panel on each of ten days. Each swab was processed and tested in a single device according to the test procedure.

The table below shows the rate of RSV detection for weakly reactive samples when the test was used by untrained intended users. The rate of RSV detection for weakly reactive samples when the test was used by trained laboratorians at one clinical laboratory site is also provided.

Sample	Untrained Users						Trained Users	
	POC 1		POC 2		POC 3		Clinical Site	
	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI	% Pos	95% CI
High negative RSV	5.0% (1/20)	0.9%, 23.6%	0.0% (0/20)	0.0%, 16.1%	5.0% (1/20)	0.9%, 23.6%	15.0% (3/20)	5.2%, 36.0%
Low positive RSV	75.0% (15/20)	53.1%, 88.8%	75.0% (15/20)	53.1%, 88.8%	70.0% (14/20)	48.1%, 85.5%	75.0% (15/20)	53.1%, 88.8%
Moderate positive RSV	100% (20/20)	83.9%, 100%	100% (20/20)	83.9%, 100%	100% (20/20)	83.9%, 100%	100% (20/20)	83.9%, 100%

Using risk analysis as a guide, analytical flex studies were conducted. The studies demonstrated that the test is insensitive to stresses due to environmental conditions or from potential user errors.

### Technical Support

For questions, or to report a problem please call **Technical Support at 1.800.638.8663**. Test system problems may also be reported to the FDA using the MedWatch reporting system (phone: 1.800.FDA.1088; fax: 1.800.FDA.1078; or <http://www.fda.gov/medwatch>).

### AVAILABILITY

#### Cat. No. Description

256038	BD Veritor™ System for Rapid Detection of Respiratory Syncytial Virus (RSV), CLIA waived kit, 30 tests
256042	BD Veritor™ System for Rapid Detection of Respiratory Syncytial Virus (RSV) Laboratory kit, 30 tests
256045	BD Veritor™ System for Rapid Detection of Flu A+B, CLIA Waived kit, 30 tests
256041	BD Veritor™ System for Rapid Detection of Flu A+B, Laboratory kit, 30 tests
256061	BD Veritor™ System RSV Control Swab Set, 10 pairs of swabs
220252	COPAN Flexible Minitip Flocked Swab, 100 swabs
256066	BD Veritor™ Plus Analyzer
256068	BD Veritor™ InfoScan Module
443907	USB Printer Cable for BD Veritor™ Analyzer

To network a BD Veritor Plus Analyzer to an LIS, contact BD Technical Services for details.

## REFERENCES

1. Hall CB, Weinberg GA, Iwane MK, et al., The Burden of Respiratory Syncytial Virus Infection in Young Children. *N Engl J Med* 2009;360:588–98.
2. Nair H, Nokes DJ, Gessner BD, et al., Global burden of acute lower respiratory infections due to respiratory syncytial virus in young children: a systematic review and meta-analysis. *Lancet* 2010;375:1545–55.
3. Falsey AR and EE Walsh, Respiratory syncytial virus infection in adults. *Clin Microbiol Rev* 2000;13:371–84.
4. Murata Y and AR Falsey, RSV Infection in Elderly Adults, In: Patricia Cane, Editor(s), *Perspectives in Medical Virology*, 2006, Elsevier, Volume 14, Pages 163–82.
5. Crowcroft NS, Cutts F and Zambon MC, Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. *Commun Dis Public Health*. 1999;2:234–41.
6. Thompson WW, Shay DK, Weintraub E, et al., Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA*. 2003;289:179–86.
7. Henrickson KJ and CB Hall, Diagnostic assays for respiratory syncytial virus disease. *Pediatr Infect Dis J* 2007;26 (Suppl):S36–40.
8. Popow-Kraupp, T and JH Aberle, Diagnosis of Respiratory Syncytial Virus Infection. *Open Microbiol J* 2011;5:128–34. Published online 2011 December 30.
9. Barenfanger J, Drake C, Leon N, et al., Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. *J Clin Microbiol* 2000;38:2824–8.
10. Clinical and Laboratory Standards Institute. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed., CLSI 2005. Wayne, PA.
11. Garner JS, Hospital Infection Control Practices Advisory Committee, U. S. Department of Health and Human Services, Center for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect Control Hospit Epidemiol* 1996;17:53–80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021–0045.

Technical Information: In the United States, contact BD Technical Service and Support at 1.800.638.8663 or bd.com.