Veritor™ System

For Rapid Detection of Flu A+B

CLIA-waived kit configured for testing nasal and nasopharyngeal swab samples freshly collected, processed and dispensed directly onto assay test device.

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Determinations
Déterminations
Bestimmungen
Determinazioni
Determinaciones

For Rapid Detection of Flu A+B

For use with nasal and nasopharyngeal swab specimens.

CLIA Complexity-WAIVED

For *in vitro* diagnostic use only.

A Certificate of Waiver is required to perform this test in a CLIA waived setting. To obtain a Certificate of Waiver, please contact your state health department. Additional CLIA waiver information is available at the Centers for Medicare and Medicaid website at www.cms.hhs.gov/CLIA or from your state health department.

Failure to follow the instructions or modification to the test system instructions will result in the test no longer meeting the requirements for waived category.

INTENDED USE

The BD Veritor™ System for Rapid Detection of Flu A+B is a rapid chromatographic immunoassay for the direct and qualitative detection of influenza A and B viral nucleoprotein antigens from nasal and nasopharyngeal swabs of symptomatic patients. The BD Veritor System for Rapid Detection of Flu A+B (also referred to as the BD Veritor System and BD Veritor System Flu A+B) is a differentiated test, such that influenza A viral antigens can be distinguished from influenza B viral antigens from a single processed sample using a single device. The test is to be used as an aid in the diagnosis of influenza A and B viral infections. A negative test is presumptive and it is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay. Outside the U.S., a negative test is presumptive and it is recommended that these results be confirmed by viral culture or a molecular assay cleared for diagnostic use in the country of use. FDA has not cleared this device for use outside of the U.S. Negative test results do not preclude influenza viral infection and should not be used as the sole basis for treatment or other patient management decisions. The test is not intended to detect influenza C antigens.

Performance characteristics for influenza A and B were established during January through March of 2011 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2010–2011 Season, and Composition of the 2011–2012 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

If infection with a novel influenza virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Virus culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

SUMMARY AND EXPLANATION

Influenza illness classically presents with sudden onset of fever, chills, headache, myalgias, and a non-productive cough. Epidemics of influenza typically occur during winter months with estimated 114,000 hospitalizations¹ and 36,000 deaths² per year in the U.S. Influenza viruses can also cause pandemics, during which rates of illness and death from influenza-related complications can increase dramatically.

Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 hours of onset of illness. It is important to rapidly distinguish influenza A from influenza B in order to allow physicians a choice in selective antiviral intervention. Moreover, it is important to determine if influenza A or B is causing symptomatic disease in a particular institution (e.g., nursing home) or community, so that appropriate preventative intervention can be taken for susceptible individuals. It is therefore important to not only rapidly determine whether influenza is present, but also which type of influenza virus is present.³

Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction (PCR), serology, and viral culture.⁴⁻¹¹ Immunofluorescence assays entail staining of specimens immobilized on microscope slides using fluorescent-labeled antibodies for observation by fluorescence microscopy.^{6,12,13} Culture methods employ initial viral isolation in cell culture, followed by hemadsorption inhibition, immunofluorescence, or neutralization assays to confirm the presence of the influenza virus.¹³⁻¹⁵

The BD Veritor System for Rapid Detection of Flu A+B (also referred to as the BD Veritor System and BD Veritor System Flu A+B) is a chromatographic immunoassay to detect influenza A or B nucleoprotein antigens from respiratory specimens of symptomatic patients with a time to result of 10 minutes. The speed and simplified workflow of the BD Veritor System for Rapid Detection of Flu A+B makes it applicable as a "STAT" influenza A and B antigen detection test providing relevant information to assist with the diagnosis of influenza.

All BD Veritor System Flu A+B test devices are interpreted by a BD Veritor System Instrument, either a BD Veritor Reader or BD Veritor Plus Analyzer (the "Analyzer"). When using the BD Veritor Plus Analyzer, workflow steps depend on the selected operational mode and the Analyzer configuration settings. In **Analyze Now mode**, the instrument evaluates assay devices after manual timing of their development. In **Walk Away mode**, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Connection of the Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the implementation of the BD Synapsys™ Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect. Please refer to the *Analyzer Instructions for Use* for details on these features and contact BD technical support for more information.

PRINCIPLES OF THE PROCEDURE

The BD Veritor System for Rapid Detection of Flu A+B is a qualitative, digital immunoassay for the detection of influenza A and B viral antigens in samples processed from respiratory specimens. When specimens are processed and added to the test device, influenza A or B viral antigens bind to anti-influenza antibodies conjugated to detector particles in the A+B test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of antibody on the membrane. A positive result for influenza A is determined by the BD Veritor System Instrument when antigen-conjugate is deposited at the Test "A" position and the Control "C" position on the BD Veritor System Flu A+B assay device. A positive result for influenza B is determined by the BD Veritor System Instrument when antigen-conjugate is deposited at the Test "B" position and the Control "C" position in the BD Veritor System Flu A+B assay device. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

REAGENTS

The following components are included in the BD Veritor System for Rapid Detection of Flu A+B kit:

BD Veritor System Flu A+B Devices	30 devices	Foil pouched device containing one reactive strip. Each strip has two test lines of monoclonal antibody specific to either influenza A or B viral antigen and murine monoclonal control line antibodies.
RV Reagent D	30 tubes with 400 µL reagent	Detergent with <0.1% sodium azide
Flexible minitip flocked swab	30 each	Swab for nasopharyngeal or nasal collection
Control A+/B- Swab	1 each	Flu A Positive and Flu B Negative Control Swab, influenza A antigen (inactive recombinant nucleoprotein) with <0.1% sodium azide
Control B+/A- Swab	1 each	Flu A Negative and Flu B Positive Control Swab, influenza B antigen (inactive recombinant nucleoprotein) with <0.1% sodium azide

Materials Required But Not Provided: BD Veritor™ Plus Analyzer (Cat. No. 256066), Timer, Tube Rack for specimen testing.

Optional Equipment: BD Veritor™ InfoScan Module (Cat. No. 256068), USB Printer cable for BD Veritor™ Analyzer (Cat. No. 443907), Epson Printer model TM-T20 II, BD Veritor Plus Connect (contact BD Technical Services for details).

Warnings and Precautions:

Warning

H302 Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

P273 Avoid release to the environment. **P264** Wash thoroughly after handling. **P270** Do not eat, drink or smoke when using this product. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P330** Rinse mouth. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.
4. Pathogenic microorganisms, including hepatitis viruses, Human Immunodeficiency Virus and novel influenza viruses, may be present in clinical specimens. "Standard Precautions"¹⁶⁻¹⁹ and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used BD Veritor System test devices as biohazardous waste in accordance with federal, state and local requirements.

6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. For optimal results, use the flocked swabs provided with the kit for specimen collection.
8. Other than the flocked swabs that are used for specimen collection, kit components should not make contact with the patient.
9. Do not use kit components beyond the expiration date.
10. Do not reuse the device.
11. Do not use the kit if the Control A+/B- swab and Control B+/A- swab do not yield appropriate results.
12. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
13. To avoid erroneous results, swab specimens must be processed as indicated in the assay procedure section.
14. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.
15. FluMist® is made from attenuated live flu virus and although the concentration tested (1%) was non-interfering, it is possible when tested with higher concentrations that an influenza A and/or influenza B false positive may occur.

Storage and Handling: Kits may be stored at 2–30 °C. DO NOT FREEZE. **Reagents and devices must be at room temperature (15–30 °C) when used for testing.**

SPECIMEN COLLECTION

Acceptable specimens for testing with the BD Veritor System Flu A+B test include nasal swabs and nasopharyngeal (NP) swabs. Freshly collected specimens should be processed within 1 hour. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers.

Inadequate specimen collection, improper specimen handling and/or transport may yield a false negative result; therefore, specimen collection requires specific training and guidance due to the importance of specimen quality to accurate test results.

Proper Nasal Swab Sample Collection


1. The BD Veritor System Kit includes swabs with a flocked tip for nasal specimen collection.

2. Insert the swab into one nostril of the patient.

3. Rotate the swab two complete 360-degree turns; pressing firmly against the nasal mucosa to help ensure the swab obtains both cells and mucus.

4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the BD Veritor System Kit.

Proper Nasopharyngeal Swab Sample Collection


1. The BD Veritor System Kit includes swabs with a flocked tip for nasopharyngeal specimen collection.

2. Insert the swab into one nostril of the patient, reaching the surface of the posterior nasopharynx.

3. Rotate the swab over the surface of the posterior nasopharynx.

4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the BD Veritor System Kit.

DOs and DON'Ts of Sample Collection

- Do collect sample as soon as possible after onset of symptoms
- Do test sample immediately
- BD recommends flocked swabs which are provided in the BD Veritor System Flu A+B Kit
- Do not use cotton tips and wood shafts
- Do not use calcium alginate swabs

TEST PROCEDURE

Nasal and Nasopharyngeal Swab Test Procedure

NOTES:


- Reagents, specimens and devices must be at room temperature (15–30 °C) prior to testing.
- The CLIA-waived BD Veritor System for Rapid Detection of Flu A+B kit is only intended for nasal and nasopharyngeal swab specimens that are collected and tested directly (i.e., dry swabs that have **NOT** been placed in transport media). The kit includes a pre-diluted process reagent in a ready to use "unitized" tube. This CLIA-waived kit IS NOT INTENDED for testing liquid samples such as wash or aspirate samples or swabs in transport media as results can be compromised by over dilution.

Prepare for Testing

The following steps assume that users of a BD Veritor Plus Analyzer have chosen and set all configuration options, and that the Analyzer is ready to use. To choose or change these settings, see the BD Veritor Plus Analyzer *Instructions for Use*, section 4.7. A printer is not necessary to display results. However, if your facility has chosen to connect the Analyzer to a printer, check that the printer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

Step 1: For each patient specimen:

- Remove one RV Reagent D tube/tip and one BD Veritor System Flu A+B device from its foil pouch immediately before testing.
- Label with patient's name or ID number.
- Place the labeled RV Reagent D tube(s) in the designated area of the tube rack.

Prepare the Sample


Step 2:

- Remove and discard the cap from the RV Reagent D tube corresponding to the sample to be tested.


Step 3:

- Insert the patient sample swab all the way into the RV Reagent D tube and swirl it against the inside wall three (3) times.


Step 4:

- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Properly discard of the swab.


Step 5:

- Press the attached tip firmly onto the RV Reagent D tube containing the processed sample (threading/twisting not required).

- Vortex or mix thoroughly by swirling or flicking the bottom of the tube before adding to assay device.
- Note: Do not use tips from any other product, including other products from BD or other manufacturers.


After step 5, choose from the model and workflow option below before continuing to step 6:

	BD Veritor Reader or Analyzer in Analyze Now mode	BD Veritor Plus Analyzer in Walk Away mode	BD Veritor Plus Analyzer with the BD Veritor InfoScan module in Analyze Now mode—or— Walk Away mode	
Instructions in section:	 A 	B 	C 	D 

Step 6A: Adding the specimen

- Invert the RV Reagent D tube and hold the tube vertically (approximately one inch above the labeled BD Veritor System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System Flu A+B device.


NOTE: Squeezing the tube too close to the tip may cause leakage.

Step 7A: Timing development

- After adding the sample, allow the test to run for 10 minutes before inserting into the BD Veritor Instrument.
- **CAUTION: incorrect results may occur if development time is less than 10 minutes.** Some lines may appear on the device sooner. Do not read device visually.
- **NOTE:** If running test under laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8A: Using the BD Veritor Instrument**

- During incubation time, turn the BD Veritor Instrument on by pressing the power button once.
- Insert assay device when the 10-minute assay development time is complete. Follow the on-screen prompt to complete the procedure.
- The status of the assay analysis process appears in the display window, followed by the result display.

**Step 9A: Record the Result**

- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

To use Walk Away mode - connect the AC power adapter to the Analyzer and a power source


Step 6B: Starting Walk Away Mode

- Turn the Analyzer on by pressing the blue power button once
- When the display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”
 - **Double-click** the blue power button.

Step 7B: Adding the specimen

- When the display window reads “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device.

NOTE: Squeezing the tube too close to the tip may cause leakage.

CAUTION: A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.


Step 8B: Starting the development and reading sequence

- Insert the test device into the slot on the right side of the Analyzer.

The test device must remain horizontal to prevent spilling the specimen out of the sample well

- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.


Step 9B: Record the Result

- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).

Step 6C: Adding the specimen

- Invert the tube, holding it vertically (approximately one inch above the BD Veritor SystemFlu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System Flu A+B device.

NOTE: Squeezing the tube close to the tip may cause leakage.

Step 7C: Timing development

- Allow the test to develop for 10 minutes.
- **CAUTION: incorrect results may occur if development time is less than 10 minutes.** Some lines may appear on the device sooner. Do not read device visually.
- If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.


Step 8C: Using the Analyzer

During the incubation time, turn the Analyzer on by pressing the blue button once.


The display window briefly shows “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration steps. Ignore this message and postpone this process when an assay is awaiting analysis.

- When assay development time is complete and the Analyzer display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”:
 - Insert the BD Veritor System Flu A+B device into the slot on the right side of the Analyzer.

Step 9C: Using the Bar Code scanner

- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER

- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new reading sequence.
- Move barcodes slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the Kit Lot Number and expiration date in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

After required scans are completed, the Analyzer displays a countdown timer and test analysis begins.

- **Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**
- When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

If a printer is not connected, record the result before removing the assay device.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

Step 10C: Remove the test device


- Remove and then discard the test device appropriately. The display will show “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE” to indicate the Analyzer is ready to perform another test.


If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

D**Using the BD Veritor Plus Analyzer In “Walk Away” mode:
with the BD Veritor InfoScan module installed****To use Walk Away mode - connect the AC power adapter to the Analyzer and a power source****Step 6D: Starting Walk Away mode**

- Turn the Analyzer on by pressing the blue power button once.
The display window will briefly show “SCAN CONFIG BARCODE.” This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration steps. Ignore this message and postpone this process when an assay is awaiting analysis.
- When the display window reads:
“INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”
 - **Double-click** the blue power button.

**Step 7D: Using the Bar Code scanner**


- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER


- Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 6D. To restart this step, double-click the power button.
- Move barcodes slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the Kit Lot Number and expiration date in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

Step 8D: Add the specimen to the test device

- When the display window reads: ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY:
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
 - **Gently** squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device. **NOTE: Squeezing the tube close to the tip may cause leakage.**
 - **CAUTION: A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.**

**Step 9D: Starting the development and reading sequence**

- Insert the test device into the slot on the right side of the Analyzer. **The test device must remain horizontal to prevent spilling the specimen out of the sample well.**
- “DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.

- When analysis is complete, a result appears in the display window. If configured to display, the Specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).

Step 10D: Removing the test device

- Remove and then discard the test device appropriately. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to Analyze Now mode at the conclusion of each read sequence.

If the Veritor Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. In the event that a network connection is not detected while the ENVELOPE symbol is still displayed, the Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

OPTIONAL TEST PROCEDURE: Testing for INFLUENZA A+B and RSV using a single NP swab

Note: The BD Veritor™ System for Rapid Detection of RSV (Cat. No. 256038) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of Flu A+B (Cat. No. 256045).

IMPORTANT NOTE: THE SAMPLE TO BE TESTED IN THE RSV KIT MUST BE FROM A PATIENT LESS THAN 6 YEARS OF AGE AS INDICATED IN THE BD VERITOR RSV POC KIT PACKAGE INSERT. THE PROCESSED SAMPLE SHOULD BE TESTED WITHIN 15 MINUTES.

This procedure allows for use of the remaining processed sample from Step 5 above to test additionally for RSV. When using this optional test procedure, the sample may be used up to 15 minutes after initial processing.

1. Collect NP swab from the patient and follow Steps 1–5 of the test procedure above as instructed for the BD Veritor System for Rapid Detection of Flu A+B.
2. Using the sample from Step 5, continue the test procedure using the test device for RSV.
3. Refer to the product insert for BD Veritor System for Rapid Detection of RSV, (Cat. No. 256038) for the test procedure and full description of the BD Veritor RSV test. Follow the Instructions in the insert and the Instrument on-screen prompts to complete the test procedure and obtain results. Refer to the product insert for the BD Veritor System RSV Kit (Cat. No. 256038) for result interpretation.

INTERPRETATION OF RESULTS

The BD Veritor System Instrument (purchased separately) must be used for interpretation of all test results. Operators should not attempt to interpret assay results directly from the test strip contained within the BD Veritor System Flu A+B assay device. With some specimens, up to four lines may be visible on the test device. The Instrument will appropriately interpret the result.

Display	Interpretation
FLU A: + FLU B: -	Positive Test for Flu A (influenza A antigen present)
FLU A: - FLU B: +	Positive Test for Flu B (influenza B antigen present)
FLU A: - FLU B: -	Negative Test for Flu A and Flu B (no antigen detected)
RESULT INVALID	Result Invalid. Repeat the test.
POSITIVE CONTROL INVALID	Test Invalid. Repeat the test.
NEGATIVE CONTROL INVALID	Test Invalid. Repeat the test.

Invalid Test – If the test is invalid, the BD Veritor System Instrument will display “RESULT INVALID” or “POSITIVE CONTROL INVALID” or “NEGATIVE CONTROL INVALID” and the test or control must then be repeated. Because true dual positives are exceptionally rare, the BD Veritor System Instrument reports dual positive influenza A and influenza B results as “Result Invalid.” Specimens generating a “Result Invalid” should be retested. Upon retesting, if the specimen produces a “Result Invalid” the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus. If either the POSITIVE or NEGATIVE “CONTROL INVALID” reading recurs, contact BD Technical Support.

REPORTING OF RESULTS

- Positive Test** Positive for the presence of influenza A or influenza B antigen. A positive result may occur in the absence of viable virus.
- Negative Test** Negative for the presence of influenza A or influenza B antigen. Infection due to influenza cannot be ruled out because the antigen present in the sample may be below the detection limit of the test. It is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay.
- Invalid Test** Test result is inconclusive. Do not report results. Repeat the test.

QUALITY CONTROL

To utilize the Analyzer's QC documentation capability, specimen barcode scanning must be enabled on an Analyzer equipped with the BD Veritor InfoScan module. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration.

Each BD Veritor System Flu A+B device contains both positive and negative internal/procedural controls:

1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
2. The membrane area surrounding test lines functions as a background check on the assay device.

The BD Veritor System Instrument evaluates the positive and negative internal/procedural controls after insertion of each BD Veritor System test device. The BD Veritor Instrument also prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result. NOTE: The internal controls do not assess proper sample collection technique.

External Positive and Negative Controls:

Flu A positive/B negative and Flu B positive/A negative control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the BD Veritor System Instrument perform as expected. Prepare kit control swabs and test using the same procedure (either **Analyze Now** or **Walk Away** mode) as used for patient specimen swabs. When using the barcode scanning feature to document QC procedures, scan the barcode on the control swab packaging when prompted for a Specimen ID.

Your laboratory's standard Quality Control procedures and applicable local, state and/or federal regulations or accreditation requirements dictate the performance of external quality control procedures.

BD recommends controls be run once for:

- each new kit lot,
- each new operator,
- each new shipment of test kits,
- as required by site quality control procedures and in accordance with local, state and federal regulations or accreditation requirements.

Control Swab Test Procedure – Preparing the sample

1. Insert the swab all the way into the appropriately labeled RV Reagent D tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
2. Continue processing the swab according to the Test Procedure for Nasal and Nasopharyngeal swabs above beginning at Step 4 **"Remove the swab."**

If the kit controls do not perform as expected, do not report patient results. Contact BD Technical Services at 1.800.638.8663.

LIMITATIONS OF THE PROCEDURE

- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of influenza type A and B antigens from nasal swab and nasopharyngeal swab specimens.
- The BD Veritor System for Rapid Detection of Flu A+B is capable of detecting both viable and non-viable influenza particles. The BD Veritor System for Rapid Detection of Flu A+B performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Results from the BD Veritor System for Rapid Detection of Flu A+B test should be correlated with the clinical history, epidemiological data and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of influenza A or influenza B infection, and should be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not identify specific influenza A virus subtypes.
- Negative test results are not intended to rule in other non-influenza viral or bacterial infections.
- Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no influenza activity when disease prevalence is low. False negative test results are more likely during peak influenza activity when prevalence of disease is high.
- This device has been evaluated for use with human specimen material only.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, influenza A viruses that have undergone minor amino acid changes in the target epitope region.

- The analytical reactivity of this device has not been established for avian or swine origin influenza strains other than those included in the "strain reactivity" tables.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection.
- The BD Veritor System Instrument reports dual positive influenza A and influenza B results as "Result Invalid." True dual positives are exceptionally rare. Specimens generating a "Result Invalid" should be retested. Upon retesting, if the specimen produces a "Result Invalid" the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus.

EXPECTED VALUES

The rate of positivity observed in respiratory testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, the time of year, age of the patient, geographic location and most importantly, local disease prevalence. The overall prevalence observed with an FDA-cleared Influenza A and B molecular assay in the U.S. during the 2010–2011 clinical study was 29.9% for Influenza A and 19.7% for influenza B. The overall prevalence observed with the same FDA-cleared Influenza A and B molecular assay in Japan during the 2010–2011 clinical study was 32.2% for Influenza A and 31.7% for influenza B.

PERFORMANCE CHARACTERISTICS

Clinical Performance:

Performance characteristics for the **BD Veritor System for Rapid Detection of Flu A+B** test were established in multi-center Point-of-Care (POC) studies conducted at five U.S. trial sites and eight Japan trial sites during the 2010–2011 respiratory season. A total of 736 prospective specimens (515 in the U.S. and 221 in Japan) were tested using the BD Veritor System for Rapid Detection of Flu A+B test. These specimens consisted of nasal and nasopharyngeal swabs from symptomatic patients. In the U.S., 54% of the samples were from females and 46% from males. 20.3% of the samples were from patients less than or equal to 5 years of age, 40.8% were from patients in the 6-21 year age group, 35.6% were from 22–59 years of age, and the remaining 3.3% were obtained from people greater than or equal to age 60. In Japan, 43.3% of the samples were from females and 56.7% from males. 27.3% of the samples were from patients less than or equal to 5 years of age, 58.4% were from patients in the 16–21 year age group, 13.1% from 22–59 years of age, and 1.3% were obtained from people greater than or equal to age 60.

The performance of the BD Veritor System for Rapid Detection of Flu A+B test at the U.S. sites were compared to an FDA-cleared Influenza A and B molecular assay (PCR).

Explanation of Terms:

- PPA: Positive Percent Agreement = $a / (a+c) \times 100\%$
 NPA: Negative Percent Agreement = $d / (b+d) \times 100\%$
 P: Positive
 N: Negative
 C.I: Confidence Interval

New Test Method	Comparator Method	
	P	N
P	a	b
N	c	d
Total	(a+c)	(b+d)

The performance is presented in Table 1 through Table 4 below.

Table 1: Summary of the Performance Data for the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for All Swabs - All Sites

POC: BD Flu A	Reference PCR		
	P	N	Total
P	189	13	202
N	37	497	534
Total	226	510	736
Reference Method: PCR PPA: 83.6% (76.1%, 89.1%) NPA: 97.5% (95.7%, 98.5%)			

POC: BD Flu B	Reference PCR		
	P	N	Total
P	139	10	149
N	32	555	587
Total	171	565	736
Reference Method: PCR PPA: 81.3% (71.1%, 88.5%) NPA: 98.2% (95.7%, 99.3%)			

Wald 95% Confidence intervals corrected for over-dispersion, where needed, due to potential variability between sites.

Table 2: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for All Swabs - U.S. Sites

POC: BD Flu A	Reference PCR		
	P	N	Total
P	122	8	130
N	33*	352	385
Total	155	360	515
Reference Method: PCR PPA: 78.7% (95% C.I.: 71.6–84.4%) NPA: 97.8% (95% C.I.: 95.7–98.9%)			

POC: BD Flu B	Reference PCR		
	P	N	Total
P	75	2	77
N	26**	412	438
Total	101	414	515
Reference Method: PCR PPA: 74.3% (95% C.I.: 65–81.8%) NPA: 99.5% (95% C.I.: 98.3–99.9%)			

*Of the 33 PCR positive, BD Veritor negative Influenza A specimens, eight were positive in the BD Veritor assay using a second swab specimen (reference method specimen) collected from the same patient.

**Of the 26 PCR positive, BD Veritor negative Influenza B specimens, six were positive in the BD Veritor assay using a second swab specimen (reference method specimen) collected from the same patient.

Table 3: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Nasopharyngeal Swabs - U.S. Sites

POC: BD Flu A	Reference PCR		
	P	N	Total
P	53	5	58
N	18	135	153
Total	71	140	211
Reference Method: PCR PPA: 74.6% (95% C.I.: 63.4–83.3%) NPA: 96.4% (95% C.I.: 91.9–98.5%)			

POC: BD Flu B	Reference PCR		
	P	N	Total
P	22	1	23
N	8	180	188
Total	30	181	211
Reference Method: PCR PPA: 73.3% (95% C.I.: 55.6–85.8%) NPA: 99.4% (95% C.I.: 96.9–99.9%)			

Table 4: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Nasal Swabs – U.S. Sites

POC: BD Flu A	Reference PCR		
	P	N	Total
P	69	3	72
N	15	217	232
Total	84	220	304
Reference Method: PCR PPA: 82.1% (95% C.I.: 72.6–88.9%) NPA: 98.6% (95% C.I.: 96.1–99.5%)			

POC: BD Flu B	Reference PCR		
	P	N	Total
P	53	1	54
N	18	232	250
Total	71	233	304
Reference Method: PCR PPA: 74.6% (95% C.I.: 63.4–83.3%) NPA: 99.6% (95% C.I.: 97.6–99.9%)			

The performance of the BD Veritor System for Rapid Detection of Flu A+B test at the Japan sites were also compared to the same FDA-cleared Influenza A and B molecular assay (PCR) and are presented in Table 5 through Table 7.

Table 5: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for all Swabs - Japan Sites

POC: BD Flu A	Reference PCR		
	P	N	Total
P	67	5	72
N	4	145	149
Total	71	150	221
Reference Method: PCR PPA: 94.4% (95% C.I.: 86.4–97.8%) NPA: 96.7% (95% C.I.: 92.4–98.6%)			

POC: BD Flu B	Reference PCR		
	P	N	Total
P	64	8	72
N	6	143	149
Total	70	151	221
Reference Method: PCR PPA: 91.4% (95% C.I.: 82.5–96%) NPA: 94.7% (95% C.I.: 89.9–97.3%)			

Table 6: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for all Nasopharyngeal Swabs - Japan Sites

POC: BD Flu A	Reference PCR		
	P	N	Total
P	30	1	31
N	2	83	85
Total	32	84	116
Reference Method: PCR PPA: 93.8% (95% C.I.: 79.9–98.3%) NPA: 98.8% (95% C.I.: 93.6–99.8%)			

POC: BD Flu B	Reference PCR		
	P	N	Total
P	38	2	40
N	1	75	76
Total	39	77	116
Reference Method: PCR PPA: 97.4% (95% C.I.: 86.8–99.5%) NPA: 97.4% (95% C.I.: 91.0–99.3%)			

Table 7: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Nasal Swabs – Japan Sites

POC: BD Flu A	Reference PCR		
	P	N	Total
P	37	4	41
N	2	62	64
Total	39	66	105
Reference Method: PCR			
PPA: 94.9% (95% C.I.: 83.1–98.6%)			
NPA: 93.9% (95% C.I.: 85.4–97.6%)			

POC: BD Flu B	Reference PCR		
	P	N	Total
P	26	6	32
N	5	68	73
Total	31	74	105
Reference Method: PCR			
PPA: 83.9% (95% C.I.: 67.4–92.9%)			
NPA: 91.9% (95% C.I.: 83.4–96.2%)			

Reproducibility

The reproducibility of the BD Veritor System for Rapid Detection of Flu A+B test was evaluated at three POC sites. The reproducibility panel was composed of 30 simulated influenza A or B samples. These included moderate positive samples, low positive samples (near the assay limit of detection), high negative samples (i.e., containing very low concentrations of virus such that positive results occur ~5% of the time) and negative samples. The panel was tested by two operators at each site for five consecutive days. The results are summarized below.

Reproducibility Results – Percent of Flu A Positives				
Sample	Site 1	Site 2	Site 3	Total
High negative H1N1 A	0.0% (0/30) (95% C.I.: 0.0–11.3%)	10% (3/30) (95% C.I.: 3.5–25.6%)	26.7% (8/30) (95% C.I.: 14.2–44.4%)	12.2% (11/90) (95% C.I.: 7.0–20.6%)
Low positive H1N1 A	86.7% (26/30) (95% C.I.: 70.3–94.7%)	96.7% (29/30) (95% C.I.: 83.3–99.4%)	100% (30/30) (95% C.I.: 88.6–100%)	94.4% (85/90) (95% C.I.: 87.6–97.6%)
Moderate positive H1N1 A	100% (30/30) (95% C.I.: 88.6–100%)	100% (30/30) (95% C.I.: 88.6–100%)	100% (30/30) (95% C.I.: 88.6–100%)	100% (90/90) (95% C.I.: 95.9–100%)
High negative H3N2 A	0.0% (0/30) (95% C.I.: 0.0–11.3%)	10% (3/30) (95% C.I.: 3.5–25.6%)	16.7% (5/30) (95% C.I.: 7.3–33.6%)	8.9% (8/90) (95% C.I.: 4.6–16.6%)
Low positive H3N2 A	100% (30/30) (95% C.I.: 88.6–100%)	93.3% (28/30) (95% C.I.: 78.7–98.2%)	96.7% (29/30) (95% C.I.: 83.3–99.4%)	96.7% (87/90) (95% C.I.: 90.7–98.9%)
Moderate positive H3N2 A	100% (30/30) (95% C.I.: 88.6–100%)	100% (30/30) (95% C.I.: 88.6–100%)	100% (30/30) (95% C.I.: 88.6–100%)	100% (90/90) (95% C.I.: 95.9–100%)
Negatives	0.0% (0/119) (95% C.I.: 0.0–3.1%)	0.8% (1/119) (95% C.I.: 0.1–4.6%)	0.0% (0/119) (95% C.I.: 0.0–3.1%)	0.3% (1/357) (95% C.I.: 0.0–1.6%)

Reproducibility Results – Percent of Flu B Positives				
Sample	Site 1	Site 2	Site 3	Total
High negative B	0.0% (0/30) (95% C.I.: 0.0–11.3%)	3.3% (1/30) (95% C.I.: 0.6–16.7%)	26.7% (8/30) (95% C.I.: 14.2–44.4%)	10% (9/90) (95% C.I.: 5.4–17.9%)
Low positive B	73.3% (22/30) (95% C.I.: 55.6–85.8%)	90% (27/30) (95% C.I.: 74.4–96.5%)	90% (27/30) (95% C.I.: 74.4–96.5%)	84.4% (76/90) (95% C.I.: 75.6–90.5%)
Moderate positive B	100% (29/29) (95% C.I.: 88.3–100%)	96.6% (28/29) (95% C.I.: 82.8–99.4%)	100% (29/29) (95% C.I.: 88.3–100%)	98.9% (86/87) (95% C.I.: 93.8–99.8%)
Negatives	0.0% (0/210) (95% C.I.: 0.0–1.8%)	1.0% (2/210) (95% C.I.: 0.3–3.4%)	0.0% (0/210) (95% C.I.: 0.0–1.8%)	0.3% (2/630) (95% C.I.: 0.1–1.2%)

Analytical Studies

Analytical Sensitivity (Limit of Detection)

The limit of detection (LOD) for the BD Veritor System for Rapid Detection of Flu A+B test was established for a total of 8 influenza strains: 5 influenza A and 3 influenza B. The LOD for each strain represents the lowest concentration producing a positivity rate of ≥95% based on testing 20 to 60 replicates.

Type	Influenza Viral Strain	Calculated LOD (TCID ₅₀ /mL)	Calculated LOD (EID ₅₀ /mL)	No. Positive / Total	% Positive
A	A/Brisbane/10/2007 H3N2	7.27 x 10 ²	N/A	57/60	95%
A	A/Brisbane/59/2007 H1N1	3.30 x 10 ²	N/A	57/60	95%
A	A/California/7/2009 H1N1	5.00 x 10 ³	N/A	57/60	95%
A	A/Victoria/3/75 H3N2	3.11 x 10 ³	N/A	59/60	98.3%
A	A/Anhui/1/2013 H7N9	N/A	5.42 x 10 ⁶	59/60	98.3%
B	B/Brisbane/60/2008	7.42 x 10 ³	N/A	58/60	96.7%
B	B/Florida/4/2006	1.30 x 10 ³	N/A	58/60	96.7%
B	B/Lee/40	4.44 x 10 ⁴	N/A	20/20	100%

TCID₅₀/mL = 50% Tissue Culture Infectious Dose
EID₅₀/mL = 50% Egg Infectious Dose

Strain Reactivity with Influenza A and B Viruses

The BD Veritor System for Rapid Detection of Flu A+B test was evaluated using a panel of influenza strains. Each strain was diluted and tested in triplicate until a point where not all of the replicates were positive. The dilution prior to that is provided in the table below as a minimal detected concentration. All influenza A strains showed positive Flu A test results and negative Flu B test results. Conversely, all of the influenza B strains showed positive Flu B test results and negative Flu A test results. Although this test has been shown to detect novel avian influenza A (H7N9) and H3N2v cultured viruses the performance characteristics of this device with clinical specimens that are positive for novel avian influenza A (H7N9) and H3N2v influenza viruses has not been established. The BD Veritor System Flu A+B assay can distinguish between influenza A and B viruses, but it cannot differentiate influenza A subtypes.

Strain	Subtype	Minimal Detected Concentration
A/Brisbane/59/2007	H1N1	3.3 x 10 ² TCID ₅₀ /mL*
A/California/7/2009	H1N1	5.0 x 10 ³ TCID ₅₀ /mL*
A/Denver/1/57	H1N1	4.45 x 10 ⁴ CEID ₅₀ /mL
A/FM/1/47	H1N1	7.91 x 10 ⁴ CEID ₅₀ /mL
A/Fujian-Gulou/1896/2009	H1N1	4.5 x 10 ⁵ CEID ₅₀ /mL
A/Mal/302/54	H1N1	2.22 x 10 ⁵ CEID ₅₀ /mL
A/New Caledonia/20/1999	H1N1	2.5 x 10 ³ TCID ₅₀ /mL
A/New Jersey/8/76	H1N1	1.58 x 10 ³ CEID ₅₀ /mL
A/NWS/33	H1N1	1.58 x 10 ⁴ CEID ₅₀ /mL
A/PR/8/34	H1N1	6.31 x 10 ² TCID ₅₀ /mL
A/Solomon Island/03/2006	H1N1	2.5 x 10 ³ TCID ₅₀ /mL
A/Washington/24/2012	H1N1	3.16 x 10 ⁴ EID ₅₀ /mL
A/Weiss/43	H1N1	7.03 x 10 ⁶ CEID ₅₀ /mL
A/WS/33	H1N1	7.91 x 10 ² CEID ₅₀ /mL
A/Aichi/2/68	H3N2	7.91 x 10 ³ CEID ₅₀ /mL
A/Brisbane/10/2007	H3N2	7.27 x 10 ² TCID ₅₀ /mL*
A/California/02/2014	H3N2	1.45 x 10 ² TCID ₅₀ /mL
A/Hong Kong/8/68	H3N2	8.89 x 10 ⁴ CEID ₅₀ /mL
A/Moscow/10/99	H3N2	5.8 x 10 ⁶ TCID ₅₀ /mL
A/Perth/16/2009	H3N2	1.0 x 10 ⁶ TCID ₅₀ /mL
A/Port Chalmers/1/73	H3N2	3.95 x 10 ⁴ CEID ₅₀ /mL
A/Switzerland/9715293/2013	H3N2	3.25 x 10 ² TCID ₅₀ /mL
A/Texas/50/2012	H3N2	1.75 x 10 ³ TCID ₅₀ /mL
A/Wisconsin/67/2005	H3N2	2.5 x 10 ⁵ TCID ₅₀ /mL
A/Victoria/3/75	H3N2	3.11 x 10 ³ TCID ₅₀ /mL*
A/Indiana/08/2011	H3N2v	1 x 10 ⁴ TCID ₅₀ /mL
A/Indiana/10/2011	H3N2v	7.9 x 10 ⁶ CEID ₅₀ /mL
A/Kansas/13/2009	H3N2v	1.0 x 10 ³ TCID ₅₀ /mL
A/Minnesota/11/2010	H3N2v	7.9 x 10 ⁵ CEID ₅₀ /mL
A/Pennsylvania/14/2010	H3N2v	1.26 x 10 ⁶ CEID ₅₀ /mL
A/West Virginia/06/2011	H3N2v	7.9 x 10 ³ TCID ₅₀ /mL
A/Anhui/01/2005	H5N1	0.512 HA
A/Vietnam/1203/2004	H5N1	0.512 HA
A/Northern Pintail/Washington/40964/2014	H5N2	6.28 x 10 ⁵ EID ₅₀ /mL
A/Pheasant/NewJersey/1355/1998	H5N2	0.256 HA
A/Gyrfalcon/Washington/41088-6/2014	H5N8	1.98 x 10 ⁶ EID ₅₀ /mL
A/Mallard/Netherlands/12/2000	H7N7	0.256 HA
A/Anhui/1/2013	H7N9	5.42 x 10 ⁶ CEID ₅₀ /mL*
A/Chicken/HongKong/G9/1997	H9N2	1.024 HA

*Values taken from preceding Analytical Limit of Detection Table

a. EID₅₀ = 50% Egg Infectious Dose

b. TCID₅₀ = 50% Tissue Culture Infectious Dose

c. CEID₅₀ = 50% Chicken Embryo Infectious Dose

d. HA = Hemagglutination Assay

Strain	Minimal Detected Concentration
B/Brazil/178/96	2.32 x 10 ⁴ TCID ₅₀ /mL
B/Brisbane/33/2008 (Victoria Lineage)	2.45 x 10 ⁵ CEID ₅₀ /mL
B/Brisbane/60/2008	7.42 x 10 ³ TCID ₅₀ /mL*
B/Brisbane/72/97	1.00 x 10 ⁴ TCID ₅₀ /mL
B/Canada/548/99	>0.64 HA
B/Egypt/393/99	>1.28 HA
B/Florida/2/2006	1.08 x 10 ⁴ TCID ₅₀ /mL
B/Florida/4/2006	1.30 x 10 ³ TCID ₅₀ /mL*
B/Fujian/93/97	3.95 x 10 ⁵ TCID ₅₀ /mL
B/Fukushima/220/99	9.33 x 10 ² TCID ₅₀ /mL
B/Guangdong-Liwan/1133/2014 (Yamagata Lineage)	9.0 x 10 ⁵ CEID ₅₀ /mL
B/GuangXi/547/98	2.32 x 10 ⁵ TCID ₅₀ /mL
B/Hawaii/01/97	>6.4 HA
B/Hong Kong/5/72	1.11 x 10 ⁴ CEID ₅₀ /mL
B/Hong Kong/219/98	>1 HA
B/Hong Kong/259/2010(Victoria Lineage)	1.35 x 10 ⁶ CEID ₅₀ /mL
B/Jiangsu/10/2003	1.16 x 10 ⁴ TCID ₅₀ /mL
B/Johannesburg/5/99	3.95 x 10 ⁴ TCID ₅₀ /mL
B/Lee/40	4.44 x 10 ⁴ CEID ₅₀ /mL*
B/Lisbon/03/96	>0.08 HA
B/Malaysia/2506/2004	5.0 x 10 ⁴ TCID ₅₀ /mL
B/Maryland/1/59	3.51 x 10 ² CEID ₅₀ /mL
B/Massachusetts/2/2012 (Yamagata Lineage)	1.25 x 10 ⁶ CEID ₅₀ /mL
B/Mass/3/66	1.58 x 10 ⁵ CEID ₅₀ /mL
B/Montana/5/2012	3.14 x 10 ⁵ EID ₅₀ /mL
B/Ohio/11/96	>0.16 HA
B/Ohio/1/05	1.34 x 10 ⁵ TCID ₅₀ /mL
B/Phuket/3073/2013	6.08 x 10 ³ TCID ₅₀ /mL
B/Puerto Mont/10427/98	0.02 HA
B/Russia/69	3.9 x 10 ² TCID ₅₀ /mL
B/Shandong/7/97	1.58 x 10 ⁶ TCID ₅₀ /mL
B/Shanghai/04/97	1.58 x 10 ⁵ TCID ₅₀ /mL
B/Shenzhen/135/97	3.16 x 10 ⁴ TCID ₅₀ /mL
B/Sichuan/116/96	0.016 HA
B/Taiwan/2/62	2.81 x 10 ² CEID ₅₀ /mL
B/Texas/06/2011 (Yamagata Lineage)	6.2 x 10 ⁵ CEID ₅₀ /mL
B/Texas/02/2013 (Victoria Lineage)	2.75 x 10 ⁴ CEID ₅₀ /mL
B/Utah/09/2014 (Yamagata Lineage)	6.3 x 10 ³ CEID ₅₀ /mL
B/Victoria/504/00	4.64 x 10 ⁴ TCID ₅₀ /mL
B/Wisconsin/01/2010 (Yamagata Lineage)	7.0 x 10 ² CEID ₅₀ /mL
B/Yamagata/16/88	9.75 x 10 ³ TCID ₅₀ /mL
B/Yamanashi/166/98	4.88 x 10 ⁴ TCID ₅₀ /mL

*Values taken from preceding Analytical Limit of Detection Table

- a. EID₅₀ = 50% Egg Infectious Dose
- b. TCID₅₀ = 50% Tissue Culture Infectious Dose
- c. CEID₅₀ = 50% Chicken Embryo Infectious Dose
- d. HA = Hemagglutination Assay

On January 12, 2017 the US FDA published an order reclassifying antigen based rapid influenza virus detection systems (RIDTs) intended to detect influenza virus antigen directly from clinical specimens from Class I to Class II with special controls. One of these special controls is a requirement for annual analytical reactivity testing of contemporary influenza strains identified by the Centers for Disease control (CDC) using a standardized dilution protocol. Results obtained using the BD Veritor System are available for viewing at bd.com/veritorsystem.

Analytical Specificity (Cross-reactivity)

The BD Veritor System for Rapid Detection of Flu A+B test was evaluated with a total of 51 microorganisms. The 37 bacteria and yeast were tested at a target concentration of approximately 10⁷ CFU/mL (CFU – Colony Forming Units) with the exception of *Staphylococcus aureus*, which was tested at a final concentration of 10⁶ CFU/mL. The 14 viruses were evaluated at concentrations of 10³ to 10¹⁰ TCID₅₀/mL. Of the 51 microorganisms tested, none showed cross-reactivity in either the Flu A or Flu B tests.

<i>Bacteroides fragilis</i>	<i>Neisseria meningitidis</i>
<i>Bordetella pertussis</i>	<i>Neisseria mucosa</i>
<i>Candida albicans</i>	<i>Neisseria sp. (Neisseria perflaus)</i>
<i>Chlamydia pneumoniae</i>	<i>Neisseria subflava</i>
<i>Corynebacterium diphtherium</i>	<i>Peptostreptococcus anaerobius</i>
<i>Escherichia coli</i>	<i>Porphyromonas asaccharolyticus</i>
<i>Fusobacterium nucleatum</i>	<i>Prevotella oralis</i>
<i>Haemophilus influenzae</i>	<i>Propionibacterium acnes</i>
<i>Haemophilus parainfluenzae</i>	<i>Proteus mirabilis</i>
<i>Kingella kingae</i>	<i>Pseudomonas aeruginosa</i>
<i>Klebsiella pneumoniae</i>	<i>Serratia marcescens</i>
<i>Lactobacillus sp.</i>	<i>Staphylococcus aureus</i>
<i>Legionella sp.</i>	<i>Staphylococcus epidermidis</i>
<i>Moraxella catarrhalis</i>	<i>Streptococcus mutans</i>
<i>Mycobacterium tuberculosis</i>	<i>Streptococcus pneumoniae</i>
<i>Mycoplasma pneumoniae</i>	<i>Streptococcus pyogenes</i>
<i>Neisseria gonorrhoeae</i>	<i>Streptococcus sp. Group C</i>

<i>Streptococcus sp. Group G</i>
<i>Streptococcus salivarius</i>
<i>Veillonella parvula</i>
Adenovirus, type 1
Adenovirus, type 7
Cytomegalovirus
Enterovirus
Epstein Barr Virus
HSV Type 1
Human Coronavirus OC43
Human Coronavirus 2229E
Human metapneumovirus (HMPV-27 A2)
Human Parainfluenza
Measles virus
Mumps virus
Respiratory syncytial virus
Rhinovirus

Interfering Substances

Various substances were evaluated with the BD Veritor System for Rapid Detection of Flu A+B test. These substances included whole blood (2%) and various medications. No interference was noted with this assay for any of the substances tested.

Substance	Concentration	Substance	Concentration
4-Acetamidophenol	10 mg/mL	Homeopathic Allergy Medicine	10 mg/mL
Acetylsalicylic acid	20 mg/mL	Ibuprofen	10 mg/mL
Albuterol	0.083 mg/mL	Loratidine	100 ng/mL
Amantadine Hydrochloride	500 ng/mL	Menthol Throat Lozenges	10 mg/mL
Ayr Saline Nasal Gel	10 mg/mL	Mometasone	500 ng/mL
Beclomethasone	500 ng/mL	Mupirocin	500 ng/mL
Budesonide	500 ng/mL	Oseltamivir	500 ng/mL
Chlorpheniramine maleate	5 mg/mL	Oxymetazoline	0.05 mg/mL
Dexamethasone	10 mg/mL	Phenylephrine	1 mg/mL
Dextromethorphan	10 mg/mL	Pseudoephedrine HCl	20 mg/mL
Diphenhydramine HCl	5 mg/mL	Purified Mucin Protein	1 mg/mL
Fexofenadine	500 ng/mL	Ribavirin	500 ng/mL
FluMist	1%	Rimantadine	500 ng/mL
Flunisolide	500 ng/mL	Three OTC mouthwashes	5%
Fluticasone	500 ng/mL	Tobramycin	500 ng/mL
Four OTC nasal sprays	10%	Triamcinolone	500 ng/mL
Four OTC throat drops	25%	Whole Blood	2%
Guaiacol Glyceryl Ether	20 mg/mL	Zanamivir	1 mg/mL

Of the 44 substances tested in this study, none exhibited interfering reactions when tested with influenza A and influenza B positive samples. Based on the data, the substances tested at the indicated concentration levels did not interfere with the BD Veritor System for Rapid Detection of Flu A+B test.

CLIA WAIVER STUDY

As part of a larger prospective study, as described in the Performance Characteristics section above, the accuracy of the BD Veritor System for Rapid Detection of Flu A+B test was evaluated at five CLIA waived testing sites. A total of 31 operators representative of CLIA waived sites (intended users) participated in the study. No training on the use of the test was provided. The study included 515 nasal/nasopharyngeal swabs prospectively collected and 80 retrospective archived specimens. The BD Veritor System results were compared with results obtained by an FDA cleared molecular influenza A and B assay, the comparator method. Three specimens were excluded due to BD Veritor invalid results. The invalid rate was 0.5% (3/598) with 95% C.I: 0.2% to 1.5%.

The positive percent agreement (PPA) and the negative percent agreement (NPA) between the BD Veritor results and the comparator method are presented in the tables below (refer to Performance Characteristics section for definition of terms).

INFLUENZA A				
Positive Percent Agreement and Negative Percent Agreement of BD Veritor Flu A+B Test with the Comparator Method				
Total Number of Samples	PPA	95% Confidence Interval	NPA	95% Confidence Interval
595	82.1% (151/184)	(75.9%, 86.9%)	98.1% (403/411)	(96.2%, 99.0%)

INFLUENZA B				
Positive Percent Agreement and Negative Percent Agreement of BD Veritor Flu A+B Test with the Comparator Method				
Total Number of Samples	PPA	95% Confidence Interval	NPA	95% Confidence Interval
595	79.7% (102/128)	(71.9%, 85.7%)	99.4% (464/467)	(98.1%, 99.8%)

Another study was designed to assess the capability of untrained users to test weakly reactive samples and deliver results with accuracy. The BD Veritor System for Rapid Detection of Flu A+B assay was evaluated at three non-laboratory CLIA waived sites using panels of simulated swab samples including two weak positives near the assay cutoff and one negative sample. The positive swab samples were formulated at two levels: a "low positive" sample targeted at the assay limit of detection; and a "high negative" sample targeted just below the assay limit of detection. The panels included two strains of Flu A viruses (A/California/7/2009 and A/Victoria 3/75) and one Flu B virus (B/Lee/40). The swab samples were randomized and masked with respect to their identity. There were two intended users at each of the CLIA waived sites (six operators in total) and each site tested the panel on each of 10 days. The same panels of simulated swab samples were also tested at three clinical laboratory sites as controls. The performance of the BD Veritor System with samples near the assay cutoff was acceptable when used by intended users.

The tables below show performance of the test with samples near the cutoff of the assay for influenza A and influenza B in the hands of untrained intended users (across all sites).

Influenza A Viral Strains		
Sample Type	Untrained Intended Users	
	Percent Detection	95% Confidence Interval
High Negative A/California/7/2009 H1N1	6.7% (4/60)	(2.6%, 15.9%)
Low Positive A/California/7/2009 H1N1	81.7% (49/60)	(70.1%, 89.4%)
High Negative A/Victoria 3/75 H3N2	6.7% (4/60)	(2.6%, 15.9%)
Low Positive A/Victoria 3/75 H3N2	80.0% (48/60)	(68.2%, 88.2%)
Negative	0.0% (0/118)*	(0.0%, 3.2%)

*Two (2) samples were excluded from the analysis due to errors in data recording.

Influenza B Viral Strain		
	Untrained Intended Users	
Sample Type	Percent Detection	95% Confidence Interval
High Negative B/Lee/40	11.7% (7/60)	(5.8%, 22.2%)
Low Positive B/Lee/40	72.4% (42/58)*	(59.8%, 82.2%)
Negative	0.0% (0/240)	(0.0%, 1.6%)

*Two (2) samples were excluded from the analysis due to errors in data recording.

Using risk analysis as a guide, analytical flex studies were conducted. The studies demonstrated that the test is insensitive to stresses of environmental conditions and potential user errors.

In support of the CLIA waiver, an additional reactivity study was performed at an independent laboratory to demonstrate reactivity of the BD Veritor System for the Rapid Detection of Flu A+B with a broad range of contemporary influenza A and influenza B viruses. The BD Veritor System yielded positive results with all 18 influenza A viruses and 7 influenza B viruses included in the test panel at acceptable viral load levels.

Technical Support

For questions, or to report a problem, please call Technical Support at 1.800.638.8663. Test system problems may also be reported to the FDA using the MedWatch reporting system (phone: 1.800.FDA.1088; fax: 1.800.FDA.1078; or <http://www.fda.gov/medwatch>).

AVAILABILITY

Cat. No. Description

256045	BD Veritor™ System for Rapid Detection of Flu A+B, 30 tests
256051	BD Veritor™ System Flu A+B Control Swab Set, 10 pairs of swabs
220252	COPAN Flexible Minitip Flocked Swab, 100 swabs
256066	BD Veritor™ Plus Analyzer
256068	BD Veritor™ InfoScan Module
443907	USB Printer Cable for BD Veritor™ Analyzer
256038	BD Veritor™ System for Rapid Detection of Respiratory Syncytial Virus (RSV)

To network a BD Veritor Plus Analyzer to an LIS, contact BD Technical Services for details.

REFERENCES

1. Simonsen L., Fukuda K, Schonberger LB, Cox NJ. Impact of influenza epidemics on hospitalizations. *J. Infect. Dis.* 2000; 181:831–7
2. Thompson WW, Shay DK, Weintraub E, et al. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA* 2003; 289:179–86
3. Treanor, J.J., Hayden, F.G., Vrooman, P.S., et al. 2000. Efficacy and safety of the oral neuraminidase inhibitor oseltamivir in treating acute influenza: a randomized controlled trial. *JAMA.* 283:1016–1024.
4. Kaiser, L., Couch, R.B., Galasso, G.J., Glezen, W.P., Webster, R.G., Wright, P.F., and Hayden, F.G. 1999. First international symposium on influenza and other respiratory viruses: summary and overview Kapalua, Maui, Hawaii, December 4–6, 1998. *Antiviral Res.*, 42:149–176
5. Cox, N.J., and Bender, C.A. 1995. The molecular epidemiology of influenza viruses. *Virology*, 6:359–370.
6. Todd, S.J., Minnich, L., and Waner, J.L. 1995. Comparison of rapid immunofluorescence procedure with TestPack RSV and Directigen Flu A for diagnosis of respiratory syncytial virus and influenza A virus. *J. Clin. Microbiol.*33:1650–1651.
7. Harris, P.O. 1989. Clinical relevance and efficient detection of seven major respiratory viruses. *ACL.* p. 15–19.
8. McElhaney, J.E., Gravenstein, S., Krause, P., Hooton, J.W., Upshaw, C.M., and Drinka, P. 1998. Assessment of markers of the cell-mediated immune response after influenza virus infection in frail older adults. *Clin. Diag. Lab. Immunol.* 5:840–844.
9. Fan, J., Henrickson, K.J., and Savatski, L.L. 1998. Rapid simultaneous diagnosis of infections with respiratory syncytial viruses A and B, influenza viruses A and B, and human parainfluenza virus types 1, 2, and 3 by multiplex quantitative reverse transcription-polymerase chain reaction-hybridization assay (hexaplex). *Clin. Infect. Disease* 26:1397–1402.
10. Wright, K.E., Wilson, G.A.R., Novosad, D., Dimock, C., Tan, D., and Weber, J.M. 1995. Typing and subtyping of influenza viruses in clinical samples by PCR. *J. Clin. Microbiol.*33:1180–1184.
11. Kendal, A.P. 1985. Influenza Viruses. p. 341–357. *Laboratory Diagnosis of Viral Infections*, *in* H. Lennette, (ed.) Marcel Dekker, Inc., New York.
12. McQuillen, J., Madeley, C.R., and Kendal, A.P. 1985. Monoclonal antibodies for the rapid diagnosis of influenza A and B virus infections by immunofluorescence. *Lancet*.ii: 911–914.

- Guenther, S.H., and Linnemann, C.C., Jr. 1988. Indirect immunofluorescence assay for rapid diagnosis of influenza virus. *Laboratory Medicine*. 19:581–583
- Minnick, L.L., and Ray, C.G. 1986. Early testing of cell cultures for detection of hemadsorbing viruses. *J. Clin. Microbiol.* 25:421–422.
- Schmidt, N.J., Ota, M., Gallo, D., and Fox, V.L. 1982. Monoclonal antibodies for rapid, strain specific identification of influenza virus isolates. *J. Clin. Microbiol.* 16:763–765.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53–80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EECP). *Office Journal L262*, 17/10/2000, p.021–0045.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or bd.com.

The logo for BD Veritor System, featuring a stylized sunburst icon to the left of the text "BD Veritor System".

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Change History

Revision	Date	Change Summary
14	2018-06	<p>Changes required after FDA review of K180438 and CR 180086. Correcting verbiage throughout, adding Caution bullets to Sections A, B, C and D. Add new Table 1 and renumber tables. Rename Analyzer "the BD Veritor Plus Analyzer" throughout. Reorganize Catalog numbers in Availability section. Add Change History Table.</p>
15	2020-04	<p>Changed name of BD Synapsys™ Microbiology Informatics Solution to BD Synapsys™ Informatics Solution. Additional result documentation capabilities made possible with the implementation of the BD Synapsys Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect (replacing the BD Veritor System Reader and the BD InfoSync Module). In sections C and D, revised Steps 10C and 10D regarding the use of the BD Veritor Plus Analyzer with the BD Veritor InfoScan. Added access information to obtain the document from bd.com/e-labeling.</p>

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US Customers only: For symbol glossary, refer to bd.com/symbols-glossary.

Manufacturer / Производител / Výrobce / Fabrikant / Hersteller / Κατασκευαστής / Fabricante / Tootja / Fabricant / Proizvođač / Gyártó /
Fabricante / Фабрышы / 제조업체 / Gamintojas / Ražotājs / Tilvirker / Producent / Producător / Производител / Výrobca / Proizvođač /
Tillverkare / Üretici / Виробник / 生产厂商

Use by / Исполняйте до / Spolfebuje do / Brug før / Verwendbar bis / Χρήση έως / Usar antes de / Kasutada enne / Date de péremption /
사용 기한 / Uputnjači do / Felhasználásig dátuma / Usare entro / Дейін пайдалануға / Naudokite iki / Izlietot līdz / Houdbaar tot / Brukes
for / Stosować do / Prazo de validade / A se utiliza până la / Исползовать до / Použite do / Använd före / Son kullanna tarini /
Використати долино / 使用截止日期

YYYY-MM-DD / YYYY-MM (MM = end of month)
ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = slutning af måneden)
JJJJ-MM-TT / JJJJ-MM (MM = Monatsende)
EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα)
AAAA-MM-DD / AAAA-MM (MM = fin del mes)
AAAA-KK-PP / AAAA-KK (KK = kuu lõpp)
AAAA-MM-JJ / AAAA-MM (MM = fin du mois)
GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
EEEE-HH-NN / EEEE-HH (HH = hónap utolsó napja)
AAAA-MM-GG / AAAA-MM (MM = finer mese)
ЖОЖК-АА-КК / ЖОЖК-АА (АА = айдын соңы)
YYYY-MM-DD/YYYY-MM (MM = 월말)
MMMM-MM-DD / MMMM-MM (MM = mēnesio pabaiga)
GGGG-MM-DD/GGGG-MM (MM = mēneša beigas)
JJJJ-MM-DD / JJJJ-MM (MM = einde maand)
AAAA-MM-DD / AAAA-MM (MM = sluttan av måneden)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
AAAA-MM-DD / AAAA-MM (MM = fim do mês)
AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii)
TTTT-MM-DD / TTTT-MM (MM = конец месяца)
RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca)
GGGG-MM-DD / GGGG-MM (MM = kraj meseca)
AAAA-MM-DD / AAAA-MM (MM = slutet av månaden)
YYYY-AA-GG / YYYY-AA (AA = ayın sonu)
PPPP-MM-DD / PPPP-MM (MM = кінець місяця)
YYYY-MM-DD / YYYY-MM (MM = 月末)

Catalog number / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Número de catálogo / Kataloginumber / Numéro
catalogue / Kataloški broj / Katalogszám / Numero di catalogo / Каталог номери / 카탈로그 번호 / Katalogo / numeris / Kataloga numurs /
Catalogus nummer / Numer katalogowy / Număr de catalog / Номер по каталогу / Katalogové číslo / Kataloški broj / Katalog numeraris / Номер по
каталогу / 目录号

EC REP

Authorized Representative in the European Community / Оторизиран представител в Европейската общност / Autorizovaný zástupce pro Evropském společenství / Autoriserer representant i De Europæiske Fællesskaber / Autorisierter Vertreter in der Europäischen Gemeinschaft / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Representante autorizado en la Comunidad Europea / Volitatud esindaja Euroopa Nõukogus / Représentant autorisé pour la Communauté européenne / Autorizuirani predstavnik u Evropskoj uniji / Meghatalmazott képviselő az Európai Közösségekben / Representante autorizzato nella Comunità Europea / Европа қауымдастығындағы уәкілетті өкіл / 유권 공동체의 위임 대표 / Įgaliojasis atstovas Europos Bendrijoje / Pilnvarotais pārstāvis Eiropas Kopienā / Bevoegde vertegenwoordiger in de Europese Gemeenschap / Autoriserer representant i EU / Autoryzowane przedstawicielstwo we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Reprezentantul autorizat pentru Comunitatea Europeană / Уполномоченный представитель в Европейском сообществе / Autorizovaný zástupca v Európskom společenstve / Autorizovano predstavništvo u Evropskoj uniji / Auktoriseret representant i Europeiska gemenskapen / Avrupa Toplumları Yetkili Temsilcisi / Уповноважений представник у країнах ЄС / 歐洲共同體授权代表

IVD

In Vitro Diagnostic Medical Device / Медицински уред за диагностика ин витро / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / Dispositivo médico para diagnóstico in vitro / In vitro diagnostiska medeltillsäppläparatur / Dispositif médical de diagnostic in vitro / Medicinska pomagala za In Vitro Dijagnostiku / In vitro diagnostická orovsi eszközök / Dispositivo medicale per diagnostica in vitro / Жасанды жағдайда жұрттегі медициналық диагностика аспабы / In Vitro Diagnostic 의료 기기 / In vitro diagnostikos prietaisai / Medicinas ierices, ko lieto in vitro diagnostika / Medisch hulpmiddel voor in-vitro diagnostiek / In vitro diagnostisk medicinsk utrustning / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Dispositif médical pentru diagnostic in vitro / Медицинский прибор для диагностики ин витро / Medicinska rökka för diagnostik in vitro / Medicinski uređaji za in vitro dijagnostiku / Medicinteknisisk produkt för in vitro-djagnostik / In Vitro Dijagnostisk Tibbi Cihaz / In vitro diagnostisk priester für diagnostik in vitro / 体外诊断医疗设备

Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrensning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperaturi pirang / Limites de température / Dozvoljena temperatura / Hömørsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperaturuulimiet / Temperaturbegrensning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenie teploty / Ograničenje temperature / Temperaturgräns / Sicaklık sınırlaması / Обмеження температури / 温度限制

LOT

Batch Code (Lot) / Код на партидата / Kód (číslo) šarže / Batch-kode (lot) / Batch-Code (Charge) / Кодικός παρτίδας (παρτίδα) / Código de lote (lote) / Partii kod / Numéro de lot / Lot (kod) / Tétel száma (Lot) / Codice batch (lotto) / Топлама коды / 배치 코드(코트) / Partijos numeris (LOT) / Partijas kods (laidiens) / Lot number / Batch-kode (part) / Kod partii (seria) / Código do lote / Cod de serie (Lot) / Код партии (lot) / Kód série (šarža) / Kod serije / Partinummer (Lot) / Parti Kodu (Lot) / Код партії / 批号 (亚批)

Contains sufficient for <=> tests / Содержително е достатъчно за <=> теста / Dostatečné množství pro <=> testů / Indeholder tilstrækkeligt til <=> tests / Ausreichend für <=> Tests / Περιέχει επαρκή ποσότητα για <=> εξετάσεις / Contenido suficiente para <=> pruebas / Küllaldane <=> testide jaoks / Contenu suffisant pour <=> tests / Sadržaj za <=> testova / <=> testzhez elegendő / Contenido suficiente per <=> test / <=> тесттери үшін жеткілікті / <=> 테스트가 충분히 포함됨 / Pakankamas kiekis atlikti <=> testų / Satur pietiekami <=> pārbaudēm / Inhoud voldoende voor "n" testen / Innholder tilstrekkelig til <=> tester / Zawiera ilość wystarczającą do <=> testów / Conteúdo suficiente para <=> testes / Conținut suficient pentru <=> teste / Достатчно до <=> тестов(а) / Obsah vystačí na <=> testov / Sadržaj dovoljan za <=> testova / Innehåller tillräckligt för <=> analyser / <=> test için yeterli miktarda içerir / Вистачить для аналізи: <=> / 足夠進行 <=> 次檢測

Consult Instructions for Use / Направете справка в инструкциите за употреба / Prosdjudite pokyny k použití / Se brugsanvisningen / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Consultar las instrucciones de uso / Lugeda kasutusjuhendit / Consulter la notice d'emploi / Koristi upute za upotrebu / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Паглядзіць інструкцыям тынасьмь апыньца / 사용 지침 참조 / Sicaiklyte naujodumo instrukcijas / Skaitli lietošanas pamācību / Raadpeel de gebruiksaanwijzing / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / См. руководство по эксплуатации / Pozri Pokyny na používanie / Pogledajte uputstvo za upotrebu / Se bruksanvisningen / Kullanim Talimatları'na başvurun / Див. інструкції з використання / 请参阅使用说明

Do not reuse / Не използвайте отново / Neopuzivajte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korraldult / Ne pas réutiliser / Ne koristiti ponovo / Egyyszer használatos / Non riutilizzare / Παιδαλανάβανζ / 재사용 금지 / Tik vienkartiniam naujodumui / Nelielet atkārtoti / Niet opnieuw gebruiken / Kun til engangsbrug / Nie stosować powtórnie / Nao reutilize / Nu refolositi / Не использовать повторно / Neopuzivajte opakovane / Ne upotrebļvajate ponovo / Får eer återanvändas / Tekrar kullanimayın / He wykorzystuwać ponownie / 请勿重复使用

SN

Serial number / Серийн номер / Sériové číslo / Serienummer / Serienummer / Σειριακό αριθμός / N° de serie / Seerianumber / Número de série / Serijski broj / Sorozatszám / Numero di serie / Топтамалық нөмірі / 일련 번호 / Serijos numeris / Sērījas numurs / Serie number / Numer serijny / Número de série / Număr de serie / Серийный номер / Seri numarası / Номер серії / 序列号

For IVD Performance evaluation only / Само за оценка качеството на работа на IVD / Pouze pro vyhodnocení výkonu IVD / Kun til evaluering af IVD ydelse / Nur für IVD-Leistungsbewertungszwecke / Móvo για αξιολόγηση απόδοσης IVD / Sólo para la evaluación del rendimiento en diagnóstico in vitro / Ainult IVD seadme hindamiseks / Réservé à l'évaluation des performances IVD / Samo u znanstvene svrhe za In Vitro Dijagnostiku / Kizárólag in vitro diagnosztikához / Solo per valutazione delle prestazioni IVD / Жасанды жағдайда «пробирка ішінде» диагностикада тек жұмысты бағалау үшін / IVD 성능 평가에 대해적인 사용 / Tik IVD prietais veikimo charakteristikoms tikrinti / Vienigi IVD darbības novērtēšanai / Uitsluitend voor doeltreffendheidsonderzoek / Kun for evaluating av IVD-ydelse / Tytko do oceny wydajności IVD / Uso exclusivo para avaliação de IVD / Numai pentru evaluarea performanței IVD / Только для оценки качества диагностики ин витро / Urcené iba na diagnostiku in vitro / Samo za procenu učinka u in vitro dijagnostici / Endast för utvärdering av diagnostisk användning in vitro / Yalnızca IVD Performans değerlendirilmesi için / Тільки для оцінювання якості діagnostiki in vitro / 仅限 IVD 性能评估

For US: "For Investigational Use Only"

Lower limit of temperature / Долен лимит на температурата / Dolní hranice teploty / Nedre temperaturgrænse / Temperaturuntergrenze / Көптірө ороу термодроқтас / Limite inferior de temperatura / Alumine temperaturuulimiet / Limite inférieure de température / Najnižja dovoljena temperatura / Alsó hőmørsékleti határ / Limite inferior de temperatura / Температураның төменгі рұқсат шері / 하한 온도 / Zemiausia laikymo temperatūra / Temperatūras zemākā robeža / Laagste temperatuurlimiet / Nedre temperaturgrænse / Dolna granica temperatury / Limite mínimo de temperatura / Limită minimă de temperatură / Нижний предел температуры / Spodná hranica teploty / Donja granica temperatury / Nedre temperaturgräns / Sicaklık alt sınırı / Мінімальна температура / 温度下限

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R_x Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

bd.com/e-labeling

Europe, CH, GB, NO:	+800 135 79 135
International:	+31 20 794 7071
AR +800 135 79 135	LT 8800 30728
AU +800 135 79 135	MT +31 20 796 5693
BR 0800 591 1055	NZ +800 135 79 135
CA +1 855 805 8539	RO 0800 895 084
CO +800 135 79 135	RU +800 135 79 135
EE 0800 0100567	SG 800 101 3366
GR 0800 161 22015 7799	SK 0800 606 287
HR 0800 804 804	TR 00800 142 064 866
IL +800 135 79 135	US +1 855 236 0910
IS 800 8996	UY +800 135 79 135
LI +31 20 796 5692	VN 122 80297

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

BD, the BD logo, Synapsys, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2020 BD. All rights reserved.

Insert ends here

DO NOT PRINT THIS PAGE

CLSI document should not be included in print file to vendor.

NOTES:

The CLSI document is to be updated at the same time as the package insert and should be released as a separate PDF document to the web.

The package insert should also be released as a separate PDF document to the web.

The package insert and CLSI document should be released as one document to SAP.

Final art provided to Buyer/print vendor should not include this page or the CLSI portion of this document.

Proofreader: please proof the CLSI document along with the package insert.

Procedure* BD Veritor™ System For Rapid Detection of Flu A+B

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Prepared by	Date Adopted	Supersedes Procedure #

Review Date	Revision Date	Signature

Distributed to	# of Copies	Distributed to	# of Copies

***Any modifications to this document are the sole responsibility of the facility. This “Sample Procedure” is not intended as a substitute for your facility procedure manual, instrument manual, or reagent labeling/package insert. This “Sample Procedure” is intended as a model for use by your facility to meet the needs of your laboratory.**

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Performance characteristics for influenza A and B were established during January through March of 2011 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2010-2011 Season, and Composition of the 2011-2012 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

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All BD Veritor System Flu A+B test devices are interpreted by a BD Veritor System Instrument, either a BD Veritor Reader or BD Veritor Plus Analyzer (the "Analyzer"). When using the BD Veritor Plus Analyzer, workflow steps depend on the selected operational mode and the Analyzer configuration settings. In **Analyze Now mode**, the instrument evaluates assay devices after manual timing of their development. In **Walk Away mode**, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Connection of the Analyzer to a printer is possible if desired. Additional result documentation capabilities are possible with the implementation of the BD Synapsys™ Informatics Solution, and with the addition of the BD Veritor InfoScan module and BD Veritor Plus Connect. Please refer to the Analyzer *Instructions for Use* for details on these features and contact BD technical support for more information.

PRINCIPLES OF THE PROCEDURE

The BD Veritor System for Rapid Detection of Flu A+B is a qualitative, digital immunoassay for the detection of influenza A and B viral antigens in samples processed from respiratory specimens. When specimens are processed and added to the test device, influenza A or B viral antigens bind to anti-influenza antibodies conjugated to detector particles in the A+B test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of antibody on the membrane. A positive result for influenza A is determined by the BD Veritor System Instrument when antigen-conjugate is deposited at the Test "A" position and the Control "C" position on the BD Veritor System Flu A+B assay device.

A positive result for influenza B is determined by the BD Veritor System Instrument when antigen-conjugate is deposited at the Test "B" position and the Control "C" position in the BD Veritor System Flu A+B assay device. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

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1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.
4. Pathogenic microorganisms, including hepatitis viruses, Human Immunodeficiency Virus and novel influenza viruses, may be present in clinical specimens. "Standard Precautions"¹⁶⁻¹⁹ and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used BD Veritor System test devices as biohazardous waste in accordance with federal, state and local requirements.
6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. For optimal results, use the flocked swabs provided with the kit for specimen collection.
8. Other than the flocked swabs that are used for specimen collection, kit components should not make contact with the patient.
9. Do not use kit components beyond the expiration date.
10. Do not reuse the device.
11. Do not use the kit if the Control A+/B- swab and Control B+/A- swab do not yield appropriate results.
12. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
13. To avoid erroneous results, swab specimens must be processed as indicated in the assay procedure section.
14. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.
15. FluMist® is made from attenuated live flu virus and although the concentration tested (1%) was non-interfering, it is possible when tested with higher concentrations that an influenza A and/or influenza B false positive may occur.

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
SPECIMEN COLLECTION

Acceptable specimens for testing with the BD Veritor System Flu A+B test include nasal swabs and nasopharyngeal (NP) swabs. Freshly collected specimens should be processed within 1 hour. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers. Inadequate specimen collection, improper specimen handling and/or transport may yield a false negative result; therefore, specimen collection requires specific training and guidance due to the importance of specimen quality to accurate test results.

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1. The BD Veritor System Kit includes swabs with a flocked tip for nasal specimen collection.
2. Insert the swab into one nostril of the patient.
3. Rotate the swab two complete 360-degree turns; pressing firmly against the nasal mucosa to help ensure the swab obtains both cells and mucus.
4. Withdraw the swab from the nasal cavity. The sample is now ready for processing using the BD Veritor System Kit.

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

Nasal and Nasopharyngeal Swab Test Procedure

NOTES:

- **Reagents, specimens and devices must be at room temperature (15–30 °C) prior to testing.**
- The CLIA-waived BD Veritor System for Rapid Detection of Flu A+B kit is only intended for nasal and nasopharyngeal swab specimens that are collected and tested directly (i.e., dry swabs that have **NOT** been placed in transport media). The kit includes a pre-diluted process reagent in a ready to use "unitized" tube. This CLIA-waived kit IS NOT INTENDED for testing liquid samples such as wash or aspirate samples or swabs in transport media as results can be compromised by over dilution.

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

**Step 9A: Record the Result**

- When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately..

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
Step 6B: Starting Walk Away Mode

- Turn the Analyzer on by pressing the blue power button once
- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
 - **Double-click** the blue power button.

Step 7B: Adding the specimen

- When the display window reads "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY"
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device.

NOTE: Squeeze the tube too close to the tip may cause leakage.

- **CAUTION:** A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.

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
- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

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
Step 6C: Adding the specimen

- Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System Flu A+B device.

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

- Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.
- When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

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

Step 6D: Starting Walk Away mode

- Turn the Analyzer on by pressing the blue power button once.
The display window will briefly show "SCAN CONFIG BARCODE." This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration steps. Ignore this message and postpone this process when an assay is awaiting analysis.
- When the display window reads:
"INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
– **Double-click** the blue power button.

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

- When the display window reads: ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY:
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
 - Gently** squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device. **NOTE: Squeezing the tube close to the tip may cause leakage.**
 - CAUTION: A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.**

Step 9D: Starting the development and reading sequence

- Insert the test device into the slot on the right side of the Analyzer. **The test device must remain horizontal to prevent spilling the specimen out of the sample well.**
- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

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1. Insert the swab all the way into the appropriately labeled RV Reagent D tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
2. Continue processing the swab according to the Test Procedure for Nasal and Nasopharyngeal swabs above beginning at Step 4 "**Remove the swab.**"

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The limit of detection (LOD) for the BD Veritor System for Rapid Detection of Flu A+B test was established for a total of 8 influenza strains: 5 influenza A and 3 influenza B. The LOD for each strain represents the lowest concentration producing a positivity rate of $\geq 95\%$ based on testing 20 to 60 replicates.

Type	Influenza Viral Strain	Calculated LOD (TCID ₅₀ /mL)	Calculated LOD (EID ₅₀ /mL)	No. Positive / Total	% Positive
A	A/Brisbane/10/2007 H3N2	7.27×10^2	N/A	57/60	95%
A	A/Brisbane/59/2007 H1N1	3.30×10^2	N/A	57/60	95%
A	A/California/7/2009 H1N1	5.00×10^3	N/A	57/60	95%
A	A/Victoria/3/75 H3N2	3.11×10^3	N/A	59/60	98.3%
A	A/Anhui/1/2013 H7N9	N/A	5.42×10^6	59/60	98.3%
B	B/Brisbane/60/2008	7.42×10^3	N/A	58/60	96.7%
B	B/Florida/4/2006	1.30×10^3	N/A	58/60	96.7%
B	B/Lee/40	4.44×10^4	N/A	20/20	100%

TCID₅₀/mL = 50% Tissue Culture Infectious Dose

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The BD Veritor System for Rapid Detection of Flu A+B test was evaluated with a total of 51 microorganisms. The 37 bacteria and yeast were tested at a target concentration of approximately 10^7 CFU/mL (CFU – Colony Forming Units) with the exception of *Staphylococcus aureus*, which was tested at a final concentration of 10^6 CFU/mL. The 14 viruses were evaluated at concentrations of 10^3 to 10^{10} TCID₅₀/mL. Of the 51 microorganisms tested, none showed cross-reactivity in either the Flu A or Flu B tests.

<i>Bacteroides fragilis</i>	<i>Neisseria meningitidis</i>
<i>Bordetella pertussis</i>	<i>Neisseria mucosa</i>
<i>Candida albicans</i>	<i>Neisseria sp. (Neisseria perflaus)</i>
<i>Chlamydia pneumoniae</i>	<i>Neisseria subflava</i>
<i>Corynebacterium diphtherium</i>	<i>Peptostreptococcus anaerobius</i>
<i>Escherichia coli</i>	<i>Porphyromonas asaccharolyticus</i>
<i>Fusobacterium nucleatum</i>	<i>Prevotella oralis</i>
<i>Haemophilus influenzae</i>	<i>Propionibacterium acnes</i>
<i>Haemophilus parainfluenzae</i>	<i>Proteus mirabilis</i>
<i>Kingella kingae</i>	<i>Pseudomonas aeruginosa</i>
<i>Klebsiella pneumoniae</i>	<i>Serratia marcescens</i>
<i>Lactobacillus sp.</i>	<i>Staphylococcus aureus</i>
<i>Legionella sp.</i>	<i>Staphylococcus epidermidis</i>
<i>Moraxella catarrhalis</i>	<i>Streptococcus mutans</i>
<i>Mycobacterium tuberculosis</i>	<i>Streptococcus pneumoniae</i>
<i>Mycoplasma pneumoniae</i>	<i>Streptococcus pyogenes</i>
<i>Neisseria gonorrhoeae</i>	<i>Streptococcus sp. Group C</i>

<i>Streptococcus sp. Group G</i>
<i>Streptococcus salivarius</i>
<i>Veillonella parvula</i>
Adenovirus, type 1
Adenovirus, type 7
Cytomegalovirus
Enterovirus
Epstein Barr Virus
HSV Type 1
Human Coronavirus OC43
Human Coronavirus 2229E
Human metapneumovirus (HMPV-27 A2)
Human Parainfluenza
Measles virus
Mumps virus
Respiratory syncytial virus
Rhinovirus

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REFERENCES

1. Simonsen L., Fukuda K, Schonberger LB, Cox NJ. Impact of influenza epidemics on hospitalizations. *J. Infect. Dis.* 2000; *181*:831–7
2. Thompson WW, Shay DK, Weintraub E, et al. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA* 2003; *289*:179–86
3. Treanor, J.J., Hayden, F.G., Vrooman, P.S., et al. 2000. Efficacy and safety of the oral neuraminidase inhibitor oseltamivir in treating acute influenza: a randomized controlled trial. *JAMA.* *283*:1016–1024.
4. Kaiser, L., Couch, R.B., Galasso, G.J., Glezen, W.P., Webster, R.G., Wright, P.F., and Hayden, F.G. 1999. First international symposium on influenza and other respiratory viruses: summary and overview Kapalua, Maui, Hawaii, December 4–6, 1998. *Antiviral Res.*, *42*:149–176
5. Cox, N.J., and Bender, C.A. 1995. The molecular epidemiology of influenza viruses. *Virology*, *6*:359–370.
6. Todd, S.J., Minnich, L., and Waner, J.L. 1995. Comparison of rapid immunofluorescence procedure with TestPack RSV and Directigen Flu A for diagnosis of respiratory syncytial virus and influenza A virus. *J. Clin. Microbiol.**33*:1650–1651.
7. Harris, P.O. 1989. Clinical relevance and efficient detection of seven major respiratory viruses. *ACL.* p. 15–19.
8. McElhaney, J.E., Gravenstein, S., Krause, P., Hooton, J.W., Upshaw, C.M., and Drinka, P. 1998. Assessment of markers of the cell-mediated immune response after influenza virus infection in frail older adults. *Clin. Diag. Lab. Immunol.* *5*:840–844.
9. Fan, J., Henrickson, K.J., and Savatski, L.L. 1998. Rapid simultaneous diagnosis of infections with respiratory syncytial viruses A and B, influenza viruses A and B, and human parainfluenza virus types 1, 2, and 3 by multiplex quantitative reverse transcription-polymerase chain reaction-hybridization assay (hexaplex). *Clin. Infect. Disease* *26*:1397–1402.
10. Wright, K.E., Wilson, G.A.R., Novosad, D., Dimock, C., Tan, D., and Weber, J.M. 1995. Typing and subtyping of influenza viruses in clinical samples by PCR. *J. Clin. Microbiol.**33*:1180–1184.
11. Kendal, A.P. 1985. Influenza Viruses. p. 341–357. *Laboratory Diagnosis of Viral Infections, In H. Lennette, (ed.) Marcel Dekker, Inc., New York.*
12. McQuillen, J., Madeley, C.R., and Kendal, A.P. 1985. Monoclonal antibodies for the rapid diagnosis of influenza A and B virus infections by immunofluorescence. *Lancet*.*ii*: 911–914.
13. Guenther, S.H., and Linnemann, C.C., Jr. 1988. Indirect immunofluorescence assay for rapid diagnosis of influenza virus. *Laboratory Medicine.**19*:581–583
14. Minnick, L.L., and Ray, C.G. 1986. Early testing of cell cultures for detection of hemadsorbing viruses. *J. Clin. Microbiol.* *25*:421–422.
15. Schmidt, N.J., Ota, M., Gallo, D., and Fox, V.L. 1982. Monoclonal antibodies for rapid, strain specific identification of influenza virus isolates. *J. Clin. Microbiol.**16*:763–765.
16. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
17. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.**17*:53–80.
18. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S Government Printing Office, Washington, D.C.
19. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EECP). *Office Journal L262, 17/10/2000, p.021–0045.*

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