# BD Veritor™ System

## For Rapid Detection of SARS-CoV-2

500058119(01)  
2021-09  
English

 256089 BD Veritor™ System for Rapid Detection of SARS-CoV-2, 30 Test

 Pages 1–13

 Pages 14 – 27

**Kit configured for testing nasal swab samples freshly collected, processed and dispensed directly onto assay test device.**

For use with the BD Veritor™ Plus Analyzer running firmware version 5.4 or later.

### INTENDED USE

The BD Veritor™ System for Rapid Detection of SARS-CoV-2 is a chromatographic digital immunoassay intended for the direct and qualitative detection of SARS-CoV-2 nucleocapsid antigens in direct anterior nasal swabs from individuals who are either suspected of COVID-19 by their healthcare provider, within the first 5 days of the onset of symptoms, or from individuals without symptoms or other epidemiological reasons to suspect COVID-19 when tested twice over two or three days with at least 24 hours, and no more than 48 hours, between tests. In the United States, testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet the requirements to perform moderate, high, or waived complexity tests. This test is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

Results are for the identification of SARS-CoV-2 nucleocapsid antigen. This antigen is generally detectable in upper respiratory samples during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. Laboratories within the United States and its territories are required to report all results to the appropriate public health authorities.

Negative results should be treated as presumptive, do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary, for patient management.

For serial testing programs, additional confirmatory testing with a molecular test for negative results may be necessary, if there is a high likelihood of SARS-CoV-2 infection, such as in an individual with a close contact with COVID-19 or with suspected exposure to COVID-19 or in communities with high prevalence of infection. Additional confirmatory testing with a molecular test for positive results may also be necessary, if there is a low likelihood of SARS-CoV-2 infection, such as in individuals without known exposures to SARS-CoV-2 or residing in communities with low prevalence of infection.

The BD Veritor™ System for Rapid Detection of SARS-CoV-2 is intended for use in point of care settings by laboratory personnel and healthcare providers appropriately trained in the use of the BD Veritor™ System.

### SUMMARY AND EXPLANATION OF THE TEST

A novel coronavirus (2019-nCoV) was identified in December 2019,<sup>1</sup> which has resulted in confirmed human infections worldwide. Cases of severe illness and deaths have been reported. On February 11, 2020 the International Committee for Taxonomy of Viruses (ICTV) renamed the virus SARS-CoV-2.

The median incubation time is estimated to be approximately 5 days<sup>2</sup> with symptoms estimated to be present within 12 days of infection. The symptoms of COVID-19 are similar to other viral respiratory diseases and include fever, cough, and shortness of breath.

The BD Veritor™ System for Rapid Detection of SARS-CoV-2 is a rapid (approximately 15 minutes) chromatographic digital immunoassay for the direct detection of the presence or absence of SARS-CoV-2 antigens in respiratory specimens taken from patients with signs and symptoms who are suspected of COVID-19. The test is intended for interpretation in both laboratory and near patient testing environments only with the BD Veritor™ Plus Analyzer Instrument. The test is not intended to be interpreted visually. Procedures to evaluate test devices depend on the BD Veritor™ Plus Analyzer workflow configuration chosen. In **Analyze Now** mode, the instrument evaluates assay devices after manual timing of their development. In **Walk Away mode**, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Additionally, connection of a BD Veritor™ Plus Analyzer to a printer or IT system is possible if desired. Additional result documentation capabilities are possible with the integration of a BD Veritor™ barcode scanning enabled module. Please refer to the BD Veritor™ Plus Analyzer Instructions for Use for details on how to implement these features.

## PRINCIPLES OF RHE PROCEDURE

The BD Veritor™ System consists of a dedicated opto-electronic interpretation instrument and immunochromatographic assays for the qualitative detection of antigens from pathogenic organisms in samples processed from respiratory specimens. When specimens are processed and added to the test device, SARS-CoV-2 antigens present in the specimen bind to antibodies conjugated to detector particles in the test strip. The antigen-conjugate complexes migrate across the test strip to the reaction area and are captured by a line of antibodies bound on the membrane. A positive result is determined by the BD Veritor™ Plus Analyzer when antigen-conjugate is deposited at the Test "T" position and the Control "C" position on the assay device. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective result.

## REAGENTS

The following components are included in the BD Veritor™ System for Rapid Detection of SARS-CoV-2 kit.

### Materials Provided:

| KIT COMPONENT                   | QUANTITY                                                                                                           | DESCRIPTION                                                                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------|--------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                 | 256089                                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                               |
| BD Veritor™ System Test Devices | 30 single use test devices                                                                                         | Foil pouched test device containing one reactive strip. Each strip has one line of murine anti-SARS coronavirus monoclonal antibody on the test line, and one of biotin coupled to bovine protein on the positive control line.<br><br>Murine and Leporine anti-SARS coronavirus and anti-biotin monoclonal antibodies conjugated to detector reagents are bound in the sample delivery area. |
| Extraction Reagent              | 30 single use reaction tubes, each with 325 µL extraction reagent and having an integral dispensing tip            | Detergent solution with less than 0.1% sodium azide (preservative).                                                                                                                                                                                                                                                                                                                           |
| Specimen sampling swabs         | 30 sterile, single use specimen sampling swabs                                                                     | For sample collection and transfer.                                                                                                                                                                                                                                                                                                                                                           |
| SARS-CoV-2 (+) Control Swab     | 1 each – individually wrapped for single use                                                                       | Non-infectious, recombinant viral protein antigen with less than 0.1% sodium azide.                                                                                                                                                                                                                                                                                                           |
| SARS-CoV-2 (–) Control Swab     | 1 each – individually wrapped for single use                                                                       | Buffer with less than 0.1% sodium azide.                                                                                                                                                                                                                                                                                                                                                      |
| Assay documentation             | 1 each – Instructions for use<br>1 each – Quick reference instruction card<br>1 each – Nasal sampling instructions |                                                                                                                                                                                                                                                                                                                                                                                               |

| MATERIALS REQUIRED BUT NOT PROVIDED                                                                                                                                                                                                                                                                 | OPTIONAL EQUIPMENT                                                                                                                                                                                                                                |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"><li>• BD Veritor™ Plus Analyzer (Catalog Number 256066)</li><li>• BD Veritor™ System Barcode Scanning Module (Catalog Number 256068 or 445010)*</li><li>• Timer</li><li>• Tube rack for specimens</li><li>• Any necessary personal protective equipment</li></ul> | <ul style="list-style-type: none"><li>• USB Printer cable for BD Veritor™ Plus Analyzer (Catalog Number 443907)</li><li>• Epson Printer model TM-T20 II</li><li>• BD Veritor™ Plus Connect (contact your BD representative for details)</li></ul> |

\* If necessary to configure instrument display language

## WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use. Do not reuse the test device or components.
2. This test has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.
3. Do not use this kit beyond the expiration date printed on the outside carton.
4. Do not use the kit to evaluate patient specimens if either the positive control swab or negative control swab fail to give expected results.
5. Test results are not meant to be visually determined. All test results must be determined using the BD Veritor™ Plus Analyzer.
6. To avoid erroneous results, specimens must be processed as indicated in the assay procedure section.
7. Do not mix components from different kits or from other BD diagnostic assays, even if they look similar.
8. Do not use extraction reagents from other BD Veritor™ assay kits as assay performance could be affected.
9. When collecting a nasal swab sample, use the nasal swab supplied in the kit.

- Other than the swabs used for specimen collection, kit components should not make contact with the patient.
- Proper specimen collection, handling, and processing are critical to the performance of this test.
- Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures. Wear protective clothing such as laboratory coats, disposable gloves, and eye protection when specimens are collected and evaluated.
- Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. Standard precautions and institutional guidelines should always be followed in handling, storing, and disposing of all specimens and all items contaminated with blood or other body fluids.
- The SARS-CoV-2 positive control swabs have been prepared from recombinant viral proteins and do not contain infectious material.
- Dispose of all BD Veritor™ kit components (test devices, reagents, and control swabs) as biohazardous waste in accordance with federal, state, and local requirements.
- Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. If there is contact with skin, wash immediately with plenty of water. Contact with acids produces very toxic gas. Do not flush reagents down the drain.
- Test devices used in a laminar flow hood or in areas with high air flow should be covered during test development to ensure proper sample flow.
- For additional information on hazard symbols, safety, handling and disposal of the components within this kit, please refer to the Safety Data Sheet (SDS) located at [bd.com](http://bd.com).

**STORAGE**

Kits must be stored at 2–30 °C. **DO NOT FREEZE. Reagents and devices must be at room temperature (15–30 °C) when used for testing.**

**SPECIMEN COLLECTION AND HANDLING**

Reagents, specimens and devices must be at room temperature (15–30 °C) for testing.

Acceptable specimens for testing with this kit only include **nasal swab specimens obtained by the dual nares collection method collected and tested directly (i.e., swabs that have NOT been placed in transport media). This kit IS NOT INTENDED for testing liquid samples such as wash or aspirate samples or swabs in transport media as results can be compromised by over dilution.** It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early during symptom onset will contain the highest viral titers; specimens obtained after five days of symptoms are more likely to produce negative results when compared to an RT-PCR assay. Inadequate specimen collection, improper specimen handling and/or transport may yield a falsely negative result; therefore, training in specimen collection is highly recommended due to the importance of specimen quality for generating accurate test results.

**Nasal Swab Specimen Collection**

The BD Veritor™ System for Rapid Detection of SARS-CoV-2 includes swabs for nasal specimen collection.

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ol style="list-style-type: none"> <li>Insert swab into one nostril of the patient. The swab tip should be inserted up to 2.5 cm (1 inch) from the edge of the nostril. Roll the swab 5 times along the mucosa inside the nostril to ensure that both mucus and cells are collected. Take approximately 15 seconds to collect the sample.</li> <li>Using the same swab, repeat this process for the other nostril to ensure that an adequate sample is collected from both nasal cavities.</li> <li>Withdraw the swab from the nasal cavity. The sample is now ready for processing using the BD Veritor™ System SARS-CoV-2 kit. The swab should be processed in the extraction reagent vial within 1 hour.</li> </ol> | <div data-bbox="429 820 979 966"> <div>1</div>  <div>2</div>  <div>3</div>  </div> <div data-bbox="429 982 979 1177"> <p><b>DO'S AND DON'TS OF SPECIMEN COLLECTION</b></p> <ul style="list-style-type: none"> <li>• Use only swabs provided with the kit.</li> <li>• Do test sample immediately and always within 1 hour of collection.</li> </ul> <p>For laboratory support for COVID-19 in the EU/EEA, visit <a href="https://www.ecdc.europa.eu/en/novel-coronavirus/laboratory-support">https://www.ecdc.europa.eu/en/novel-coronavirus/laboratory-support</a>.<br/>Outside the United States, refer to applicable guidelines from other national or local authorities.</p> </div> |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|


TEST PROCEDURE

Getting ready to test

The following steps assume that the BD Veritor™ Plus Analyzer is ready to use. To choose or change any BD Veritor™ Plus Analyzer settings, see the BD Veritor™ Plus Analyzer Instructions for Use. A printer is not necessary to display results. However, if your facility has chosen to connect the BD Veritor™ Plus Analyzer to a printer, check that the BD Veritor™ Plus Analyzer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

Once the nasal swab has been collected from the nostrils, the swab must be processed within 1 hour.

Procedural steps for Nasal Swabs or control swabs:

|   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|---|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | <ul style="list-style-type: none"><li>Remove one extraction reagent tube/tip and one BD Veritor™ System for Rapid Detection of SARS-CoV-2 test device from its foil pouch immediately before testing.</li><li>Label one test device and one extraction reagent tube for each specimen or control to be tested.</li><li>Place the labeled extraction reagent tube(s) in a rack in the designated area of the workspace.</li></ul>                                                                                  |
| 2 | Remove and discard the cap from the extraction reagent tube. The extraction reagent is formulated for use with this kit. Do not use reagent tubes from other BD Veritor™ kits as assay performance could be impacted.                                                                                                                                                                                                                                                                                             |
| 3 | Insert the swab into the tube and plunge the swab up and down in the fluid for a minimum of 15 seconds, taking care not to splash contents out of the tube.                                                                                                                                                                                                                                                                                                                                                       |
| 4 | Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.                                                                                                                                                                                                                                                                                                                                                                                                                        |
| 5 | <p>Press the attached tip firmly onto the extraction reagent tube containing the processed sample (threading or twisting is not required). Mix thoroughly by swirling or flicking the bottom of the tube.</p> <p><b>NOTE:</b> Do not use tubes or tips from any other product, including other products from BD or other manufacturers.</p> <p><b>Once the swab has been processed in the extraction reagent and the tube has been capped, the sample must be added to the test device within 30 minutes.</b></p> |


| After step 5, choose from the BD Veritor™ Plus Analyzer workflow option below before continuing to step 6: |                                                      |                                                    |                                                                                |                          |
|------------------------------------------------------------------------------------------------------------|------------------------------------------------------|----------------------------------------------------|--------------------------------------------------------------------------------|--------------------------|
|                                                                                                            | BD Veritor™ Plus Analyzer in <b>Analyze Now</b> mode | BD Veritor™ Plus Analyzer in <b>Walk Away</b> mode | BD Veritor™ Plus Analyzer with the BD Veritor™ Barcode Scanning Enabled Module |                          |
|                                                                                                            |                                                      |                                                    | in <b>Analyze Now</b> mode                                                     | in <b>Walk Away</b> mode |
| Instructions in section:                                                                                   | <b>A</b>                                             | <b>B</b>                                           | <b>C</b>                                                                       | <b>D</b>                 |

For batch testing instructions, please refer to [https://bdveritor.bd.com/content/dam/bdveritor/pdfs/728\\_US\\_0121\\_BD\\_Veritor\\_COVID\\_Analyze\\_Now\\_Batch\\_Testing\\_Guide\\_WEB.pdf](https://bdveritor.bd.com/content/dam/bdveritor/pdfs/728_US_0121_BD_Veritor_COVID_Analyze_Now_Batch_Testing_Guide_WEB.pdf)

## A Using a BD Veritor™ Plus Analyzer in “Analyze Now” mode\*:

- 6A Adding the specimen to the test device**
- Invert the extraction reagent tube and hold it vertically (approximately 1 inch above the sample well).
  - Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well.
  - Excess volume remains for retesting if necessary.

**NOTE:** Squeezing the tube too close to the tip may cause leakage.

- 7A Timing test development**
- After adding the sample, allow the test to run for **15 minutes but no longer than 20 minutes** before inserting the test device into the BD Veritor™ Plus Analyzer.
  - During incubation time, turn the BD Veritor™ Plus Analyzer on by pressing the blue power button once.

**NOTE:** If running test under laminar flow hood, cover test device to avoid inconsistent flow.

- 8A Using the BD Veritor™ Plus Analyzer**
- The BD Veritor™ Plus Analyzer will complete a self-test before it is ready for use. After the self-test the display window shows “INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE”. During incubation time, turn the BD Veritor™ Plus Analyzer on by pressing the blue power button once.
  - INSERT THE TEST DEVICE when the 15-minute assay development time is complete.
  - The status of the assay analysis process appears in the display window. Follow the on-screen prompts to complete the procedure. Do not touch the instrument or remove the test device until the result appears.
  - When analysis is complete, the test result appears in the display window.

- 9A Record the result before removing the test device.**

**\*ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the BD Veritor™ Plus Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

## B Using the BD Veritor™ Plus Analyzer in “Walk Away” mode\*: with no barcode scanning module installed


To use Walk Away mode - connect the AC power adapter to the Analyzer and a power source

- 6B Starting Walk Away Mode**
- Turn the BD Veritor™ Plus Analyzer on by pressing the blue power button once
  - When the display window reads: “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”, Double-click the blue power button.
  - The display window reads “ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY”.

- 7B Adding the specimen to the test device**
- Invert the extraction reagent tube and hold it vertically (approximately 1 inch above the sample well).
  - Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well.
  - Excess volume remains for retesting if necessary.


**NOTE:** Squeezing the tube too close to the tip may cause leakage.

**CAUTION:** A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.

|                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                  |
|-------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------|
| 8B                                                                                                                                              | <p><b>Starting the development and reading sequence</b></p> <ul style="list-style-type: none"><li>Insert the test device into the slot on the right side of the BD Veritor™ Plus Analyzer.</li></ul> <p><b>The test device must remain horizontal to prevent spilling the specimen out of the sample well.</b></p> <ul style="list-style-type: none"><li>“DO NOT DISTURB TEST IN PROGRESS” appears in the display window. Automatic timing of the assay development, image processing and result analysis begins. The status of the assay analysis process appears in the display window. Follow the on-screen prompts to complete the procedure. Do not touch the instrument or remove the test device until the result appears.</li><li>The display window shows the remaining analysis time.</li></ul> |  |
| <p><b>Do not touch the BD Veritor™ Plus Analyzer or remove the test device during this process. Doing so will abort the assay analysis.</b></p> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                  |
| 9B                                                                                                                                              | <p><b>Record the result before removing the test device.</b></p> <ul style="list-style-type: none"><li>When analysis is complete, the test result appears in the display window. Record the result and discard the test device appropriately.</li></ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                  |

**\*ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the BD Veritor™ Plus Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

|           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                                                                                   |
|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| <b>C</b>  | <b>Using the BD Veritor™ Plus Analyzer In “Analyze Now” mode with a barcode scanning module installed</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                   |
| <b>6C</b> | <b>Adding the specimen to the test device</b> <ul style="list-style-type: none"> <li>Invert the extraction reagent tube and hold it vertically (approximately 1 inch above the sample well).</li> <li>Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well.</li> <li>Excess volume remains for retesting if necessary.</li> </ul> <p><b>NOTE:</b> Squeezing the tube too close to the tip may cause leakage.</p>                                                                                                                                                                                                                                                                                                                                   |  |
| <b>7C</b> | <b>Timing test development</b> <ul style="list-style-type: none"> <li>Allow the test to develop for <b>15 minutes</b>.</li> </ul> <p><b>Caution:</b> Incorrect results may occur if development time is less than 15 minutes. Some lines may appear on the device sooner. Do not read device visually.</p> <ul style="list-style-type: none"> <li>If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.</li> </ul>                                                                                                                                                                                                                                                                                                                         |  |
| <b>8C</b> | <b>Using the BD Veritor™ Plus Analyzer</b> <p>During the incubation time, turn on the BD Veritor™ Plus Analyzer by pressing the blue button once.</p> <p>The display window briefly shows “SCAN CONFIG BARCODE”. This is an opportunity to change the configuration of the BD Veritor™ Plus Analyzer. Ignore this message and postpone this process when an assay is awaiting analysis. Please refer to the BD Veritor™ Plus Analyzer Instructions for Use for configuration steps.</p> <ul style="list-style-type: none"> <li>When assay development time is complete and the BD Veritor™ Plus Analyzer display window reads “INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE”, insert the BD Veritor™ System SARS-CoV-2 device into the slot on the right side of the BD Veritor™ Plus Analyzer.</li> </ul> |                                                                                   |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                 |                                                                                  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------|
| <b>9C</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | <b>Using the barcode scanner</b> <ul style="list-style-type: none"> <li>Follow the prompts on the display screen to complete any required barcode scans of: <ul style="list-style-type: none"> <li>— OPERATOR ID</li> <li>— SPECIMEN ID and/or</li> <li>— KIT LOT NUMBER</li> </ul> </li> </ul> |  |
| <ul style="list-style-type: none"> <li>Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the BD Veritor™ Plus Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new reading sequence.</li> <li>Move barcodes slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.</li> <li>The BD Veritor™ Plus Analyzer can record the Kit Lot Number and expiration date in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user.</li> </ul> |                                                                                                                                                                                                                                                                                                 |                                                                                  |
| <p>After required scans are completed, the BD Veritor™ Plus Analyzer displays a countdown timer and test analysis begins.</p> <ul style="list-style-type: none"> <li><b>Do not touch the BD Veritor™ Plus Analyzer or remove the test device during this process. Doing so will abort the assay analysis.</b></li> <li>When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.</li> </ul> <p><b>If the printer is not connected, record the result before removing the assay device.</b></p>                                                                                                       |                                                                                                                                                                                                                                                                                                 |                                                                                  |

**ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the BD Veritor™ Plus Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

|                                                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>10C</b><br><br> | <b>Remove the test device</b><br>Remove and then discard the test device appropriately. The display will show "INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE" to indicate the BD Veritor™ Plus Analyzer is ready to perform another test. If the BD Veritor™ Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. If a network connection is not detected while the ENVELOPE symbol is still displayed, the BD Veritor™ Plus Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted. |
|----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## **D Using the BD Veritor™ Plus Analyzer In "Walk Away" mode with a barcode scanning module installed**


**To use Walk Away mode - connect the AC power adapter to the BD Veritor™ Plus Analyzer and a power source**

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                     |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>6D</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | <b>Starting Walk Away Mode</b> <ul style="list-style-type: none"> <li>Turn the BD Veritor™ Plus Analyzer on by pressing the blue power button once. The display window will briefly show "SCAN CONFIG BARCODE". This is an opportunity to change the configuration of the BD Veritor™ Plus Analyzer. Please refer to the BD Veritor™ Plus Analyzer Instructions for Use for configuration steps. Ignore this message and postpone this process when an assay is awaiting analysis. When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE", Double-click the blue power button.</li> <li>When the display window reads "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE", double-click the blue power button.</li> </ul> |   |
| <b>7D</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | <b>Using the barcode scanner</b> <ul style="list-style-type: none"> <li>Follow the prompts on the display screen to complete any required barcode scans of: <ul style="list-style-type: none"> <li>— OPERATOR ID</li> <li>— SPECIMEN ID and/or</li> <li>— KIT LOT NUMBER</li> </ul> </li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                               |  |
| <ul style="list-style-type: none"> <li>Prompts for each scanning step appear in the display window for only 30 seconds. Failure to complete scans during that time will cause the BD Veritor™ Plus Analyzer to default to the beginning of step 6D. To restart this step, remove and reinsert the test device to initiate a new reading sequence.</li> <li>Move barcodes slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.</li> <li>The BD Veritor™ Plus Analyzer can record the Kit Lot Number and expiration date in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user.</li> </ul> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                     |

**8D Adding the specimen to the test device**

- When the display window reads: "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY": "WALK AWAY MODE", double-click the blue power button.
- Invert the tube, holding it vertically (approximately 1 inch above the BD Veritor™ System SARS-CoV-2 device sample well).
- Gently squeeze the ridged portion of the tube, dispensing three (3) drops of the processed specimen into the sample well.
- Excess volume remains for retesting if necessary.

**NOTE:** Squeezing the tube too close to the tip may cause leakage.  
**CAUTION:** A countdown timer displays the time remaining for test insertion. Walk Away mode must be activated again when this timer expires. Confirm timer is visible and Walk Away mode is activated before inserting test device.

**9D Starting the development and reading sequence**

- Insert the test device into the slot on the right side of the BD Veritor™ Plus Analyzer. The test device must remain horizontal to prevent spilling the specimen out of the sample well.
- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- The display window shows the remaining analysis time.

**Do not touch the BD Veritor™ Plus Analyzer or remove the test device during this process. Doing so will abort the assay analysis.**

- When analysis is complete, a result appears in the display window. If configured to display, the Specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

**If the printer is not connected, record the result before removing the assay device.**

**\*ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the BD Veritor™ Plus Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).**

**10D Remove the test device**

Remove and then discard the test device appropriately. The display will show "INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE" to indicate the BD Veritor™ Plus Analyzer is ready to perform another test.

If the BD Veritor™ Plus Analyzer is connected to an LIS, a steady ENVELOPE symbol will appear to indicate that results are awaiting transmission. If a network connection is not detected while the ENVELOPE symbol is still displayed, the BD Veritor™ Plus Analyzer will queue all untransmitted results and attempt to transmit them when reconnected. If it is powered off during this time, it will attempt to transmit as soon as power is restored, and connection is re-established. A flashing envelope indicates that data are in the process of being transmitted.

## INTERPRETATION OF RESULTS

The BD Veritor™ Plus Analyzer (provided separately) must be used for interpretation of all test results. Do not read results visually.

| Display         | Interpretation                                                 |
|-----------------|----------------------------------------------------------------|
| CoV2: +         | Positive Test for SARS-CoV-2 (antigen detected)                |
| CoV2: -         | Presumptive Negative Test for SARS-CoV-2 (no antigen detected) |
| CONTROL INVALID | Test Invalid.* Repeat the test.                                |

\*Invalid Test – If the test is invalid, the BD Veritor™ System Instrument will display "CONTROL INVALID" and the test or control must then be repeated. If the "CONTROL INVALID" reading recurs, contact BD.

## REPORTING OF RESULTS

**Positive Test** – Positive for the presence of SARS-CoV-2 antigen. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

**Negative Test** – Negative results are presumptive. Negative test results do not preclude infection and should not be used as the sole basis for treatment or other patient management decisions, including infection control decisions, particularly in the presence of clinical signs and symptoms consistent with COVID-19, or in those who have been in contact with the virus. It is recommended that these results be confirmed by a molecular testing method, if necessary, for patient management.

For serial testing programs, additional confirmatory testing with a molecular test for negative results may be necessary after second negative result for asymptomatic patients, if there is a high likelihood of SARS-CoV-2 infection, such as in an individual with a close contact with COVID-19 or with suspected exposure to COVID-19 or in communities with high prevalence of infection. Additional confirmatory testing with a molecular test for positive results may also be necessary, if there is a low likelihood of SARS-CoV-2 infection, such as in individuals without known exposures to SARS-CoV-2 or residing in communities with low prevalence of infection.

**Control Invalid** – Do not report results. Repeat the test. It may be necessary to collect a fresh patient specimen.

## QUALITY CONTROL

Each BD Veritor™ System SARS-CoV-2 test device contains both positive and negative internal/procedural controls:

- The internal positive control line validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
- The membrane area surrounding test lines functions as a background check on the assay device.

The BD Veritor™ System Instrument evaluates the internal/procedural controls after insertion of each test device. The BD Veritor™ Plus Analyzer prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result. NOTE: The internal controls do not assess proper sample collection technique.

## EXTERNAL POSITIVE AND NEGATIVE CONTROLS

Positive and Negative control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the BD Veritor™ System Instrument perform as expected. Prepare kit control swabs and test using the same procedure as used for patient specimens.

BD recommends controls be run once for:

- each new kit lot,
- each new operator,
- as required by internal quality control procedures and in accordance with local and federal regulations or accreditation requirements.

**If the kit controls do not perform as expected, do not report patient results. Contact your local BD representative.**

## LIMITATIONS OF THE PROCEDURE

- The performance of this test has not yet been clinically validated for use in patients without signs and symptoms of respiratory infection, or for serial screening applications and performance may differ in these populations.
- Clinical performance was evaluated with frozen samples, and test performance may be different with fresh samples.
- Users should test specimens as quickly as possible after specimen collection.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other non-SARS-CoV-2 viral or bacterial infections.
- Results from the BD Veritor™ System for Rapid Detection of SARS-CoV-2 test should be correlated with the clinical history, epidemiological data, and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or processed improperly; therefore, a negative test result does not eliminate the possibility of SARS-CoV-2 infection.
- The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5 of illness are more likely to be negative compared to a RT-PCR assay.
- Failure to follow the test procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of SARS-CoV-2 antigens from nasal swab specimens only.
- The BD Veritor™ System for Rapid Detection of SARS-CoV-2 can detect both viable and non-viable SARS-CoV-2 material. The BD Veritor™ System for Rapid Detection of SARS-CoV-2 performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
- The performance of this device has not been assessed in a population vaccinated against COVID-19.
- The performance of the device has not been assessed on specimens from individuals who have been infected with emerging variants of SARS-CoV-2 of public health concern.
- This device has been evaluated for use with human specimen material only.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no SARS-CoV-2 activity when disease prevalence is low. False negative test results are more likely when prevalence of disease caused by SARS-CoV-2 is high.
- Sensitivity of the test after the first five days of the onset of symptoms has been demonstrated to decrease as compared to a RT-PCR SARS-CoV-2 assay.
- Negative results should be treated as presumptive and confirmed with an authorized molecular assay cleared for diagnostic use in region of use.
- Users should test specimens as quickly as possible after specimen collection, within 1 hour after specimen collection and within 30 minutes of placing the swab into the extraction reagent.
- The validity of the BD Veritor™ System for Rapid Detection of SARS-CoV-2 test has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
- Based on *in vitro* testing, false positive or invalid results cannot be ruled out in patients with rheumatoid factor higher than 12.5 IU/mL in nasal fluid, although it is unknown if such concentrations are clinically relevant.

CLINICAL PERFORMANCE

The performance of the BD Veritor™ System for Rapid Detection of SARS-CoV-2 was evaluated in symptomatic individuals. Supplemental clinical data was submitted for the BD Veritor™ System for the Rapid Detection of SARS-CoV-2 and Flu A+B, in which the clinical performance did exceed the 95% confidence interval lower bound of 70% to support serial testing.

As with all antigen tests, performance may decrease as days since symptom onset increases, due to lower viral loads later in the patient's disease course. The inability to synchronize asymptomatic individuals with onset of infection may impact performance as specimens may be tested when viral loads are below the assay's limit of detection.

The performance of this test has not yet been clinically validated for use in patients without signs and symptoms of respiratory infection or for serial screening applications, and performance may differ in these populations. Clinical studies in asymptomatic patients.

Clinical studies in asymptomatic patients undergoing serial testing are ongoing to establish the clinical performance. Samples were collected by qualified personnel in 21 geographically diverse areas across the United States. Nasal swabs were collected following the dual nares method and handled as described in the package insert of the collection device. Specimens were frozen within 30 minutes of collection and stored until tested. All specimens within a pre-specified date range were selected and then sequentially tested in a blinded fashion. The performance of the BD Veritor™ System Assay was compared to results of a nasopharyngeal or oropharyngeal swab stored in 3 mL viral transport media tested with an Emergency Use Authorized molecular (RT-PCR) test for detection of SARS-CoV-2.

Performance was established with 319 direct nasal swabs prospectively collected and enrolled from individual symptomatic patients who were suspected of COVID-19 (within 5 days of onset of two or more self-reported symptoms).<sup>a</sup>

<sup>a</sup> Symptoms included new loss of taste or smell, fever, shortness of breath or difficulty breathing, headache, cough, sore throat, muscle pain, chills and repeated shaking with chills.

Summary of the Performance of the BD Veritor™ System for Rapid Detection of SARS-CoV-2 Compared to RT-PCR for Nasal Swabs in Symptomatic Individuals.

| BD Veritor™ Analyzer Results | Reference PCR Results |     |       |
|------------------------------|-----------------------|-----|-------|
|                              | POS                   | NEG | Total |
| POS                          | 51                    | 1   | 52    |
| NEG                          | 5                     | 262 | 267   |
| Total                        | 56                    | 263 | 319   |

PPA: 91.1% (C.I. 80.7%–96.1%)

NPA: 99.6% (C.I. 97.9%–99.9%)

OPA: 98.1% (C.I. 96.0%–99.1%)

PPV: 98.1% (C.I. 90.7%–99.9%)

NPV: 98.1% (C.I. 96.0%–99.4%)

EXPLANATION OF TERMS:

C.I.: Confidence Interval

PPA: Positive Percent Agreement = True Positives / (True Positives + False Negatives)

NPA: Negative Percent Agreement = True Negatives / (True Negatives + False Positives)

OPA: Overall Percent Agreement = (True Positives + True Negatives) / Total Samples

PPV: Positive Predictive Value = True Positives / (True Positives + False Positives)

NPV: Negative Predictive Value = True Negatives / (True Negatives + False Negatives)

ANALYTICAL PERFORMANCE

LIMIT OF DETECTION (ANALYTICAL SENSITIVITY)

The LOD for the BD Veritor™ System for Rapid Detection of SARS-CoV-2 was established using limiting dilutions of a chemically inactivated viral sample. The material was supplied at a concentration of 6.3 x 10<sup>6</sup> TCID<sub>50</sub>/mL. In this study, designed to estimate the LOD of the assay when using a direct nasal swab, the starting material was spiked into a volume of pooled human nasal matrix obtained from healthy donors and confirmed negative for SARS-CoV-2. At each dilution, 50 µL samples were added to swabs and then tested in the BD Veritor™ assay using the procedure appropriate for patient nasal swab specimens. The lowest concentration that demonstrated ≥95% positivity was set as the LOD.

| Starting Material Concentration              | Estimated LOD                                | No. Positive/Total | % Positive |
|----------------------------------------------|----------------------------------------------|--------------------|------------|
| 6.3 x 10 <sup>6</sup> TCID <sub>50</sub> /mL | 3.9 x 10 <sup>3</sup> TCID <sub>50</sub> /mL | 20/20              | 100%       |

## CROSS-REACTIVITY (ANALYTICAL SPECIFICITY)

Cross-reactivity of the BD Veritor™ System for Rapid Detection of SARS-CoV-2 was evaluated by testing a panel of high prevalence respiratory pathogens that could potentially cross-react with the BD Veritor™ System for Rapid Detection of SARS-CoV-2. Each organism and virus was tested in triplicate. The final concentration of each organism is documented in the following table.

| Potential Cross-Reacting Organism         | Concentration Tested                         | Cross-Reactivity (Yes/No) |
|-------------------------------------------|----------------------------------------------|---------------------------|
| Human coronavirus 229E (heat inactivated) | 1.0 x 10 <sup>5</sup> U/mL                   | No                        |
| Human coronavirus OC43                    | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Human coronavirus NL63                    | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Adenovirus                                | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Human Metapneumovirus                     | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Parainfluenza virus 1                     | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Parainfluenza virus 2                     | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Parainfluenza virus 3                     | 5.2 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Parainfluenza virus 4                     | 1.6 x 10 <sup>4</sup> TCID <sub>50</sub> /mL | No                        |
| Influenza A                               | 2.5 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Influenza B                               | 2.9 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Enterovirus                               | 4.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Respiratory syncytial virus               | 4.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| Rhinovirus                                | 1.1 x 10 <sup>5</sup> PFU/mL                 | No                        |
| SARS-coronavirus                          | 4.5 x 10 <sup>5</sup> PFU/mL                 | No                        |
| MERS-coronavirus                          | 1.5 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                        |
| <i>Haemophilus influenzae</i>             | 1.4 x 10 <sup>6</sup> CFU/mL                 | No                        |
| <i>Streptococcus pneumoniae</i>           | 1.0 x 10 <sup>6</sup> CFU/mL                 | No                        |
| <i>Streptococcus pyogenes</i>             | 1.6 x 10 <sup>6</sup> CFU/mL                 | No                        |
| <i>Candida albicans</i>                   | 1.8 x 10 <sup>6</sup> CFU/mL                 | No                        |
| Pooled human nasal wash                   | 100%                                         | No                        |
| <i>Bordetella pertussis</i>               | 1.4 x 10 <sup>6</sup> CFU/mL                 | No                        |
| <i>Mycoplasma pneumoniae</i>              | 1.0 x 10 <sup>6</sup> CFU/mL                 | No                        |
| <i>Chlamydia pneumoniae</i>               | 1.0 x 10 <sup>6</sup> IFU/mL                 | No                        |
| <i>Legionella pneumophila</i>             | 1.0 x 10 <sup>6</sup> CFU/mL                 | No                        |

To estimate the likelihood of cross-reactivity with SARS-CoV-2 of organisms that were not available for wet testing, *in silico* analysis using the Basic Local Alignment Search Tool (BLAST) managed by the National Center for Biotechnology Information (NCBI) was used to assess the degree of protein sequence homology.

- For *P. jirovecii* one area of sequence similarity shows 45.4% homology across 9% of the sequence, making cross-reactivity in the BD Veritor™ sandwich immunoassay highly unlikely.
- No protein sequence homology was found between SARS-CoV-2 and *M. tuberculosis*, and thus homology-based cross-reactivity can be ruled out.
- The comparison between SARS-CoV-2 nucleocapsid protein and human coronavirus HKU1 revealed that the only potential for homology is with the HKU1 nucleocapsid phosphoprotein. Homology is relatively low, at 36.7% across 82% of sequences, but cross-reactivity cannot be ruled out.

## ENDOGENOUS INTERFERING SUBSTANCES

Various substances were evaluated with the BD Veritor™ System for Rapid Detection of SARS-CoV-2. The substances tested included whole blood 4%, mucin and various medications. No interference was noted with this assay for any of the substances tested.

| Substance                                                         | Concentration Tested | Interference (Yes/No) |
|-------------------------------------------------------------------|----------------------|-----------------------|
| Afrin Nasal Spray (Oxymetazoline)                                 | 5% v/v               | No                    |
| Flonase (Fluticasone)                                             | 5% v/v               | No                    |
| Nasacort (Triamcinolone)                                          | 5% v/v               | No                    |
| Neo-Synephrine (Phenylephrine hydrochloride)                      | 5% v/v               | No                    |
| Oseltamivir                                                       | 2.2 µg/mL            | No                    |
| Mucin protein                                                     | 2.5 mg/mL            | No                    |
| Rhinocort (Budesonide)                                            | 5% v/v               | No                    |
| Saline nasal spray                                                | 15% v/v              | No                    |
| Zanamivir                                                         | 282 ng/mL            | No                    |
| Zicam Cold Remedy (Galphimia glauca, Luffa operculata, Sabadilla) | 5% v/v               | No                    |
| Whole blood                                                       | 4% v/v               | No                    |
| Cepacol (Menthol/Benzocaine)                                      | 1.5 mg/mL            | No                    |
| Ricola (menthol)                                                  | 1.5 mg/mL            | No                    |
| Tobramycin                                                        | 4 µg/mL              | No                    |
| Sucrets (Dyclonine/Menthol)                                       | 1.5 mg/mL            | No                    |
| NeilMed Naso Gel                                                  | 5% v/v               | No                    |
| Zicam nasal spray (Oxymetazoline)                                 | 10% v/v              | No                    |
| Alkalol nasal wash                                                | 10% v/v              | No                    |
| Fisherman's Friend (menthol)                                      | 1.5 mg/mL            | No                    |
| Chloraseptic (Phenol Spray)                                       | 15% v/v              | No                    |
| Mupirocin                                                         | 10 mg/mL             | No                    |

## MICROBIAL INTERFERENCE

The BD Veritor™ System for Rapid Detection of SARS-CoV-2 assay was evaluated with various organisms at the concentrations indicated below. No interference was noted.

| Potential Microbial Interferent | Concentration Tested                         | Interference (Yes/No) |
|---------------------------------|----------------------------------------------|-----------------------|
| Human coronavirus 229E          | 1.0 x 10 <sup>5</sup> U/mL                   | No                    |
| Human coronavirus OC43          | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Human coronavirus NL63          | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Adenovirus                      | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Human Metapneumovirus           | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Parainfluenza virus 1           | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Parainfluenza virus 2           | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Parainfluenza virus 3           | 5.2 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Parainfluenza virus 4a          | 1.5 x 10 <sup>4</sup> TCID <sub>50</sub> /mL | No                    |
| Influenza A                     | 2.5 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Influenza B                     | 2.9 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Enterovirus D68                 | 4.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| Respiratory syncytial virus     | 4.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |

| Potential Microbial Interferent | Concentration Tested                         | Interference (Yes/No) |
|---------------------------------|----------------------------------------------|-----------------------|
| Rhinovirus 3                    | 1.1 x 10 <sup>5</sup> PFU/mL                 | No                    |
| SARS-coronavirus                | 4.5 x 10 <sup>5</sup> PFU/mL                 | No                    |
| MERS-coronavirus                | 1.5 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No                    |
| <i>Haemophilus influenzae</i>   | 1.4 x 10 <sup>6</sup> CFU/mL                 | No                    |
| <i>Streptococcus pneumoniae</i> | 1.0 x 10 <sup>6</sup> CFU/mL                 | No                    |
| <i>Streptococcus pyogenes</i>   | 1.6 x 10 <sup>6</sup> CFU/mL                 | No                    |
| <i>Bordetella pertussis</i>     | 1.4 x 10 <sup>6</sup> CFU/mL                 | No                    |
| <i>Mycoplasma pneumoniae</i>    | 1.0 x 10 <sup>6</sup> CFU/mL                 | No                    |
| <i>Chlamydia pneumoniae</i>     | 1.0 x 10 <sup>6</sup> CFU/mL                 | No                    |
| <i>Legionella pneumophila</i>   | 1.0 x 10 <sup>6</sup> CFU/mL                 | No                    |
| Pooled human nasal wash         | N/A                                          | No                    |
| <i>Candida albicans</i>         | 1.8 x 10 <sup>6</sup> CFU/mL                 | No                    |

## REPRODUCIBILITY

Another study was designed to assess the capability of users to test seeded swab samples across the range of the assay with three (3) users, over three (3) days, with three (3) lots of devices. The following table shows the performance.

| Sample                            | Operator #1     |                    | Operator #2     |                    | Operator #3     |                    | Total           |                 |
|-----------------------------------|-----------------|--------------------|-----------------|--------------------|-----------------|--------------------|-----------------|-----------------|
|                                   | % Positive      | 95 % C.I.          | % Positive      | 95 % C.I.          | % Positive      | 95 % C.I.          | % Positive      | 95 % C.I.       |
| Negative                          | 0%<br>(0/27)    | (0.0%, 12.5%)      | 0%<br>(0/27)    | (0.0%, 12.5%)      | 0%<br>(0/27)    | (0.0%, 12.5%)      | 0%<br>(0/81)    | (0.0%, 4.5%)    |
| Low Positive<br>(3x LOD)          | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(81/81) | (95.5%, 100.0%) |
| Low Positive<br>(5x LOD)          | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(81/81) | (95.5%, 100.0%) |
| Moderate<br>Positive<br>(10x LOD) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(81/81) | (95.5%, 100.0%) |
| High Positive<br>(40x LOD)        | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(27/27) | (87.5%,<br>100.0%) | 100%<br>(81/81) | (95.5%, 100.0%) |

## HIGH DOSE HOOK EFFECT

No high dose hook effect was observed up to 2.8 x 10<sup>5</sup> TCID<sub>50</sub>/mL of gamma-inactivated SARS-CoV-2 with the BD Veritor™ System for Rapid Detection of SARS-CoV-2 test.

## TECHNICAL SUPPORT

Technical Service and Support: In the United States contact BD at 1.800.638.8663 or bd.com.

For regions outside of the United States, contact your local BD representative or bd.com.

## REFERENCES

- Centers for Disease Control and Prevention. <https://www.cdc.gov/coronavirus/2019-ncov/index.html>. Accessed March 30, 2020.
- <https://www.cdc.gov/flu/symptoms/flu-vs-covid19.htm>.
- CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
- CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
- Bullard, et al. Predicting Infectious Severe Acute Respiratory Syndrome Coronavirus 2 From Diagnostic Samples. CID. 2020; Nov 15;71 (10); DOI:10.1093/cid/ciaa638.

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**RÉFÉRENCES**

1. Centers for Disease Control and Prevention. <https://www.cdc.gov/coronavirus/2019-ncov/index.html>. Accessed March 30, 2020.
2. <https://www.cdc.gov/flu/symptoms/flu-vs-covid19.htm>.
3. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
4. CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
5. Bullard, et al. Predicting Infectious Severe Acute Respiratory Syndrome Coronavirus 2 From Diagnostic Samples. CID. 2020; Nov 15;71 (10); DOI:10.1093/cid/ciaa638.

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Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, Maryland 21152 USA

BD, the BD Logo, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2021 BD. All rights reserved.