

# BD Vacutainer® Blood Collection Tubes

The below table lists alternative substitute options by SKU in the event your desired product is not available.

Please follow established laboratory procedures prior to implementing any new blood collection tube.

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1	Sub Option #2
367587	GLUCOSE TUBE	Plastic	13x75	2.0	Hemo	Paper	Sodium Fluoride 3mg, Na2EDTA 6mg	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367921	GLUCOSE TUBE	Plastic	13x75	2.0	Hemo	Paper	Sodium Fluoride 5mg, Potassium Oxalate 4mg	367922, Glucose Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Fluoride 10mg, Potassium Oxalate 8mg	
367922	GLUCOSE TUBE	Plastic	13x75	4.0	Hemo	Paper	Sodium Fluoride 10mg, Potassium Oxalate 8mg	367921, Glucose Tube, Plastic, 13x75, 2.0, Hemo, Paper, Sodium Fluoride 5mg, Potassium Oxalate 4mg	368587, Glucose Tube, Plastic, 13x75, 4.0, Conv, Paper, Sodium Fluoride 10mg, Potassium Oxalate 8mg
368587	GLUCOSE TUBE	Plastic	13x75	4.0	Conv	Paper	Sodium Fluoride 10mg, Potassium Oxalate 8mg	367922, Glucose Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Fluoride 10mg, Potassium Oxalate 8mg	
367925	GLUCOSE TUBE	Plastic	13x100	6.0	Hemo	Paper	Sodium Fluoride 15mg, Potassium Oxalate 12mg	367922, Glucose Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Fluoride 10mg, Potassium Oxalate 8mg	
367729	GLUCOSE TUBE	Glass	13x100	7.0	Hemo	Paper	Sodium Fluoride 30mg	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367001	GLUCOSE TUBE	Glass	16x100	10.0	Conv	Paper	Sodium Fluoride 100mg, Potassium Oxalate 20mg	367925, Glucose Tube, Plastic, 13x100, 6.0, Hemo, Paper, Sodium Fluoride 15mg, Potassium Oxalate 12mg	
367671	SODIUM HEPARIN TUBE	Glass	13x75	2.0	Hemo	Paper	Sodium Heparin 33 USP Units	367871, Sodium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Heparin 75 USP Units	
367871	SODIUM HEPARIN TUBE	Plastic	13x75	4.0	Hemo	Paper	Sodium Heparin 75 USP Units	367878, Sodium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Heparin 75 USP Units	
367878	SODIUM HEPARIN TUBE	Plastic	13x100	6.0	Hemo	Paper	Sodium Heparin 95 USP Units	367871, Sodium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Heparin 75 USP Units	
366664	LITHIUM HEPARIN TUBE	Plastic	13x75	2.0	Hemo	Paper	Lithium Heparin 37 USP Units	366667, Lithium Heparin Tube, Plastic, 13x75, 3.0, Conv, Paper, Lithium Heparin 56 USP Units	
366667	LITHIUM HEPARIN TUBE	Plastic	13x75	3.0	Conv	Paper	Lithium Heparin 56 USP Units	366664, Lithium Heparin Tube, Plastic, 13x75, 2.0, Hemo, Paper, Lithium Heparin 37 USP Units	
366480	SODIUM HEPARIN TUBE	Glass	16x100	10.0	Conv	Paper	Sodium Heparin 158 USP Units	367874, Sodium Heparin Tube, Plastic, 16x100, 10.0, Conv, Paper, Sodium Heparin 158 USP Units	
367874	SODIUM HEPARIN TUBE	Plastic	16x100	10.0	Conv	Paper	Sodium Heparin 158 USP Units	366480, Sodium Heparin Tube, Glass, 16x100, 10.0, Conv, Paper, Sodium Heparin 158 USP Units	
367884	LITHIUM HEPARIN TUBE	Plastic	13x75	4.0	Hemo	Paper	Lithium Heparin 75 USP Units	366667, Lithium Heparin Tube, Plastic, 13x75, 3.0, Conv, Paper, Lithium Heparin 56 USP Units	
367886	LITHIUM HEPARIN TUBE	Plastic	13x100	6.0	Hemo	Paper	Lithium Heparin 95 USP Units	367884, Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units	
367880	LITHIUM HEPARIN TUBE	Plastic	16x100	10.0	Conv	Paper	Lithium Heparin 158 USP Units	367886, Lithium Heparin Tube, Plastic, 13x100, 6.0, Hemo, Paper, Lithium Heparin 95 USP Units	

(continued next page)



## BD Vacutainer® Blood Collection Tubes (continued)

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1	Sub Option #2
366668	SE RUM TUBE	Plastic	13x75	3.0	Conv	Paper	Clot Activator, Silicone Coated	367812, Serum Tube, Plastic, 13x75, 4.0, Hemo, Paper, Clot Activator/ Silicone Coated	
367812	SERUM TUBE	Plastic	13x75	4.0	Hemo	Paper	Clot Activator, Silicone Coated	366668, Serum Tube, Plastic, 13x75, 3.0, Conv, Paper, Clot Activator/ Silicone Coated	
367814	SERUM TUBE	Plastic	13x100	5.0	Hemo	Paper	Clot Activator, Silicone Coated	367815, Serum Tube, Plastic, 13x100, 6.0, Hemo, Paper, Clot Activator/ Silicone Coated	
367815	SERUM TUBE	Plastic	13x100	6.0	Hemo	Paper	Clot Activator, Silicone Coated	367814, Serum Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/ Silicone Coated	
366703	SPECIALTY TUBE	Plastic	13x75	3.0	Hemo	Paper	No Additive	366704, Specialty Tube, Plastic, 13x75, 3.0, Conv, Paper, No Additive	
368660	SERUM TUBE	Plastic	13x100	6.0	Conv	Paper	Clot Activator, Silicone Coated	367815, Serum Tube, Plastic, 13x100, 6.0, Hemo, Paper, Clot Activator/ Silicone Coated	
366430	SERUM TUBE	Glass	16x100	10.0	Conv	Paper	Silicone Coated	367820, Serum Tube, Plastic, 16x100, 10.0, Conv, Paper, Clot Activator/ Silicone Coated	
367820	SERUM TUBE	Plastic	16x100	10.0	Conv	Paper	Clot Activator, Silicone Coated	367815, Serum Tube, Plastic, 13x100, 6.0, Hemo, Paper, Clot Activator/ Silicone Coated	
366704	SPECIALTY TUBE	Plastic	13x75	3.0	Conv	Paper	No Additive	366703, Specialty Tube, Plastic, 13x75, 3.0, Hemo, Paper, No Additive	
364606	SPECIALTY TUBE	Glass	16x100	8.5	Conv	Paper	Acid Citric Dextrose (ACD) Solution A consists of trisodium citrate, 22.0 g/L, Citric Acid, 8.0 g/L, and Dextrose, 24.5 g/L	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
364816	SPECIALTY TUBE	Glass	13x100	6.0	Conv	Paper	Acid Citric Dextrose (ACD) Solution B consists of trisodium citrate, 13.2 g/L, Citric Acid, 4.8 g/L, and Dextrose, 14.7 g/L	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366401	SPECIALTY TUBE	Glass	13x100, 16x100	7.0 & 10.0	Conv	Paper	Silicone Coated K3EDTA (12.15mg). Packaged in a sterile pouch.	366643, EDTA Tube, 16x100, 10.0, Hemo, See Through, K2EDTA 18mg	
364960	SPECIALTY TUBE	Glass	16x100	8.3	Conv	See Through	Sodium Polyanethol Sulfonate (SPS) (0.35% in 0.85% Sodium Chloride) 1.7mL	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
368380	SPECIALTY TUBE	Plastic	13x100	6.0	Hemo	Paper	Serum Clot Activator (Silicone Coated) for Trace Element	368381, Specialty Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg for Trace Element	
368381	SPECIALTY TUBE	Plastic	13x100	6.0	Hemo	Paper	K2EDTA 10.8mg for Trace Element	368380, Specialty Tube, Plastic, 13x100, 6.0, Hemo, Paper, Serum Clot Activator/Silicone Coated for Trace Element	
368774	RST TUBE	Plastic	13x100	5.0	Hemo	Paper	Thrombin-Based Clot Activator	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/ Polymer Gel	
367842	EDTA TUBE	Plastic	13x75	2.0	Hemo	Block Label	K2EDTA 3.6mg	367841, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Paper, K2EDTA 3.6mg	367835, EDTA Tube, Plastic, 13x75, 3.0, Conv, Paper, K2EDTA 5.4mg
367841	EDTA TUBE	Plastic	13x75	2.0	Hemo	Paper	K2EDTA 3.6mg	367842, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Block Label, K2EDTA 3.6mg	367835, EDTA Tube, Plastic, 13x75, 3.0, Conv, Paper, K2EDTA 5.4mg
367835	EDTA TUBE	Plastic	13x75	3.0	Conv	Paper	K2EDTA 5.4mg	367841, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Paper, K2EDTA 3.6mg	367842, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Block Label, K2EDTA 3.6mg

(continued next page)



## BD Vacutainer® Blood Collection Tubes (continued)

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1	Sub Option #2
368056	PST TUBE	Plastic	13x75	3.0	Conv	Paper	Lithium Heparin 56 USP Units w/ Polymer Gel	367960, PST Tube, Plastic, 13x75, 3.0, Hemo, Paper, Lithium Heparin 56 USP Units w/ Polymer Gel	366667, Lithium Heparin Tube, Plastic, 13x75, 3.0, Conv, Paper, Lithium Heparin 56 USP Units
367960	PST TUBE	Plastic	13x75	3.0	Hemo	Paper	Lithium Heparin 56 USP Units w/ Polymer Gel	368056, PST Tube, Plastic, 13x75, 3.0, Conv, Paper, Lithium Heparin 56 USP Units w/ Polymer Gel	367884 Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units
367856	EDTA TUBE	Plastic	13x75	3.0	Hemo	Paper	K2EDTA 5.4mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg	367841, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Paper, K2EDTA 3.6mg
367981	SST TUBE	Plastic	13x75	3.5	Conv	Paper	Clot Activator/Polymer Gel	367983, SST Tube, Plastic, 13x75, 3.5, Hemo, Paper, Clot Activator/Polymer Gel	
367961	PST TUBE	Plastic	13x100	3.5	Hemo	Paper	Lithium Heparin 65 USP Units w/ Polymer Gel	367962, PST Tube, Plastic, 13x100, 4.5, Hemo, Paper, Lithium Heparin 84 USP Units w/ Polymer Gel	367884 Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units
367983	SST TUBE	Plastic	13x75	3.5	Hemo	Paper	Clot Activator/Polymer Gel	367981, SST Tube, Plastic, 13x75, 3.5, Conv, Paper, Clot Activator/Polymer Gel	
367862	EDTA TUBE	Plastic	13x75	4.0	Hemo	See Through	K2EDTA 7.2mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg	
367861	EDTA TUBE	Plastic	13x75	4.0	Hemo	Paper	K2EDTA 7.2mg	367844, EDTA Tube, Plastic, 13x75, 4.0, Conv, Paper, K2EDTA 7.2mg	367856, EDTA Tube, Plastic, 13x75, 3.0, Hemo, Paper, K2EDTA 5.4mg
367977	SST TUBE	Plastic	13x100	4.0	Hemo	Paper	Clot Activator/Polymer Gel	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/Polymer Gel	367983, SST Tube, Plastic, 13x75, 3.5, Hemo, Paper, Clot Activator/Polymer Gel
367844	EDTA TUBE	Plastic	13x75	4.0	Conv	Paper	K2EDTA 7.2mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg	368661, EDTA Tube, Plastic, 13x100, 6.0, Conv, Paper, K2EDTA 10.8mg
368661	EDTA TUBE	Plastic	13x100	6.0	Conv	Paper	K2EDTA 10.8mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg	
367989	SST TUBE	Plastic	13x100	5.0	Hemo	See Through	Clot Activator/Polymer Gel	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/Polymer Gel	
363080	CITRATE TUBE	Plastic	13x75	1.8	Hemo	Paper	Buffered Sodium Citrate (0.109M, 3.2%)	363083, Citrate Tube, Plastic, 13x75, 2.7, Hemo, Paper, Buffered Sodium Citrate (0.109M, 3.2%)	
363083	CITRATE TUBE	Plastic	13x75	2.7	Hemo	Paper	Buffered Sodium Citrate (0.109M, 3.2%)	363080, Citrate Tube, Plastic, 13x75, 1.8, Hemo, Paper, Buffered Sodium Citrate (0.109M, 3.2%)	
367899	EDTA TUBE	Plastic	13x100	6.0	Hemo	Block Label	K2EDTA 10.8mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg
367963	PST TUBE	Plastic	13x100	4.5	Hemo	See Through	Lithium Heparin 84 USP Units w/ Polymer Gel	367962, PST Tube, Plastic, 13x100, 4.5, Hemo, Paper, Lithium Heparin 84 USP Units w/ Polymer Gel	367886, Lithium Heparin Tube, Plastic, 13x100, 6.0, Hemo, Paper, Lithium Heparin 95 USP Units
367986	SST TUBE	Plastic	13x100	5.0	Hemo	Paper	Clot Activator/Polymer Gel	367977, SST Tube, Plastic, 13x100, 4.0, Hemo, Paper, Clot Activator/Polymer Gel	367989, SST Tube, Plastic, 13x100, 5.0, Hemo, See Through, Clot Activator/Polymer Gel

(continued next page)

## BD Vacutainer® Blood Collection Tubes (continued)

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1	Sub Option #2
367863	EDTA TUBE	Plastic	13x100	6.0	Hemo	Paper	K2EDTA 10.8mg	367899, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Block Label, K2EDTA 10.8mg	368661, EDTA Tube, Plastic, 13x100, 6.0, Conv, Paper, K2EDTA 10.8mg
367962	PST TUBE	Plastic	13x100	4.5	Hemo	Paper	Lithium Heparin 84 USP Units w/ Polymer Gel	367961, PST Tube, Plastic, 13x100, 3.5, Hemo, Paper, Lithium Heparin 65 USP Units of Polymer Gel	367884 Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units
366450	EDTA TUBE	Glass	13x100	7.0	Conv	Paper	K3EDTA 12.15mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg	
367987	SST TUBE	Plastic	16x100	7.5	Conv	Paper	Transport Tube Clot Activator/ Double Polymer Gel	367988, SST Tube, Plastic, 16x100, 8.5, Conv, Paper, Clot Activator/ Polymer Gel	367985, SST Tube, Plastic, 16x125, 10.0, Conv, Paper, Transport Tube Clot Activator/ Double Polymer Gel
368589	EDTA TUBE	Plastic	16x100	10.0	Conv	Block Label	K2EDTA 18mg	367899, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Block Label, K2EDTA 10.8mg	366643, EDTA Tube, Plastic, 16x100, 10.0, Hemo, See Through, K2EDTA 18mg
367988	SST TUBE	Plastic	16x100	8.5	Conv	Paper	Clot Activator/Polymer Gel	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/ Polymer Gel	
366643	EDTA TUBE	Plastic	16x100	10.0	Hemo	See Through	K2EDTA 18mg	368589, EDTA Tube, Plastic, 16x100, 10.0, Conv, Block Label, K2EDTA 18mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg
367964	PST TUBE	Plastic	16x100	8.0	Conv	Paper	Lithium Heparin 126 USP Units w/ Polymer Gel	367962, PST Tube, Plastic, 13x100, 4.5, Hemo, Paper, Lithium Heparin 84 USP Units w/ Polymer Gel	367880 Lithium Heparin Tube, Plastic, 16x100, 10.0, Conv, Paper, Lithium Heparin 158 USP Units
367985	SST TUBE	Plastic	16x125	10.0	Conv	Paper	Transport Tube Clot Activator/ Double Polymer Gel	367988, SST Tube, Plastic, 16x100, 8.5, Conv, Paper, Clot Activator/ Polymer Gel	367987, SST Tube, Plastic, 16x100, 7.5, Conv, Paper, Transport Tube Clot Activator/ Double Polymer Gel
369741	SPECIALTY TUBE	Glass	13x75	2.4	Hemo	Paper	Buffered Citrate (32.0 mg sodium citrate, 4.2 mg Citric Acid/mL) 0.6mL	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366408	SPECIALTY TUBE	Plastic	13x100	6	Hemo	Paper	No Additive	366703, Specialty Tube, Plastic, 13x75, 3.0, Hemo, Paper, No Additive	

BD Life Sciences, 7 Loveton Circle, Sparks, MD, 21152, USA  
800.638.8663

[bd.com](http://bd.com)

