

BD Vacutainer® Blood Collection Tubes

Last Updated April 2022

SKU Listing with alternative substitute options in the event that a desired product is not available.

It is recommended that customers refer to their institution's policies and procedures for product validation prior to use.

 = Temporarily halted production, must select an alternative

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1*	Sub Option #2*
363080	CITRATE TUBE	Plastic	13x75	1.8	Hemo	Paper	Buffered Sodium Citrate (0.109M, 3.2%)	363083, Citrate Tube, Plastic, 13x75, 2.7, Hemo, Paper, Buffered Sodium Citrate (0.109M, 3.2%)	
363083	CITRATE TUBE	Plastic	13x75	2.7	Hemo	Paper	Buffered Sodium Citrate (0.109M, 3.2%)	363080, Citrate Tube, Plastic, 13x75, 1.8, Hemo, Paper, Buffered Sodium Citrate (0.109M, 3.2%)	
367835	EDTA TUBE	Plastic	13x75	3.0	Conv	Paper	K2EDTA 5.4mg	367841, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Paper, K2EDTA 3.6mg	367842, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Block Label, K2EDTA 3.6mg
367842	EDTA TUBE	Plastic	13x75	2.0	Hemo	Block Label	K2EDTA 3.6mg	367841, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Paper, K2EDTA 3.6mg	
367841	EDTA TUBE	Plastic	13x75	2.0	Hemo	Paper	K2EDTA 3.6mg	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367844	EDTA TUBE	Plastic	13x75	4.0	Conv	Paper	K2EDTA 7.2mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg	
367856	EDTA TUBE	Plastic	13x75	3.0	Hemo	Paper	K2EDTA 5.4mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg	367841, EDTA Tube, Plastic, 13x75, 2.0, Hemo, Paper, K2EDTA 3.6mg
367862	EDTA TUBE	Plastic	13x75	4.0	Hemo	See Through	K2EDTA 7.2mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg	
367861	EDTA TUBE	Plastic	13x75	4.0	Hemo	Paper	K2EDTA 7.2mg	367856, EDTA Tube, Plastic, 13x75, 3.0, Hemo, Paper, K2EDTA 5.4mg	
368661	EDTA TUBE	Plastic	13x100	6.0	Conv	Paper	K2EDTA 10.8mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg	
367899	EDTA TUBE	Plastic	13x100	6.0	Hemo	Block Label	K2EDTA 10.8mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg	367861, EDTA Tube, Plastic, 13x75, 4.0, Hemo, Paper, K2EDTA 7.2mg
367863	EDTA TUBE	Plastic	13x100	6.0	Hemo	Paper	K2EDTA 10.8mg	367899, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Block Label, K2EDTA 10.8mg	
366450	EDTA TUBE	Glass	13x100	7.0	Conv	Paper	K3EDTA 12.15mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg	
368589	EDTA TUBE	Plastic	16x100	10.0	Conv	Block Label	K2EDTA 18mg	367899, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Block Label, K2EDTA 10.8mg	366643, EDTA Tube, Plastic, 16x100, 10.0, Hemo, See Through, K2EDTA 18mg
366643	EDTA TUBE	Plastic	16x100	10.0	Hemo	See Through	K2EDTA 18mg	367863, EDTA Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg	
367587	GLUCOSE TUBE	Plastic	13x75	2.0	Hemo	Paper	Sodium Fluoride 3mg, Na2EDTA 6mg	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
368587	GLUCOSE TUBE	Plastic	13x75	4.0	Conv	Paper	Sodium Fluoride 10mg, Potassium Oxalate 8mg	367922, Glucose Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Fluoride 10mg, Potassium Oxalate 8mg	

(continued next page)



BD Vacutainer® Blood Collection Tubes

 = Temporarily halted production, must select an alternative

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1*	Sub Option #2*
367921	GLUCOSE TUBE	Plastic	13x75	2.0	Hemo	Paper	Sodium Fluoride 5mg, Potassium Oxalate 4mg	367922, Glucose Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Fluoride 10mg, Potassium Oxalate 8mg	
367922	GLUCOSE TUBE	Plastic	13x75	4.0	Hemo	Paper	Sodium Fluoride 10mg, Potassium Oxalate 8mg	367921, Glucose Tube, Plastic, 13x75, 2.0, Hemo, Paper, Sodium Fluoride 5mg, Potassium Oxalate 4mg	
367925	GLUCOSE TUBE	Plastic	13x100	6.0	Hemo	Paper	Sodium Fluoride 15mg, Potassium Oxalate 12mg	367922, Glucose Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Fluoride 10mg, Potassium Oxalate 8mg	
367729	GLUCOSE TUBE	Glass	13x100	7.0	Hemo	Paper	Sodium Fluoride 30mg	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366664	LITHIUM HEPARIN TUBE	Plastic	13x75	2.0	Hemo	Paper	Lithium Heparin 37 USP Units	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366667	LITHIUM HEPARIN TUBE	Plastic	13x75	3.0	Conv	Paper	Lithium Heparin 56 USP Units	366664, Lithium Heparin Tube, Plastic, 13x75, 2.0, Hemo Conv, Paper, Lithium Heparin 37 USP Units	
367884	LITHIUM HEPARIN TUBE	Plastic	13x75	4.0	Hemo	Paper	Lithium Heparin 75 USP Units	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367886	LITHIUM HEPARIN TUBE	Plastic	13x100	6.0	Hemo	Paper	Lithium Heparin 95 USP Units	367884, Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units	
367880	LITHIUM HEPARIN TUBE	Plastic	16x100	10.0	Conv	Paper	Lithium Heparin 158 USP Units	367884, Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units	
368056	PST TUBE	Plastic	13x75	3.0	Conv	Paper	Lithium Heparin 56 USP Units w/ Polymer Gel	367960, PST Tube, Plastic, 13x75, 3.0, Hemo, Paper, Lithium Heparin 56 USP Units w/ Polymer Gel	
367960	PST TUBE	Plastic	13x75	3.0	Hemo	Paper	Lithium Heparin 56 USP Units w/ Polymer Gel	367884 Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units	
367961	PST TUBE	Plastic	13x100	3.5	Hemo	Paper	Lithium Heparin 65 USP Units w/ Polymer Gel	367962, PST Tube, Plastic, 13x100, 4.5, Hemo, Paper, Lithium Heparin 84 USP Units w/ Polymer Gel	367884 Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units
367963	PST TUBE	Plastic	13x100	4.5	Hemo	See Through	Lithium Heparin 84 USP Units w/ Polymer Gel	367962, PST Tube, Plastic, 13x100, 4.5, Hemo, Paper, Lithium Heparin 84 USP Units w/ Polymer Gel	
367962	PST TUBE	Plastic	13x100	4.5	Hemo	Paper	Lithium Heparin 84 USP Units w/ Polymer Gel	367884 Lithium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Lithium Heparin 75 USP Units	
367964	PST TUBE	Plastic	16x100	8.0	Conv	Paper	Lithium Heparin 126 USP Units w/ Polymer Gel	367962, PST Tube, Plastic, 13x100, 4.5, Hemo, Paper, Lithium Heparin 84 USP Units w/ Polymer Gel	
368774	RST TUBE	Plastic	13x100	5.0	Hemo	Paper	Thrombin-Based Clot Activator	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/ Polymer Gel	
366668	SERUM TUBE	Plastic	13x75	3.0	Conv	Paper	Clot Activator, Silicone Coated	367812, Serum Tube, Plastic, 13x75, 4.0, Hemo, Paper, Clot Activator/ Silicone Coated	
367812	SERUM TUBE	Plastic	13x75	4.0	Hemo	Paper	Clot Activator, Silicone Coated	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367814	SERUM TUBE	Plastic	13x100	5.0	Hemo	Paper	Clot Activator, Silicone Coated	367815, Serum Tube, Plastic, 13x100, 6.0, Hemo, Paper, Clot Activator/ Silicone Coated	
367815	SERUM TUBE	Plastic	13x100	6.0	Hemo	Paper	Clot Activator, Silicone Coated	367812, Serum Tube, Plastic, 13x75, 4.0, Hemo, Paper, Clot Activator/ Silicone Coated	
368660	SERUM TUBE	Plastic	13x100	6.0	Conv	Paper	Clot Activator, Silicone Coated	367815, Serum Tube, Plastic, 13x100, 6.0, Hemo, Paper, Clot Activator/ Silicone Coated	
366430	SERUM TUBE	Glass	16x100	10.0	Conv	Paper	Silicone Coated	367820, Serum Tube, Plastic, 16x100, 10.0, Conv, Paper, Clot Activator/ Silicone Coated	

(continued next page)



BD Vacutainer® Blood Collection Tubes

■ = Temporarily halted production, must select an alternative

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1*	Sub Option #2*
367820	SERUM TUBE	Plastic	16x100	10.0	Conv	Paper	Clot Activator, Silicone Coated	367815, Serum Tube, Plastic, 13x100, 6.0, Hemo, Paper, Clot Activator/ Silicone Coated	
367671	SODIUM HEPARIN TUBE	Glass	13x75	2.0	Hemo	Paper	Sodium Heparin 33 USP Units	367871, Sodium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Heparin 75 USP Units	
367871	SODIUM HEPARIN TUBE	Plastic	13x75	4.0	Hemo	Paper	Sodium Heparin 75 USP Units	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367878	SODIUM HEPARIN TUBE	Plastic	13x100	6.0	Hemo	Paper	Sodium Heparin 95 USP Units	367871, Sodium Heparin Tube, Plastic, 13x75, 4.0, Hemo, Paper, Sodium Heparin 75 USP Units	
366480	SODIUM HEPARIN TUBE	Glass	16x100	10.0	Conv	Paper	Sodium Heparin 158 USP Units	367874, Sodium Heparin Tube, Plastic, 16x100, 10.0, Conv, Paper, Sodium Heparin 158 USP Units	
367874	SODIUM HEPARIN TUBE	Plastic	16x100	10.0	Conv	Paper	Sodium Heparin 158 USP Units	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367001	SPECIALTY TUBE	Glass	16x100	10.0	Conv	Paper	Sodium Fluoride 100mg, Potassium Oxalate 20mg	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366703	SPECIALTY TUBE	Plastic	13x75	3.0	Hemo	Paper	No Additive	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366704	SPECIALTY TUBE	Plastic	13x75	3.0	Conv	Paper	No Additive	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
364606	SPECIALTY TUBE	Glass	16x100	8.5	Conv	Paper	Acid Citric Dextrose (ACD) Solution A consists of trisodium citrate, 22.0 g/L, Citric Acid, 8.0 g/L, and Dextrose, 24.5 g/L	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
364816	SPECIALTY TUBE	Glass	13x100	6.0	Conv	Paper	Acid Citric Dextrose (ACD) Solution B consists of trisodium citrate, 13.2 g/L, Citric Acid, 4.8 g/L, and Dextrose, 14.7 g/L	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366401	SPECIALTY TUBE	Glass	13x100, 16x100	7.0 & 10.0	Conv	Paper	Silicone Coated K3EDTA (12.15mg). Packaged in a sterile pouch.	366643, EDTA Tube, 16x100, 10.0, Hemo, See Through, K2EDTA 18mg	
364960	SPECIALTY TUBE	Glass	16x100	8.3	Conv	See Through	Sodium Polyanethol Sulfonate (SPS) (0.35% in 0.85% Sodium Chloride) 1.7mL	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
368380	SPECIALTY TUBE	Plastic	13x100	6.0	Hemo	Paper	Serum Clot Activator (Silicone Coated) for Trace Element	368381, Specialty Tube, Plastic, 13x100, 6.0, Hemo, Paper, K2EDTA 10.8mg for Trace Element	
368381	SPECIALTY TUBE	Plastic	13x100	6.0	Hemo	Paper	K2EDTA 10.8mg for Trace Element	368380, Specialty Tube, Plastic, 13x100, 6.0, Hemo, Paper, Serum Clot Activator/Silicone Coated for Trace Element	
369741	SPECIALTY TUBE	Glass	13x75	2.4	Hemo	Paper	Buffered Citrate (32.0 mg sodium citrate, 4.2 mg Citric Acid/mL) 0.6mL	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
366408	SPECIALTY TUBE	Plastic	13x100	6	Hemo	Paper	No Additive	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	
367983	SST TUBE	Plastic	13x75	3.5	Hemo	Paper	Clot Activator/Polymer Gel	Please contact BD Technical Services for the best sub option. Call 1.800.631.0174	

(continued next page)



BD Vacutainer® Blood Collection Tubes

 = Temporarily halted production, must select an alternative

SKU	Family	Glass/ Plastic	Tube Size	Draw Vol (mL)	Closure	Label Type	Additive	Sub Option #1*	Sub Option #2*
367977	SST TUBE	Plastic	13x100	4.0	Hemo	Paper	Clot Activator/Polymer Gel	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/ Polymer Gel	367983, SST Tube, Plastic, 13x75, 3.5, Hemo, Paper, Clot Activator/Polymer Gel
367981	SST TUBE	Plastic	13x75	3.5	Conv	Paper	Clot Activator/Polymer Gel	367983, SST Tube, Plastic, 13x75, 3.5, Hemo, Paper, Clot Activator/ Polymer Gel	
367989	SST TUBE	Plastic	13x100	5.0	Hemo	See Through	Clot Activator/Polymer Gel	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/ Polymer Gel	
367986	SST TUBE	Plastic	13x100	5.0	Hemo	Paper	Clot Activator/Polymer Gel	367988, SST Tube, Plastic, 16x100, 8.5, Conv, Paper, Clot Activator/ Polymer Gel	
367987	SST TUBE	Plastic	16x100	7.5	Conv	Paper	Transport Tube Clot Activator/ Double Polymer Gel	367988, SST Tube, Plastic, 16x100, 8.5, Conv, Paper, Clot Activator/ Polymer Gel	
367988	SST TUBE	Plastic	16x100	8.5	Conv	Paper	Clot Activator/Polymer Gel	367986, SST Tube, Plastic, 13x100, 5.0, Hemo, Paper, Clot Activator/ Polymer Gel	
367985	SST TUBE	Plastic	16x125	10.0	Conv	Paper	Transport Tube Clot Activator/ Double Polymer Gel	367988, SST Tube, Plastic, 16x100, 8.5, Conv, Paper, Clot Activator/ Polymer Gel	

*The potential alternative products are BD's suggestions based on its own comparison of the form, features and functions of the original products and the potential alternative products. Depending on a customer's circumstances, clinical requirements and intended use of the proposed alternative product, such product may not be a clinically appropriate alternative to the original product, and any such determination shall be made by the customer.

BD Life Sciences, 7 Loveton Circle, Sparks, MD, 21152, USA
800.638.8663

bd.com

