

There's a life behind every test

Simple **RSV** testing for those at highest risk

The portable, easy-to-use **BD Veritor™ Plus System** provides reliable RSV results in 10 minutes.



Be prepared for respiratory illness

RSV

When your pediatric patients[†] are in need of fast, reliable RSV testing, turn to the **BD Veritor™ Plus System**. Offering reliable results at point of care in a simple-to-operate, handheld instrument.



Simplifies the testing process

- Easy operation and single-button functionality may help reduce potential manual errors
- Provides intuitive sample processing with pre-filled, unitized tubes



Delivers workflow efficiency

- Adapts easily to your workflow by offering 2 operational modes
 - » **Walk Away:** the test device is inserted into the Analyzer immediately after preparation, allowing staff to multitask while the sample incubates
 - » **Analyze Now:** the test device is inserted after the incubation time is complete, allowing batches of samples to be tested at once
- Option to test the same sample with a BD Veritor™ Flu A+B test device.¹



Achieves reliable, rapid results

- Displays easy-to-read digital results for RSV in 10 minutes¹
- Provides clear, objective results by correcting for nonspecific binding and detecting positives not recognized by the unaided eye¹



Provides result traceability

- Downloadable test records matched to kit lot number, operator, patient or specimen ID with the BD Veritor™ InfoWiFi Module or via the companion ImageMover application
- Offers result printing capabilities via USB port

¹The BD Veritor™ RSV point-of-care kit is indicated for use in patients under 6 years of age. The BD Veritor™ RSV laboratory kit is indicated for use in patients under 20 years of age.

Reference: 1. BD Veritor™ System for Rapid Detection of Respiratory Syncytial Virus (RSV), package insert, 8086098. Franklin Lakes, NJ: Becton, Dickinson and Company

BD Life Sciences, 2100 Derry Road West, Mississauga, ON L5N 0B3 CA
866-979-9408

BDVeritor.com

BD, the BD Logo and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved.

Offers proven performance (point-of-care kit)

	PPA [†]	NPA [‡]	Results in
Compared to viral culture	91.8% (95% CI: 85.9%, 95.4%)	93.3% (95% CI: 90.4%, 94.4%)	10
Compared to PCR	81.6% (95% CI: 75.2%, 86.6%)	99.1% (95% CI: 97.5%, 99.7%)	10

[†]positive percent agreement
[‡]negative percent agreement

Ordering information	Cat no.	Qty
BD Veritor™ Plus System RSV kit	256038	30 tests
BD Veritor™ Plus System RSV laboratory kit	256042	30 tests
BD Veritor™ Plus Analyzer	256066	1 unit

*BD Veritor™ System for Rapid Detection of SARS-CoV-2:

- This test has been authorized for sale in Canada by Health Canada under Interim Order;
- This test is only authorized for the duration of the Interim Order Respecting the Importation and Sale of Medical Devices for Use in Relation to COVID-19, unless the authorization is terminated or revoked sooner;
- This test has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.

